

# <sup>13</sup>C Baseline Correction on Probes with a High Q-Factor


Written by Markus Voehler, January 17, 2013

All our cryoprobes have cold <sup>13</sup>C preamplifiers, which means they have the highest possible sensitivity for direct detection with a very high Q-factor. Because of that, baseline distortion due to first point imperfections are common and need to be accounted for.

You can either do this with an AU-program or manually. Both ways are described in this tutorial.

## 1) Using the AU-program

- a) Go to *ProcPars* → *Automation*, and set the AU-program “c13cryo” in the AUNMP field



- b) In the command line type “xaup”. This program will modify your initial FID in the current expno, but save the original FID in the expno+100,000.

## 2) Manual adjustments

- a) Change NZP = 4  
b) Execute ZP  
c) Execute CONVDTA → this will ask for a new expno to modify the FID in  
d) Change the following Processing Parameters in ProcPars  
i) ME\_mod LPbc  
ii) NCOEFF 32  
iii) TDoff 32  
e) Do an efp to process the spectrum in the new expno.