

# ***Solvent Suppression***

# Why solvent suppression?



**There are two major problems arising from an intense solvent peak:**

- The un-deuterated solvent contributes the major intensity of the recorded FID. Therefore the receiver gain must be lowered and becomes very insensitive. Consequently the small signal contributions of the molecule of interest are monitored imperfectly.  
(Dynamic range problem)
- The solvent peak becomes very broad and overlays with peaks arising from the molecule of interest.

# Strategies for different methods



- **Selective dephasing of the solvent resonance**

- Purge spin-lock pulses
- WATERGATE, excitation sculpting, WET

- **Realigning solvent magnetization along z-axis**

- Water flip back (often together with WATERGATE)
- Jump and Return
- Manipulation of radiation damping



Radiation damping = Exotisch

- **Other strategies**

- Presaturation (Chemical shift reagents will also be presaturated)
- Coherence selection (e.g. gradient HSQC with EA)

ist auch dephasing aber nicht selective dephasing, gehört also so bisschen zu i

# Presaturation



- On- or off-resonance presaturation during recycle or other delays (e.g. during NOE - mixing time  $d8$  in `noesy1dpr`)
- Problem: saturation transfer to exchanging protons
- **zgpr**: simple one solvent presaturation:  
power: **pl9** (55 - 65dB, calculate in *edprosol*) , duration: **d1** (2s)  
offset: on-resonant on solvent peak.  
(Optimize **o1** on minimal FID in `gs` mode)
- **noesypr1d**: 1D NOESY presaturation for single solvent suppression:  
power: **pl9**, duration: **d1** and **d8** (NOESY mixing time)  
offset: on-resonant on solvent peak.  
(Optimize **o1** on minimal FID in `gs` mode)

# Gradient coherence selection



**Trick:** only specific coherences are selected, all others are dephased

**Problem:** exchanging protons are also dephased

- **Pulse programs:** large number of 2D pulse programs (e.g. hsqcetgpsi, cosygpmfph etc.)
- Optimization:
  - Use exact gradient ratios based on gyromagnetic ratios (e.g.: 80 : 8 instead of 80 : 8.1 may lead to 30% sensitivity loss)
  - Use lowest possible gradient power to avoid sensitivity losses due to diffusion effects
  - For samples with short T2 relaxation times, use short gradient pulses
  - Avoid static gradients due to large z5 or z6 shim values. They could act as refocussing gradients and compromise the solvent suppression.
    - in hsqcetgpsi, use 80 : 8.1 : -80 : 8.1 or -80 : -8.1 : 80 : -8.1
    - optimize values of z5 and z6 in the gs mode.

kti:

Artifacts:

Static gradient - which could be created by the usage of large values for high-order shims (Z5, Z6) could act as refocussing gradients and thereby compromise the quality of solvent suppression. Some tips to solve this problem,

- deconvolution of the solvent signal
- removes all signals within the defined range around the solvent frequency

e.g. for `1m1etgss1`, use gradient ratios `80 : 8.1 : -80 : 8.1` instead of

`80 : 8.1 : 80 : -8.1`


Parameters

- `bc_mod` = `qpol` subtracts a 5th order polynomial from the FID
- `bc_mod` = `qfil` applies a filter according to

3. A shim perfect for presaturation might not be the best shim for a gradient solvent suppression experiment. Change Z5, Z6 within the GS-mode and observe the changes in the quality of water suppression.

`bcfw` range for deconvolution, in ppm  
`corroffs` offset from spectrum midpoint, in Hz

# Selective dephasing of solvent signal by rf pulses



## Purge pulses (a.k.a spin-lock- or trim pulses)

selectively dephase coherences orthogonal to the RF-field ( $I_x, I_z$ ) while preserving coherences locked along the RF-field ( $I_y$ )

### pulse program

many HSQC-type inverse and triple resonance experiments.

### purge pulse

**p28** (500 $\mu$ sec - 2ms)

### power level for purge pulse

**p11** (same as hard proton pulses)

**Note: make sure that the probe can stand the rather high power level used for the trim pulse. Otherwise the probe is damaged.**


### Optimizing conditions:

*Change the trim pulse p28 within the GS-mode for minimum solvent signal.*

### Warning:

*When a simple HSQC experiment does not give the expected sensitivity or even no*


# Selective dephasing of solvent signal by rf pulses



Purge pulses (a.k.a spin-lock- or trim pulses)


$p28 = 500\mu\text{s} - 2\text{ms}$  @  $p11$  !  
use  $p28$  only if you have a solvent signal to suppress

Example: hsqcetgpsi





# Selective dephasing of solvent signal by gradients



selectively dephase solvent magnetization

dephasing gradients are also called purge or crush gradients

Pulse programs

- many inverse and triple resonance experiments
- Schemes: Watergate, WET, Excitation Sculpting


Optimization

- use lowest possible gradient power strength
- in Watergate use shortest possible gradients to minimize J-coupling artifacts.

Note:

dephasing gradients (or dephasing schemes, watergate etc.) can be preceded by flip-back pulses to reduce the gradient strength

# WATERGATE



## WATER suppression by GrAdent Tailored Excitation


selectively dephase solvent magnetization by Gradients

<b>pulse program</b>	<b>p3919gp</b> <b>19 or wg</b>	(1D-NMR with binomial WATERGATE) (nD-experiments using WATERGATE)
<b>watergate pulse (on-resonant)</b>	<b>p27</b>	(= hard pulse or $\sim 35\mu\text{sec}$ )
<b>power level for selective pulse</b>	<b>pl18</b>	(= value of hard pulse or $\sim 35\mu\text{sec}$ pulse)
<b>delay for binomial suppression</b>	<b>d19</b>	(2 * dwell time, 2 * <b>dw</b> , depends on <b>TD</b> and <b>SW</b> )

Optimizing conditions:


- Use lowest gradient power level as possible.
- Keep the length of the gradient pulse as short as possible to

# WATERGATE



## Original WATERGATE Sequence:

- Based on gradient spin echo experiment to refocus chemical shift evolution.
- Coherence orders of water and other signals are changed differently:
  - On resonant coherence unchanged
  - Off resonant coherence inverted
- Gradients discriminate on- and off-resonance coherences by CTP selection.





$$E_{tot1} = +20 (\gamma \cdot 1) = +20$$

$$E_{tot2} = +20 (\gamma \cdot -1) = -20$$

0


## WATERGATE properties:

- Excitation profile: Notch (band stop-) filter
- in **p3919gp** a binomial pulse sequence is used which is enclosed by the gradients: (no further pulse adjustments)
  - The binomial pulse sequence changes the coherence order of all resonances except that of the on-resonant solvent.
  - The second field gradient dephases the water signal and rephases all other resonances.




- Gradients and selective pulses: *BEHIND* excitation
  - ⇒ total length of pulse sequence is 4-6ms
  - ⇒ evolution of  $^3J(\text{HH})$  scalar couplings leads to phase errors

# WATERGATE phase errors



Evolution of  $^3J(\text{HH})$  scalar couplings leads to phase errors



# Water-flip-back



Common building block in heteronuclear correlation experiments to improve WATERGATE

- flip-back pulses prior to WATERGATE to align transverse water magnetization along z
- Pulse sequence names: *...fp...* with the following parameters:


**selective on-resonant flip-back pulse:** **p11** (2 - 3ms the shorter the broader)

**shape of the flip-back pulse:** **spnam1** (sinc or gauss pulse)

**power of the flip-back pulse:** **sp1** (calculate selective 90° pulse in

HSQC with water-flip-back improved WATERGATE: **stdisp** and optimize in **gs-mode**)

**invifpf3gp phwg**



# Excitation Sculpting



same principle as WATERGATE:

- selectively dephase the solvent magnetization

Difference:

- two identical spin echoes flanked by gradients (instead of one spin echo)
- the flanking gradients are of different strength

**pulse program**

**zgesgp**

(1D-NMR with selective pulses) or

**zggpw5**

*binomial W5 pulse train*  $\equiv$  *binomial WATERGATE*)

**es**

(nD-experiments)

**Parameters for **zgesgp** and **...es...**:**

**selective 180° pulse (on-resonant)**

**p12**

(1 - 5ms, the shorter the broader )

**power level for selective pulse**

**sp1**

(calculate with “*pulse*” for **p12** length subtract 6dB for 180° pulse)

**shape of selective 180° pulse**

**spnam1**



(Squa100.1000)

Optimization:

# Excitation Sculpting: no phase distortion




zgesgp





# Multiple Solvent Suppression



- **zgps:** multiple solvent or off resonant presaturation:  
power: **sp6** (= **p19** +6dB per additional solvent),  
duration: **l6** \* **p18** (**p18** = 10ms, **l6** = 200)  
**spnam6:** square (manipulate with solvent freqlist in *stdisp*)



- **zgesgp:** multiple solvent suppression using excitation sculpting  
**p12** (1 - 5ms)  
**spnam1** (square, manipulate with solvent freqlist in *stdisp*)  
**sp1** (calculate  $90^\circ$  in pulse, subtract 6dB to get  $180^\circ$  and subtract 6dB for each additional signal that has to be suppressed).

# Shapes for Multiple Solvent Suppression




- 1 phase
- calibrate
- integrate
- utilities
- dual

utilities



Please enter type of list (f1, f2, f3) : f1  
Please enter name of f1 list: frqlist  
Write name of f1 list to acqu parameters ?

2  
frqlist



4 start the shape tool  
type: **stdisp**

# Shapes for Multiple Solvent Suppression



5

select Shapes → e.g. Rectangle

**Editor for Shape Parameters**

Parameters for Shape Rectangle:

Size of Shape: 2000

Amplitude: 6 100

OK Apply Cancel

**XWIN-NMR Shape Tool**

File Split Shapes Analyze Manipulate Options **7**

\*2 \*8 [up arrow] Selecti  
/2 /8 [checkbox] Cal

Phase Modulation acc. to Offset Freq. [offs] (circled in red)  
Single Sine Modulation [sinm2]  
Single Cosine Modulation [cosm2]  
Modulation acc. to Freq. Sweep  
Power of Amplitude [power]  
Scale Amplitude [scale]  
Add constant Phase [addphase]  
Time Reversal [trev]  
Calc. Shape from Excitation Region [region]  
Add Shapes

Rectangle

Amp

50

0

Save Ctrl+S

Save As... **9**

Save Gradient

**Editor for Command Parameters**

Parameters for Manipulate Command offs: **8**

- Beginning at Phase 0 (ly -> lz)
- Phase = 0 at Middle of Shape (ly -> -ly, lz -> -lz)
- Ending at Phase 0 (lz -> -ly)

Reference Frequency:

- No Reference Frequency specified
- Reference = O1 from current Data Set
- Reference = First Frequency in List

Frequencies taken from Frequency List

With additional Phase Setting

With additional Scaling

Length of Pulse (in us): 2000

Name of Frequency List: freqlist

OK Apply Cancel

# Multiple Solvent Suppression: WET




## Water suppression Enhanced through T1 effects

- transferring solvent magnetization to the transverse plane and
  - dephasing this solvent magnetization prior to the excitation of the sample magnetization
- ⇒ 4 *pulse-gradient* units to reduce recovered solvent magnetization (recovery based on T1 relaxation during dephasing gradients).
- Multiple solvent suppression by using off-resonance phase shifted pulses (shifted laminar pulses, SLP)

## Gradients and selective pulses: **BEFORE** excitation

⇒ no evolution of  $^3J(\text{HH})$  scalar couplings, no phase errors.



# Multiple Solvent Suppression: WET



WET is commonly used in LC-NMR, often with  $^{13}\text{C}$  decoupling during WET and during acquisition.

pulse program	<b>lc1wet,</b> <b>lc1wetdc</b> or <b>lc1wetdw</b> <b>wt</b>	(1D-NMR with WET, dc = decoupling)  (nD-experiments using WET)
WET shape pulse length	<b>p11</b>	(= 15 - 25 ms)
power level for WET pulses	<b>sp7</b> <b>sp8</b> <b>sp9</b> <b>sp10</b>	$\text{sp1} + 0.87 \text{ dB}$ $\text{sp1} - 1.04 \text{ dB}$ $\text{sp1} + 2.27 \text{ dB}$ $\text{sp1} - 5.05 \text{ dB}$
shape for WET pulse	<b>spnam7-10</b>	gauss or sinc pulse, manipulate with solvent freqlist in <b>stdisp</b>
reference power level for WET	<b>sp1</b>	$90^\circ$ gauss or sinc pulse, calculate power

Numbers originate from numeric solutions of Bloch equations with average values for T1

# Multiple Solvent Suppression: WET



WET can be set up automatically, start with **xaua**

Parameter set

**LC1DWTDC**

(1D-NMR with WET using **lc1wetdc**)

Number of solvents to suppress

**130**

default value in parameter set: 2

parameters for decoupling

**pcpd2**

decoupling pulse length 80 - 100us

**pl12**

power level for GARP decoupling

**cpdprg2**

GARP

if no decoupling is desired: **pl12** = 120

parameters for WET pulses **sp1** **spnam1**

Standard soft pulses for 1H on channel F1 routed to amplifier A2:

Description:	Pulses:	P.Level:	Alignm.:	Name:	
90/270 excitation	PSH1	0	120	calc. 0.5	Gaus1.1000
180 refocussing	PSH2	0	120	calc. 0.5	Gaus1.1000
psh3	PSH3	14	54.25	calc. 0.5	Gaus1.1000
psh4	PSH4	0	120	calc. 0.5	Gaus1.1000

set parameters in **edprosol** for PSH3[F1]

in *standard soft pulses*: use ca. 20ms

90° gauss or sinc pulse, calculate power

in **stdisp**, use analyze → integrate

(older instruments: chose a pulse length

in a way that **sp1** < 55dB)

# Multiple Solvent Suppression: WET



Gradients and selective pulses: *BEFORE* excitation

- ⇒ no evolution of  $^3J(\text{HH})$  scalar couplings
- ⇒ no phase errors.

