



***New Features for Pulse Programs
and Xwin-NMR 3.5***

Detlef Moskau

March 2005

New Features for Pulse Programs

Frequency Lists in Pulse Programs

Define frequencies via constants



- XWIN-NMR 3.5 allows to set frequencies within a pulse program without the need of external frequency lists.
- Example:

.....

```
3 d11 fq=cnst23(bf ppm):f2
```

....


fq: take BF2 (as :f2 is used),
 add a frequency $\text{cnst23} \cdot \text{BF2} \cdot 10\text{e-6}$

$\text{SFO2} = \text{BF2} + \text{cnst23} \cdot \text{BF2} \cdot 10\text{e-6}$

cnst23: could be C-alpha chemical shift, e.g.

*Calculate Excitation Offset for
Shaped Pulses within a Pulse
Programs*

spoff calculation



- Excitation offsets for shaped pulses can be defined directly in a pulse program:

example 1:

```
;hbhaconhgp3d
```

```
;avance version (02/05/31)
```

```
;HBHACONH
```

```
"spoffs2=0"
```


```
"spoffs3=0"
```

```
"spoffs5=bf2*((cnst21-cnst23)/1000000)"
```

```
;cnst21: CO chemical shift (offset, in ppm)
```

```
;cnst23: Caliphatic chemical shift (offset, in ppm)
```

spoff calculation



- Excitation offsets for shaped pulses can be defined directly in a pulse program:

example 2:


```
;hbhanhgp3d
;avance version (02/05/31)
;HBHANH

"spoffs2=0"
"spoffs3=0"
"spoffs5=bf2*(cnst21/1000000)-o2"

;cnst21: CO chemical shift (offset, in ppm)
```

Alignment of Pulses

Pulse Alignment



- XWIN-NMR 3.5 allows alignment of pulses in a very easy way
- Example:

(**center** (p2 ph1):f1 (p14:sp3 ph1):f2)

center pulses

(**ralign** (p1 ph2) (p3 ph9):f2)

**right alignment of
shorter pulse**


(**lalign** (p1 ph2) (p3 ph9):f2)

**left alignment of
shorter pulse**

(**reference** (p1 ph2) **center** (p3 ph9):f2 **ralign** (p21 ph2):f3)

**use p1 as reference
for alignment**

Pulse Alignment



- Simultaneous pulses do not need to be coded on the same line
- Example:

(**reference** (p1 ph2) **center** (p3 ph9):f2 **ralign** (p21 ph2):f3)


OR:

(
 reference (p1 ph2)
 center (p3 ph9):f2
 ralign (p21 ph2):f3
)


Pulse Alignment



```
(  
  reference (d0 p1 ph2 d0) :f1  
  center (p2 ph3) :f2  
  ralign (p3 ph3) :f3  
)
```




$p1 > p2$
 $p1 > p3$



$p1 < p2$
 $p1 > p3$

Pulse Alignment




- New syntax for delays
- Example:

"DELTA=d0*2+larger(p14,p22)-p14+36u"

*New syntax for
pulsed field gradients*

Pulse Programming: PFG's



XWINNMR3.1 options for programming gradients

GRADIENT(cnst21)

gradient of length *p16* '(GS-syntax')
gradient sequence and assignment of gradient
strength for individual pulses are defined in
GRDPROG

2u:ngrad

set next value according to GRDPROG

gron1, groff

gradient ON / OFF according to GPZ1


p16:gp1

gradient of length *p16*,
gradient strength = *gpz1* (*gpx1, gpy1*)
gradient shape defined by *gpnam1*

*p16:gp1*EA*


multiplier for gradient strength; allow to invert
gradient
sign (useful for echo-antiecho gradient selection)
only XWIN-NMR version ≥ 2.5

Pulse Programming: PFG's



- Syntax *GRADIENT(cnst21)* not used anymore
- Syntax *ngrad* not used anymore
- Pulse programs **gs** substituted by **gp**
- Content of *GRDPROG* is now coded within pulse program

Pulse Programming: PFG's



```
;imggp3d
```

```
.....
```

```
lgrad r2d<2d> = td1
```

```
lgrad r3d<3d> = td2
```

```
.....
```

```
5u          grad{ (0)+r3d(5.057) | (0)+r2d(5.057) | (3.112) }
```

```
d27
```

```
5u          grad{ (0)          | (0)          | (0)          }
```

```
d15
```

```
5u          grad{ (0)          | (0)          | (-6.739) }
```

```
d21
```

```
ACQ_START(ph30,ph31)
```

```
aq DWELL_GEN:f1
```

```
5u          grad{ (0)          | (0)          | (0)          }
```

```
rcyc=2
```

```
5m st lgrad r2d
```


```
lo to 2 times td1
```

```
5m zgrad r2d lgrad r3d
```

```
lo to 3 times td2
```



Setting precompiler options

Setting precompiler options



- Precompiler options allow parts of the pulse program to be used or to be neglected, depending on precompiler options
- This allows to have one pulse program, which can be used for solvent suppression either by presaturation or by WATERGATE, e.g.
- Optioned can be
 - defined in the pulse program
 - disadvantage: pulse program has to be modified when the option is changed
 - Acquisition parameter ZGOPTNS
 - several option flags can be set together

Setting precompiler options



Example: setting precompiler options within the pulse program:

presaturation or WATERGATE?

```
define PRESAT
;undef PRESAT
;define WATERGATE
undef WATERGATE
```

```
#ifdef PRESAT
```

```
  d12 pl9:f1
```

```
  d1 cw:f1
```

```
  d13 do:f1
```

```
  d12 pl1:f1
```

```
#endif PRESAT
```


```
#ifdef WATERGATE
```

```
  d1 pl1:f1
```

```
#endif WATERGATE
```

.....

Setting precompiler options



Example: setting precompiler options with ZGOPTNS:

presaturation or WATERGATE?

A) ZGOPTNS	-DPRESAT	used presaturation
B) ZGOPTNS	-DWATERGATE	used WATERGATE
C) ZGOPTNS	-DPRESAT -DWATERGATE	used presaturation and WATERGATE

```
#ifdef PRESAT
```

```
  d12 pl9:f1
```

```
  d1 cw:f1
```

```
  d13 do:f1
```

```
  d12 pl1:f1
```

```
#endif PRESAT
```

```
#ifdef WATERGATE
```

```
  d1 pl1:f1
```


```
#endif WATERGATE
```

.....

*The **mc** macro statement*

(available since XWIN-NMR 3.1)

Pulse Programming: example HNCO



p16:gp4*EA

echo-antiecho gradient

go=2 ph31 cpd3:f3

d11 do:f3 wr #0 if #0 zd

6m ip6*2 igrad EA


increment phase, invert echo-antiech gradient

if "d30 < d21-p16-d16-4u" goto 98 conditional loop with runtime check


d11 id10
d11 dd30
d11
goto 99

98 d11 id10
d11 dd30
d11 id28

99 d12
lo to 4 times l13
d11 rd10 ip4
d11 rd30
d11 rd28
lo to 5 times 2
d11 id0
lo to 6 times l3



Pulse Programming: the mc macro statement




Advantages of the mc statement:

simplified pulse programming:

One pulse program can be used for different 2D phase modes (TPPI, States or States-TPPI)

Pulse Programming: the mc macro statement




mc clause	possible mode	values of F1-FnMODE
F1QF	phase insensitive	QF
F1PH	phase sensitive	QSEQ, States, TPPI, States-TPPI
F1EA	Echo-Antiecho	Echo-Antiecho

The acquisition and processing parameter FnMODE define the mode for F1 (and F2, for 3D-experiments) quadrature

For details, see XWINNMR help -> Other topics -> Writing pulse program

Pulse Programming: the mc macro statement



up to XWINNMR 2.6

new: XWINNMR 3.0

```
;noesytp
```

```
#include <Avance.incl>
```

```
"d0=3u"
```

```
1 ze
```

```
2 d1
```

```
3 p1 ph1
```

```
d0
```

```
p1 ph2
```

```
d8
```

```
p1 ph3
```

```
go=2 ph31
```

```
d1 wr #0 if #0 ip1 id0 zd
```

```
lo to 3 times td1
```

```
exit
```

```
;noesyph
```

```
#include <Avance.incl>
```

```
"d0=3u"
```

```
1 ze
```

```
2 d1
```

```
3 p1 ph1
```

```
d0
```

```
p1 ph2
```

```
d8
```


```
p1 ph3
```

```
go=2 ph31
```

```
d1 mc #0 to 2 F1PH(ip1, id0)
```

```
exit
```


Pulse Programming: the mc macro statement



Examples:

MC clause	t1 quadrature mode	action after	
		odd increment	even increment
F1PH(ip1, id0) F1PH(ip1, id0)	TPPI States-TPPI	ip1+ id0 ip1	ip1 + id0 again id0
F1PH(rd10 & rd30 & ip4, id0) F1PH(rd10 & rd30 & ip4, id0)	TPPI States-TPPI	rd10+rd30+ip4+id0 rd10+rd30+ip4	id0

Pulse Programming: the mc macro statement



Details

within XWINNMR, nD-pulse sequences will be expanded:

example for NOESY, F1-FnMODE = States-TPPI

```
define loopcounter ST1CNT
```

```
"ST1CNT = td1 / ( 2 ) "
```

```
"MCWRK = 0.500000 * d1"
```

```
"MCRES = d1- d1"
```

```
1 ze
```

```
2
```

```
LBLSTS1,
```

```
LBLF1,
```

```
3 p1 ph1
```

```
d0
```

```
p1 ph2
```

```
d8
```

```
p1 ph3
```

```
go=2 ph31
```

```
MCWRK wr #0 if #0 zd ip1
```

```
lo to LBLSTS1 times 2
```

```
MCWRK id0
```


```
lo to LBLF1 times ST1CNT
```

loop counter

delays to ensure steady-state

jump address

Pulse Programming: the mc macro statement




Example 1: **zg**

```
1 ze  
2 30m  
  d1  
  p1 ph1  
  go=2 ph31  
  30m mc #0 to 2 F0(zd)  
exit
```

New jump address added here

Jump address 2: just a delay

Pulse Programming: the mc macro statement



Example 2: **zgdc**


```
1 ze
  d11 pl12:f2
2 30m do:f2
  d11 cpd2:f2
  d1
  p1 ph1
  go=2 ph31
  30m do:f2 mc #0 to 2 F0(zd)
exit
```

 New jump address added here

Jump address 2: switch **OFF** the decoupler!!!

go=2 loop goes to jump address 2

Pulse Programming: the mc macro statement




Example 3: **zgpr**

```
1 ze
2 30m
  d12 pl9:f1
  d1 cw:f1 ph29
  4u do:f1
  d12 pl1:f1
  p1 ph1
  go=2 ph31
  30m mc #0 to 2 F0(zd)
exit
```

New jump address added here

Jump address 2: just a delay
power level setting done directly after

Pulse Programming: the mc macro statement




Example 4: **hsqcetgp**

```
1 ze
  d11 pl12:f2
2 d1 do:f2
3 (p1 ph1)
  d4 pl2:f2
.....
go=2 ph31 cpd2:f2
  d1 do:f2 mc #0 to 2
    F1EA(igrad EA, id0 & ip3*2 & ip6*2 & ip31*2)
exit
```

A red arrow points from the text "New jump address added here" to the line "2 d1 do:f2".

Jump address 2: switch OFF the decoupler!!!

Pulse Programming: the mc macro statement



Example 5: **hncogp3d**

```
1 d11 ze
  d11 pl16:f3
2 d11 do:f3
3 d1 pl1:f1
  p1 ph1
  d26 pl3:f3
.....
go=2 ph31 cpd3:f3
d11 do:f3 mc #0 to 2
  F1PH(rd10 & rd29 & rd30 & ip4, id0)
  F2EA(igrad EA & ip6*2, id10 & id29 & dd30)
exit
```

A red arrow points from the text "New jump address added here" to the line "2 d11 do:f3".

Jump address 2: switch **OFF** the decoupler!!!