

A large, semi-circular graphic on the left side of the slide. It has a blue-to-purple gradient and contains a white, stylized waveform that resembles a series of pulses or a modulated signal.

# ***Introduction to shaped pulses***

## ***Definition:***

***a 'shaped pulse' is an amplitude modulated RF-pulse***

- *they are pulses with a customized bandwidth and behavior*

## *Bandwidth considerations:*

- *selective excitation/inversion/refocussing of a **narrow** bandwidth  
(e.g. 270° Gaussian pulse in 1D-COSY, 1D-TOCSY)*
- *selective excitation / inversion / refocussing of a **larger** bandwidth  
(e.g. G4 pulses for <sup>13</sup>C bandselective excitation in HSQC)*
- ***wideband** excitation at lowest possible power level  
(e.g. CHIRP and WURST decoupling pulses)*

# Examples of shaped pulses



## Classical

Square  
Gauss  
Sinc  
Burp

e-Burp1  
e-Burp2  
i-Burp1  
i-Burp2  
re-Burp  
u-Burp

## Gaussian Cascades

G4  
G3  
Q5  
Q3

## Snob

e-Snob  
i-Snob2

## r-Snob

## Vega

EVega1  
IVega

## Adiabatic

Hyperbolic Secant  
Sine/Cosine

Smoothed Chirp  
Composite smoothed Chirp  
Wurst

## Decoupling

Swirl

i=inversion 180°


e=excitation 90°

r=refocusing 180°

u=universal 90°, 180°




# Common pulse shapes and their FT



## Ideal selective pulse:


- **Selects a rectangular narrow region**
- **Uniform excitation (flat top)**
- **Negligible sidebands**
- **Uniform phase response (self refocussing)**
- **Short time duration (relaxation, J-coupling)**
- **Suitable for  $M_z$ ,  $M_{xy}$ , 90 and 180 pulses**

**Some are highly specialized pulses (single application)**




# Selective 90 degree pulses


### Square, 10 ms




### Sinc, 10ms




### Gauss, 10ms




### G4, 10 ms



### u-BURP, 10 ms




### e-BURP, 10 ms




# Selective 180 degree pulses


### Square, 10 ms




### G3, 10 ms



### i-Burp, 10ms



### re-Burp, 10ms





# *How to choose a shaped pulse*

## 1) *Selectivity*

The region of interest should be irradiated as selectively as possible

## 2) *Length*

Relaxation and J-evolution might take place during shaped pulses

## 3) *Power*

Low integral power or low peak power

## 4) *Phase*


Pulse should give pure phase within the region of interest

# Selecting a shape

- according to its 'top hat' behaviour (simulated profile):

<u>Excitation</u>	e-Burp2 G4 EVega1	↑	higher priority
<u>Refocussing</u>	re-Burp Q3 Hermite r-Snob Sinc (3) Gauss Hyperbolic Secant composite smoothed CHIRP (wideband refocussing)	↑	higher priority
<u>Inversion</u>	G3 Hyperbolic Secant i-Burp2 i-Snob3 Hermite Sinc (3) CHIRP (wideband inversion)	↑	higher priority

# The bandwidth parameter $\Delta \omega_e / \Delta \omega_t$



$\Delta \omega_e / \Delta \omega_t = \text{excitation region} / \text{transition region}$

width of excitation

$\Delta \omega_e$

width of transition region

$\Delta \omega_t$

bandwidth

$\Delta \omega$

pulse length

$\Delta T$

peak amplitude

$\gamma B_{1\text{max}}$

Values for refocussing:


Gauss	0.16
Square	0.38
Sinc (3)	0.41
Sinc (7)	0.92
Hermite	0.97
r-Snob	0.99
G3	1.69
Q3	2.10
re-Burp	3.33
IVega	6.88



higher priority

larger value for  $\Delta \omega_e / \Delta \omega_t$  means better selectivity

# The bandwidth factor $\Delta \omega * \Delta T$



*How long the pulse must be to give a certain excitation bandwidth?*

width of excitation

$\Delta \omega_e$

width of transition region

$\Delta \omega_t$

bandwidth

$\Delta \omega$

pulse length

$\Delta T$

peak amplitude

$\gamma B_{1,max}$


# The bandwidth factor $\Delta \omega^* \Delta T$

Values for refocussing:

Square		0.75	
Gauss		0.88	
Hermite	2.24		
r-Snob		2.33	
G3		3.42	↑ higher priority
Q3		3.45	
Sinc (3)		4.52	
re-Burp		5.81	
Sinc (7)		8.66	
IVega		9.25	

*smaller value for  $\Delta \omega^* \Delta T$  means shorter pulse length required for a certain bandwidth*

*maximum pulse power =  $\gamma B_1 \max$  (peak amplitude)*



Shaped pulse

$$I = (A_i) * \gamma B_1 \max / 2\pi * \tau_p$$

Square pulse

$$I = \gamma B_1 \max / 2\pi * \tau_p$$




Values for refocussing 25Hz region (high selectivity,  $\Delta\omega$  is constant):

Shape	$\Delta\omega_e / \Delta\omega_t$	pulse length [ms]	$\gamma B_1 \max / 2\pi$ [Hz]	
Sinc (3)	0.41	180.8	15.6	
IVega	6.88	369.6	16.6	
Square	0.38	30.0	16.7	
Sinc (7)	0.92	346.0	19.7	↑ Less power
Q3	2.10	138.0	23.9	
r-Snob	0.99	93.2	25.1	
G3	1.69	136.8	26.2	
re-Burp	3.33	232.4	26.9	
Hermite	0.97	89.6	30.4	
Gauss	0.16	35.2	34.5	

# Maximum pulse power $\gamma B_1 \max$

Values for refocussing 25Hz region (high selectivity,  $\Delta\omega$  is constant):

Shape	$\Delta\omega_e / \Delta\omega_t$	pulse length [ms]	$\gamma B_1 \max / 2\pi$ [Hz]
Sinc (3)	0.41	180.8	15.6
IVega	6.88	369.6	16.6
Square	0.38	30.0	16.7
Sinc (7)	0.92	346.0	19.7
Q3	2.10	138.0	23.9
r-Snob	0.99	93.2	25.1
G3	1.69	136.8	26.2
re-Burp	3.33	232.4	26.9
Hermite	0.97	89.6	30.4
Gauss	0.16	35.2	34.5

 Less power

# Comparison of the best pulses



*for band-selective excitation*

Shape	$\Delta \omega_e / \Delta \omega_t$	$\Delta \omega * \Delta T$	$\gamma B_1 \text{max} / 2\pi$ [Hz]
	<i>selectivity</i>	<i>length</i>	<i>power</i>
re-Burp	3.33	5.81	26.9
Q3	2.10	3.45	23.9
Hermite	0.97	2.24	30.4
r-Snob	0.99	2.33	25.1
Sinc (3)	0.41	4.52	15.6

*Choose either re-Burp or Q3 as a good refocussing pulse:*

**re-Burp:**

***advantage: more selective***

***disadvantage: longer pulse length, more power***

**Q3**

***advantage: shorter pulse, less peak power***

***disadvantage: less selective***





# ***3. Introduction to adiabatic pulses***

## 1) *In insensitive to $B_1$ inhomogeneity*

- *magnetization can be 'picked up' , excited or inverted in coil regions of reduced  $B_1$  homogeneity*
- *increased signal-to-noise can be expected*

## 2) *Wide inversion bandwidth*

- *HSQC / DEPT experiments: solving the problem of  $^{13}\text{C}$  pulse offset effects*
- *CHIRP95 and WURST decoupling: less power or wider decoupling range*


## 3) *In insensitive to power missettings*

- *threshold of minimum and maximum RF field where pulses are adiabatic*


*Note: Adiabatic 180-degree pulses discussed only*

# The adiabatic rotation

*Schematic representation of movement of magnetization along the effective  $B_1$  field during an adiabatic pulse*



a) magnetization  $M$  during an adiabatic pulse, no apodization on amplitude



no apodization on amplitude  
with apodization on amplitude

b) trajectory of  $M$  on a sphere during an adiabatic pulse

# The adiabatic condition


The adiabatic condition:

$$\omega / \gamma B_0 \ll 1$$

The adiabatic pulse:

- *'rapid passage' of frequency*
- *quasi parabolic phase change*
- *magnetization 'spin-locked' during pulse*

amplitude      frequency      phase




**BIR-4 pulse:**

*'apodization' of amplitude*

*off-resonance: fast frequency change*

*on-resonance: slow frequency change*

amplitude      frequency      phase



**CHIRP pulse:**

*optional 'apodization' of amplitude*

*constant frequency sweep*

*quasi-parabolic phase change*

**NOTE: SLP-pulses have  
constant phase change**



Shape	$\Delta \omega_e / \Delta \omega_t$ <b>selectivity</b>	$\Delta \omega * \Delta T$ <b>length</b>	$\gamma B_1 \text{max} * \Delta T$ <b>power*length</b>	
sm. CHIRP (30%)	2.81	29.16	6.18	
sm. CHIRP (10%)	1.49	29.52	6.18	↑ <b>higher priority</b>
wurst (n=20)	2.16	28.26	6.18	
power h. secant (n=8)	2.33	34.10	6.94	
sincos	4.57	35.12	7.71	
power wurst (n=2)	5.48	36.70	10.1	

*∴ All these pulses are quite similar*

## *Which shape?*

**CHIRP and WURST are the most common adiabatic pulses**

## *What are the typical parameters?*

for inversion: **smoothed CHIRP**

**0.5ms, 20% smoothing, 60-80kHz sweep**

for refocussing: **composite smoothed CHIRP**

**2ms, 20% smoothing, 60-80kHz sweep**

CHIRP-decoupling: **smoothed CHIRP**

**1.5ms, 20% smoothing, 18-40kHz sweep**

**use decoupling program `p5m4sp180`**

*Should all 180° pulses in HSQC be substituted by adiabatic pulses?*

organic molecules: **CHIRP's** for refocussing and inversion

<sup>13</sup>C labeled proteins: **CHIRP's** for inversion only

*What is the difference between smoothed CHIRP and composite smoothed CHIRP?*

smoothed CHIRP: use for inversion only

comp. Sm. CHIRP: can be used for inversion and refocussing

## *How to calibrate an adiabatic pulse?*

ShapeTool -> analyse -> „integrate adiabatic shape“

## *Power level for $^{13}\text{C}$ 0.5ms smoothed CHIRP and*

*$^{13}\text{C}$  2ms composite smoothed CHIRP ?*

corresponds to 25  $\mu\text{sec}$   $^{13}\text{C}$  pulse

adiabaticity Q-factor : 5 full adiabaticity

## *Power level for 1.5ms smoothed CHIRP 18-40kHz in CHIRP95 decoupling?*

use GARP power level + 2-6 dB

adiabaticity Q-factor : 2-3 is sufficient (additional averaging due to the supercycles in the decoupling program)

**1) Adiabatic pulses are included in XWIN-NMR, e.g.**

- *Crp60,0.5.20.1* for inversion
- *Crp60comp.4* for refocussing

**2) Many experiments are written with adiabatic pulses**

**3) Power of adiabatic pulses directly from ShapeTool**

*„integrate adiabatic pulse“ - option*

**Crp60,0.5.20.1** *size of shape* 1000  
*total sweep-width* 60000  
*length of pulse* 500  
*% to be smoothed* 20  
*low to high field* -1

*power = same as 25us hard pulse*




**ShapeTool -> analyse-> integrate adiabatic shape**  
**select Q-factor=5**

# Simulation of inversion behaviour



**ShapeTool -> Analyse -> Simulate**

**power level=same as 25us hard pulse**  
**power:  $1/(4*0.000025s)=10000\text{Hz}$**



Check the experiment parameters

SPNAME0

Mx (0)  initial magnetization

My (0)  initial magnetization

Mz (0)  initial magnetization

P 0  us fixed pulse

SP 0  Hz power for shaped pulses

N  number of calculated points

Start  Hz first RF offset

Step  Hz offset increment

OK Cancel

**low adiabaticity**  
**(Q-factor~2-3)**

**high adiabaticity**  
**(Q-factor=5)**




# Adiabatic refocussing

**Crp60comp.4**

*size of shape*                    1000  
*total sweep-width*            60000  
*length of basic unit*          500  
*% to be smoothed*            20  
*low to high field*              -1

*length of pulse*                2000us  
*power level= 25us hard pulse*



**ShapeTool -> analyse-> integrate adiabatic shape**  
**select Q-factor = 5**




# Create a pulse for adiabatic decoupling



**Crp40,1.5.20.1** size of shape 1000  
total sweep-width 40000  
length of pulse 1500  
% to be smoothed 20  
low to high field -1

**cpd-program p5m4sp180**  
sweep-width (decouples ~80% of 40 000Hz)  
length of pulse 1500us  
power level corresponding to ~80us hard pulse  
=> ca. 3000Hz



**Results**

integradia:  
Sweep Rate on Resonance (in Hz/sec): 2.66667e+07  
gammaB1(max)/2pi / sqrt(Q) (in Hz): 2060.13  
Corresponding 90 degree square pulse (in usec) 76.7495  
Change of power level compared to level of hard pulse in dB: 14.1797


Calculator:  
Q:   
GammaB1(max)/2pi (in Hz):

Update Parameters      Seen

**ShapeTool: analyse-> integrate adiabatic shape**  
**select Q-factor= 2-3 sufficient (supercycles!)**




# Adiabatic refocusing pulse in DEPT



yellow =  
standard  
DEPT45

black =  
DEPT45 with  
adiabatic  
inversion



*Refocusing pulse*  
*Crp60comp.4*

## Adiabatic refocussing pulse in DEPT experiments:

pulprog  
spnam2  
sp2  
p12

deptsp45 / deptsp90 / deptsp135


Crp60comp.4

power level corresponding to 25us  $90^\circ$ - $^{13}\text{C}$ -pulse

2000us

# Adiabatic pulses in $^{13}\text{C}$ -HSQC

**Advantage: higher sensitivity over a larger excitation bandwidth**



green =  
standard HSQC

red = HSQC  
with adiabatic  
inversion


# Adiabatic pulses in $^{13}\text{C}$ -HSQC

**Two options:**

1) *hsqcetgpsp* with adiabatic inversion  
(proteins,  $^{13}\text{C}$  labeled)

2) *hsqcetgpsp.2* with adiabatic inversion and refocussing  
( $^{13}\text{C}$  at natural abundance)

HSQC with gradients



## $^{13}\text{C}$ labeled - only adiabatic inversion:

<b>pulprog</b>	<b>hsqcetgpsisp</b>
<b>spnam3</b>	<b>Crp60,0.5.20.1</b>
<b>sp3</b>	<b>power level corresponding to 25us <math>90^\circ</math>-<math>^{13}\text{C}</math>-pulse</b>
<b>p14</b>	<b>500us</b>

## non labeled - adiabatic inversion und refocussing:

<b>pulprog</b>	<b>hsqcetgpsisp.2</b>
<b>spnam3</b>	<b>Crp60,0.5.20.1</b>
<b>sp3</b>	<b>power level corresponding to 25us <math>90^\circ</math>-<math>^{13}\text{C}</math>-pulse</b>
<b>p14</b>	<b>500us</b>

<b>spnam7</b>	<b>Crp60comp.4</b>
<b>sp7</b>	<b>power level corresponding to 25us <math>90^\circ</math>-<math>^{13}\text{C}</math>-pulse</b>
<b>p24</b>	<b>2000us</b>

**/NMRhome/exp/stan/nmr/lists/cpd/p5m4sp180**

**Note: this cpd-program uses sp15**

<b><i>cpdprg2</i></b>	<b><i>p5m4180</i></b>
<b><i>pcpd2</i></b>	<b><i>1500us</i></b>
<b><i>spnam15</i></b>	<b><i>Crp40,0.5.20.1</i></b>
<b><i>sp15</i></b>	<b><i>power (with Q-factor 2-3)</i></b>
<b><i>pl12</i></b>	<b><i>= sp15</i></b>

```
#setphase  
1 pcpd:sp15:0  
pcpd:sp15:150  
pcpd:sp15:60  
pcpd:sp15:150  
pcpd:sp15:0  
pcpd:sp15:0  
pcpd:sp15:150  
pcpd:sp15:60  
pcpd:sp15:150  
pcpd:sp15:0  
pcpd:sp15:180  
pcpd:sp15:330  
pcpd:sp15:240  
pcpd:sp15:330  
pcpd:sp15:180  
pcpd:sp15:180  
pcpd:sp15:330  
pcpd:sp15:240  
pcpd:sp15:330  
pcpd:sp15:180  
jump to 1
```

# Adiabatic bilevel decoupling programs



**/NMRhome/exp/stan/nmr/lists/cpd/  
bi\_p5m4sp\_4pl**

- this cpd-program uses sp15

<b>cpdprg2</b>	<b>bi_p5m4sp_4pl</b>
<b>pcpd2</b>	<b>1500us</b>
<b>spnam15</b>	<b>Crp40,0.5.20.1</b>
<b>sp15</b>	<b>power at Q-factor 2-3</b>
<b>pl30</b>	<b>= sp15</b>

**pl31 = sp15 minus 6dB**

**Note - use synchronous mode, cpd2s:f2**

**#setphase**

**bilev "l31=(nsdone+ds)%4+1"**

**1 pcpd\*0.5:0 pl=pl31**

**lo to 1 times l31**

**2 pcpd:sp15:0 pl=pl30**

**pcpd:sp15:150**

**pcpd:sp15:60**

**pcpd:sp15:150**

**pcpd:sp15:0**

**pcpd:sp15:0**

**pcpd:sp15:150**

**pcpd:sp15:60**

**pcpd:sp15:150**

**pcpd:sp15:0**

**pcpd:sp15:180**

**pcpd:sp15:330**

**pcpd:sp15:240**

**pcpd:sp15:330**

**pcpd:sp15:180**

**pcpd:sp15:180**

**pcpd:sp15:330**

**pcpd:sp15:240**

**pcpd:sp15:330**

**pcpd:sp15:180**

