

4D-NMR Experiments with Avance Spectrometer

Topspin 1.3

General procedure to setup a 4D-experiment



GETTING STARTED:

simply read standard parameter set with 'rpar HNCOCAGP3D' and select the corresponding 4D pulse program, then load pulses with 'getprosol'

DEFINING THE SPECTRAL WINDOWS FOR THE INDIRECT DIMENSIONS:

1. use CNST20 to define spectral width for CO-dimension
2. set SW{F1} and SW{F2} in the eda-window to define spectral window, e.g. for N and CA

DEFINING THE NUMBER OF INCREMENTS FOR THE INDIVIDUAL DIMENSIONS:

3. Use TD{F1} and TD{F2} to define number of increments for the F1- and F2-dimensions.
Use 'parmode' = 3D.
4. Use L3 to define number of increments for F3-dimension.
Use 'parmode' = 2D.

SELECTING DIMENSIONS:

5. use the acquisition parameter ZGOPTN to define which dimensions should be measured:
examples
a) 'ZGOPTN = -DNH_EVOL' would give the NH plane
b) 'ZGOPTN = -DCA_EVOL -DNH_EVOL' will give a regular 3D HNCOCA spectrum
c) 'ZGOPTN = -DCA_EVOL -DNH_EVOL -DCO_EVOL' will give a 4D HNCOCA spectrum



Example of a 4D Pulse Program



1. Defining the 4th dimension:

Parameters like spectral window and number of time increments for the first three frequency dimensions can be directly set within the ,eda' window of XWIN-NMR.

Parameters for the 4th dimension cannot be set with the ,eda' window, and therefore will be coded in the pulse program

2. Writing a 4D pulse program:

Starting from a corresponding 3D-pulse program, just a few modifications have to be done:

- b) Automated calculation of the time increment for the evolution delay of the 4th dimension
- c) Adding the evolution for the 4th dimension
- d) Adding actors for the evolution in the 4th dimension, e.g. TPPI phase incrementation, incrementation of evolution delay

Example of a 4D Pulse Program



Automated calculation of the time increment for the evolution delay of the 4th dimension:

```
"in20=1/(4*bf2*(cnst20/1000000))"
```

Adding evolution period for the 4th dimension (here, the CO evolution):

```
d20 ;start t3 evolution  
(center (p14:sp5 ph1):f2 (p22 ph1):f3 ) ;refocuss J(CO,N) and  
J(CaCO)  
d20 ;end t3 evolution
```

Example of a 4D Pulse Program



Factors for quadrature and evolution;

/ Definitions for individual planes*/*

ifdef CA_EVOL

d11 do:f3 wr #0 if #0 ip4 zd

lo to 3 times 2

d11 id0

lo to 4 times I1

endif

ifdef NH_EVOL

*d11 do:f3 wr #0 if #0 igrad EA ip6*2 zd*

lo to 5 times 2

d11 id10

d11 id29

d11 dd30

lo to 6 times I2

endif

ifdef CO_EVOL

d11 do:f3 wr #0 if #0 id20 ip3 zd

lo to 7 times I3


endif

individual or combinations of dimensions can be selected by setting the acquisition parameter ZGOPTNS.

NH plane of 4D-HNCOCA experiment




hncocagp4d 15N plane



CA plane of 4D-HNCOCA experiment



hncocagp4d CA plane



CO plane of 4D-HNCOCA experiment



hncocagp4d CO plane

