



AU-Programs

Detlef Moskau

Bruker BioSpin

AU programs

Pulse program: performed at run time on the spectrometer
real time
edpul

AU-program: performed on data station (INDY, ASPECT Station, X32)
no run-time access to spectrometer
indirect action on pulse program execution by changing
acquisition parameters prior acquisition start
slow
edau

Macros: Macros simply contain a series of Topspin-commands
slow
edmac

Commands used for AU-programs:
pure C code
pure Topspin standard commands (via library)
mixture of C-code and Topspin commands



AU programs

General structure:

GETCURDATA
FT
QUIT

Program loops:

TIMES(10)
<command sequence>
END

Additional loop commands:

TIMES2(i)
TIMES3(i)

TIMES(n)
...
TIMES2(m)
...
TIMES3(p)
...
END
END

AU programs

Termination of an AU program:

QUIT
QUITMSG ("end of AU program")
STOP
STOPMSQ("AU program stopped")

Nested AU programs:

XAU("AU-program")



AU programs

Predefined variables:

```
int i1, i2, i3;  
float f1, f2, f3;  
double d1, d2, d3;  
char [text100];  
char name[15], disk[15], user[15];  
int expno, procno;
```

AU programs

Setting parameters:

STOREPAR("parameter name", parameter value) F2-dimension

example:

```
STOREPAR("SW",12.5)
STOREPAR("D 5",0.001)
STOREPAR("SI", 16*1024)
```

STOREPARS("parameter name", parameter value) F2-dimension status parameter

STOREPAR1("parameter name", parameter value) F1-dimension

STOREPAR1S("parameter name", parameter value) F1-dimension status parameter



AU programs

Reading parameters:

FETCHPAR("parameter name", &variable) F2-dimension

FETCHPAR1("parameter name", &variable) F1-dimension

example: FETCHPAR("P 1", &f1)
 f1+=10;
 STOREPAR("P 1", f1)
 ZG



AU programs

Setting and reading plot parameters:

STOREPLPAR()

FETCHPLPAR()

syntax for content within brackets corresponds to STOREPAR and FETCHPAR



AU programs

Reading and writing complete parameter sets:

WPAR("parameter set name","parameter class")

RPAR("parameter set name","parameter class")

example: RPAR("COSY","all")

 RPAR("COSY","proc")

 WPAR("COSY","proc plot outd")



AU programs


Simple and useful C commands:

entering a variable:

```
char *Gets();  
...  
d1=atoi(Gets("enter number of experiments","10"));  
f1=atof(Gets("enter evolution delay 1/2J(HC)", "0.0036"));
```

Display error window, select OK or Cancel:

```
AUERR=Proc_err(ERRORP_AK_CAN|ERRORP_BEEP_DEF,  
                 "text");  
if(AUERR==ERR_CANCEL) abort
```



www.bruker-biospin.com