


DOSY

Diffusion Ordered SpectroscopY (DOSY)

An Introduction



- Separation of signals from different components
- Separation is based on differences in diffusion rates of the compounds
- A reasonable difference in the diffusion rates is required for successful separation
- Relative diffusion rates are determined
- Careful calibration of gradient strength of the spectrometer is **NOT** required

In contrast to the DOSY experiment:

Diffusion Experiments

- Absolute diffusion rates are determined
- Careful calibration of gradient strength of the spectrometer **IS** required

DOSY and Diffusion




DOSY and diffusion experiments have in common:

- identical experiments
- same pulse programs can be used

DOSY and diffusion experiments differ in:

- processing
- diffusion experiments require careful calibration of the gradient constant



Translational Diffusion




$$J = dn/dt/A$$

J = flux

dn/dt = number molecules transported/sec

A = sampling area of the reference plane



Translational Diffusion



$$I(q) = I_0 e^{-Dq^2 \Delta'}$$

D : diffusion coefficient

Δ' : $\Delta - \delta/3$

Δ : diffusion time and δ is a correction factor for finite gradients

q : $\gamma g \delta$

γ : gyromagnetic ratio


g : is the amplitude of the applied gradient

δ : duration of the applied gradient


I(q) depends on:

size, shape, aggregation etc. of molecule
temperature
viscosity

Translational Diffusion



Simple Pulse Sequence and Parameters



Acquisition parameters:

G_z : gradient strength

Δ : ,big delta'

diffusion takes place during this delay


δ : ,little delta'

duration of the gradient pulse

diffusion also takes place during gradient pulses

As the name says, the duration for ,big delta' is longer than ,little delta'

Simple Pulse Sequence and Parameters




The DOSY / diffusion experiment:

- In principle, the parameters G_z , Δ or δ can be varied for a gradient experiment
- Varying Δ or δ will vary the length of the pulse sequence and T_2 relaxation effects will add to diffusion effects
- To avoid relaxation effects the gradient strength G_z is varied

Translational Diffusion



Decay of magnetization as function of gradient intensity



Setup DOSY / Diffusion Experiment



1. Start from a ^1H spectrum:

- Switch OFF sample rotation
- Create new data set, change to 2D-mode
- TD{F1}: **24 – 32** (gradient ramp with 24 – 32 different gradient values)
- PULPROG: **stebpgp1s** or **ledbpgp2s**
- TD: **16k - 32k**
- D20: **100ms** (,big DELTA': needs to be adjusted according to diffusion rates of compounds)
- D21: **5ms** (eddy current delay)
- GPNAM*: **SINE.100** (square gradient can also be used, e.g. **SMSQ10.100**)
- P19: **600us** (spoil gradient)
- P30: **2 – 5ms** (,little DELTA': needs to be adjusted according to diffusion rates of compounds)

Note: for CryoProbes, max p30=4ms, GPNAM* = SINE.SHAPE

Setup DOSY / Diffusion Experiment



2. Calculate gradient ramp and start experiment:

- Use AU program „dosy“ to calculate the gradient ramp:

```
xau dosy 2 95 32 I n
```

(first gradient: 2%, last gradient 95%, 32 gradient values using a linear ramp, do not start the acquisition)

- Adjust receiver gain: “rga”
- Start experiment: “zg”

3. Optimize values for little / big delta:

- Process data with „xf2“ and phase correction
- Check signal decay:

Signals should decay to 5-10% of the original intensity

Decay incomplete: increase D20 and/or P30

Decay too fast: decrease D20 and/or P30


- Repeat acquisition

Setup DOSY / Diffusion Experiment



The decay can be easily checked from the stacked view after 'xf2'


Stacked plot of a diffusion spectrum after F2 transformation. On the F2 axis some of the strong signals of a peptide mixture are visible. The F1 axis shows the decay of the signals with increased gradient strength due to the diffusion of the molecules in solution.




Ex. the pulse calibration sample



optimal decay
d20=100ms



2D-DOSY with bipolar gradients pulprog=stebpgp1s
p30=5ms, sine shape gradients
d20=20ms (big delta delay)
d1=4sec, aq=3 sec
processing: LWF=1 SpiSup=1
sample: pulse calibration sample




Note: for CryoProbes, max p30=4ms, GPNAM* = SINE.SHAPE


Ex. the pulse calibration sample



decay incomplete
 $d_{20}=50\text{ms}$



decay too fast
 $d_{20}=200\text{ms}$



1. Setup processing parameters:

- `xau setdiffparm` copy acquisition parameters to DOSY processing parameters (little delta, big DELTA, gamma,..)
- `edp` edit parameters for FT in F2-dimension and number of points in the diffusion (F1) dimension


Note: DOSY processing requires good S/N, therefore window functions like EM with high values for LB, e.g. LB=5 might be required.

- `eddosy` edit processing parameters

2. Processing:

- `xf2` F2 Fourier transformation
- Phase correction
- `abs2` F2 automatic base line correction
- `dosy2d, dosy3d` calculate the 2D or 3D DOSY spectrum

eddosy parameters



eddosy parameters



Spectrum	ProcPars	AcquPars	Title	PulseProg	Peaks	Integrals	Sample	Structure	Fid	
<div style="display: flex; justify-content: space-between; align-items: center;"> ↶ P G 🔍 📄 📊 📈 📉 📌 <input style="width: 100px;" type="text"/> </div>										
General	General									
First	FITMODE =	exponential								Type of fit
Second	ExpVar =	Gradient								Variable parameter
Third	Xlist =	difflist								Variable parameter values file name
Baseline	Nstart =	0								Start of input points
	Ndata =	16								Number of input points (TD)
	Maxiter =	100								Maximum number of iterations
	EPS =	1								Tolerance
	Nexp =	1								Number of components to fit
	Noise =	63637.00								Noise level (S_DEV)
	PC =	4								Noise sensitivity factor
	SpiSup =	1								Spike suppression factor
	F1mode =	Peaks								F1 output data mode
	lmode =	Integral								Fitted intensity meaning
	Scale =	Linear								Scaling
	LWF =	1								Line width factor
	DISPmin =	1e-010								Lower display limit
	DISPmax =	1e-008								Upper display limit
	Npars =	7								Number of parameters
	Nvar =	2								Number of parameters to fit
	Gamma [Hz/G] =	4257.64000								Gamma
	Grad [G/cm] =	0.00000								Diffusion gradient
	Gdist [ms] =	0.00000								Gradient distance, big delta
	Glen [ms] =	0.00000								Gradient length, little delta
	First component									
	I1vary =	Yes								Fit intensity?
	I1 =	1000000000								Intensity
	I1min =	-2147483647								Minimum intensity
	I1max =	2147483647								Maximum intensity
	D1vary =	Yes								Fit diffusion coefficient?
	D1 [m²/s] =	1e-009								Diffusion coefficient
	D1min [m²/s] =	0								Minimum diffusion coefficient
	D1max [m²/s] =	1e-008								Maximum diffusion coefficient

Parameters will be set with *xau dosy*

Parameters will be set with *xau setdiffparm*

eddosy parameters



Spectrum	ProcPars	AcquPars	Title	PulseProg	Peaks	Integrals	Sample	Structure	Fid
<div style="display: flex; justify-content: space-between; align-items: center;"> General General </div>									
First	FITMODE =	exponential	Type of fit						
Second	ExpVar =	Gradient	Variable parameter						
Third	Xlist =	difflist	Variable parameter values file name						
Baseline	Nstart =	0	Start of input points						
	Ndata =	16	Number of input points (TD)						
	Maxiter =	100	Maximum number of iterations						
	EPS =	1	Tolerance						
	Nexp =	1	Number of components to fit						
	Noise =	63637.00	Noise level (S_DEV)						
	PC =	4	Noise sensitivity factor						
	SpiSup =	1	Spike suppression factor						
	F1mode =	Peaks	F1 output data mode						
	lmode =	Integral	Fitted intensity meaning						
	Scale =	Linear	Scaling						
	LWF =	1	Line width factor						
	DISPmin =	1e-010	Lower display limit						
	DISPmax =	1e-008	Upper display limit						
	Npars =	7	Number of parameters						
	Nvar =	2	Number of parameters to fit						
	Gamma [Hz/G] =	4257.64000	Gamma						
	Grad [G/cm] =	0.00000	Diffusion gradient						
	Gdist [ms] =	0.00000	Gradient distance, big delta						
	Glen [ms] =	0.00000	Gradient length, little delta						
	First component								
	I1vary =	Yes	Fit intensity?						
	I1 =	1000000000	Intensity						
	I1min =	-2147483647	Minimum intensity						
	I1max =	2147483647	Maximum intensity						
	D1vary =	Yes	Fit diffusion coefficient?						
	D1 [m ² /s] =	1e-009	Diffusion coefficient						
	D1min [m ² /s] =	0	Minimum diffusion coefficient						
	D1max [m ² /s] =	1e-008	Maximum diffusion coefficient						

PC: same function as for 'pp'

Suppression of 'tails' on peaks

Line broadening in diffusion dimension

DOSY peaks will be calculated and displayed within those values only!!!

eddosy in ProcPars



eddosy setdiffparm xf2 dosy2d

The screenshot shows a software interface with a toolbar and a parameter settings panel. The toolbar contains icons for undo, a 'P' button, a 'G' button, a plot icon, a fit icon, a zoom icon, a data icon, and a multi-view icon. A red arrow points from 'eddosy' to the 'G' button. A blue arrow points from 'setdiffparm' to the plot icon. A black arrow points from 'xf2' to the fit icon. A green arrow points from 'dosy2d' to the data icon. The parameter settings panel has a 'General' tab selected, with a text box containing 'Copy parameters from experiment (setdiffparm)'. Below this are two rows: 'First' with 'FITMODE =' and a dropdown menu set to 'exponential' (labeled 'Type of fit'), and 'Second' with 'ExpVar =' and a dropdown menu set to 'Gradient' (labeled 'Variable parameter').

1. No peaks at all in the DOSY spectrum:

- Are the correct values for big / little delta (*D20 / P30*) entered in eddosy?
- Do the display limits *DISPmin / DISPmax* in eddosy correspond to the expected diffusion rates?
- Check peak picking parameter PC

2. Poor separation of peaks:

- Effect not related to processing parameters:

Difference in diffusion rates is too small

Poor S/N, record spectrum again with more scans

- Change eddosy parameters:

SpiSup

larger values will *reduce tails* on peaks, but separation in F2 dimension will be reduced

LWF


values <1 will result in *better resolution* in diffusion dimension, but truncation artifacts might appear

DOSY processing: example for LWF




LWF: broadening in diffusion dimension, reduces truncation artifacts

LWF 0.3, SpiSup 1



LWF 4, SpiSup 1




DOSY processing: example for SpiSup




SpiSup: broadening in chemical shift dimension, reduces tails

LWF 4, SpiSup 1



LWF 4, SpiSup 100



DOSY: data



Data can also be analysed with the Topspin T_1/T_2 tool

F1 axis shows diffusion constant:

- logarithmic: $\log D$
- linear: $D * 1e^9$
- The unit for the diffusion constant is: m^2/sec


DOSY: sources of artifacts




1. Low signal-to-noise ratio

–Poor description of signal decay and bad separation of signals in diffusion dimension

Signal decay: high S/N



Signal decay: low S/N




DOSY: sources of artifacts



2. Convection

self diffusion measurements of a peptide in DMSO at different temperatures and with different pulse sequences. The signal amplitude is plotted versus the square of the gradient strength.




DOSY: effects of convection



2. Convection

- Reduce convection by:
 - High flow rates for temperature gas
 - Capillary tube instead of 5 mm tube
- Use pulse sequences with bipolar gradient to reduce the influence of convection in the DOSY / diffusion experiment



Pulse sequences using bipolar gradients:

step1s (ste=stimulated echo)


ledgp2s (led=longitudinal eddy current delay)

Note: for CryoProbes, max p30=4ms, GPNAM* = SINE.SHAPE

Example: two water species




Intra- and inter cellular water



Example: mixture of four esters




P-Hydroxybenzoic acid ester
R= -Me, -Et, -Prop, -But



Example: mixture of four esters




P-Hydroxybenzoic acid ester
R= -Me, -Et, -Prop, -But



Example: 3D DOSY-TOCSY




Projection view of the diffusion separation of compounds



Example: 3D DOSY-TOCSY



3D view shows TOCSY correlation peaks



AU Programs

- dosy, setdiffparm

Pulse Programs (1D-versions not shown)

Stimulated Echo Sequences:

- stegp1s, ste**bp**gp1s **bp**: bipolar gradients
- ste**bp**gp1s**19** **19**: 3919 WATERGATE
- ste**bp**gp**in**1s **in**: INEPT, for X-nuclei

Double Stimulated Echo Sequences:

- dstegp3s

Stimulated Echo plus LED Sequences:

- ledgp2s
- led**bp**gp2s **bp**: bipolar gradients
- ledbpgp**ml**2s**2d** **ml**: diffusion filtered **2D** TOCSY
- ledbpgp**ml**2s**19****2d** **19**: additional WATERGATE
- ledbpgp**co**2s**3d** **co**: **3D** DOSY-COSY
- ledbpgp**ml**2s**3d** **ml**: **3D** DOSY-TOCSY
- ledbpgp**no**3s**3d** **no**: **3D** DOSY-NOESY

Files used for DOSY and diffusion



Gradient calibration constant is stored in:

XWINNMRHOME/conf/instr/gradient_calib

The gradient ramp (created with 'xau dosy') is stored in:

XWINNMRHOME/exp/stan/nmr/lists/gp/Difframp

For calculation of diffusion values the following file is used:

/DISK/data/USER/nmr/NAME/EXPNO/proc/PROCNO/difflist

A decorative graphic on the left side of the slide, consisting of a semi-circular shape with a blue-to-purple gradient. Inside this shape is a white sine wave that oscillates across the width of the semi-circle. To the right of the semi-circle, there are faint, light blue horizontal lines that resemble a scale or a grid.


Processing Diffusion Experiments

Diffusion experiments: processing




Processing of data from diffusion experiments:

- Diffusion experiments can be performed using exactly the same setup procedure and pulse programs like a 2D DOSY experiment
- Processing will be done using a guided processing tool
- First, process data with *xf2* and phase correction




Diffusion experiments: processing



Diffusion experiments: processing




Go through the individual points step-by-step:



Diffusion experiments: processing




Extract FID:



Diffusion experiments: processing




FT and phase correction:



Diffusion experiments: processing




Define peaks to be analysed using integral ranges:



Diffusion experiments: processing




Define peaks to be analysed using integral ranges:



Diffusion experiments: processing



Define peaks to be analysed using integral ranges:



Diffusion experiments: processing



Define parameters for fitting:

1.

Relaxation parameters

General Parameters

1 FID # for phase determination

1000.0 Left limit for baseline correction

-1000.0 Right limit for baseline correction

5 Number of drift points

1.0E-5 Convergence limit

16 Number of points

1 First slice

1 Slice increment

Fitting Function

vargrad Function Type

expdec Number of components

gaussdec List file name

lorgauss Increment (auto)

linear

varbigdel

varlitdel

vargrad

raddamp

Additional Parameters

10000.0 GAMMA(Hz/G)

10.0 LITDEL(msec)

100.0 BIGDEL(msec)

1.0 GRADIEN(G/cm)

OK Cancel

Fitting Function

Start Calculation

Display Report

Print/Export

2.

Relaxation parameters

General Parameters

1 FID # for phase determination

1000.0 Left limit for baseline correction

-1000.0 Right limit for baseline correction

5 Number of drift points

1.0E-5 Convergence limit

16 Number of points

1 First slice

1 Slice increment

Fitting Function

vargrad Function Type

1 Number of components

vdlist List file name

vdlist Increment (auto)

auto

vplist

vclist

difflist

dw

10000.0 GAMMA(Hz/G)

10.0 LITDEL(msec)

100.0 BIGDEL(msec)

1.0 GRADIEN(G/cm)

OK Cancel

Fitting Function

Start Calculation


Display Report

Print/Export

Diffusion experiments: processing



Start the calculation



A decorative graphic on the left side of the slide, consisting of a semi-circular shape with a blue-to-purple gradient. Inside the semi-circle is a white sine wave that oscillates across the width of the semi-circle. To the right of the semi-circle are several horizontal lines of varying lengths, resembling a scale or a set of data points.

Gradient Calibration for Diffusion Experiments

Gradient Calibration



1. Calibrate the temperature:

- Record a ^1H -spectrum of methanol or glycol:
 - Methanol: Low temperature 175 to 330 K (-98 to +57 °C)
 - Glycol: High temperature 310 to 410 K (+37 to +137 °C)
- Define plot region and use *xau calctemp* to calculate the temperature
- Methanol sample: 4% MeOH in Methanol-d4
- Glycol sample: 80% Glycol in DMSO-d6

2. Gradient calibration:

1. A diffusion experiment will be done for a sample with known diffusion value

Sample	Temperature [C]	D [$10^{-9}\text{m}^2/\text{s}$]
H ₂ O	20	2.031
H ₂ O	25	2.299
D ₂ O	25	1.872

Data by courtesy of Dr. M. Holz (Inst. of. Phys. Chem., Uni Karlsruhe, FRG)

Gradient Calibration



2. Gradient calibration:

2. From the difference $D_{measured}$ and $D_{literature}$ the gradient calibration constant **GCC** will be calculated:

$$GCC_{New} = GCC_{old} * \sqrt{\frac{D_{measured}}{D_{literature}}}$$

3. The new gradient calibration constant has to be entered using the TopSpin command **setpre**:

Parameter	Value	Default
Grad. calib. const. [Hz/cm]	9.998	0.100
Grad. calib. const. [G/mm]	0.00	1.00
Rate to measure temperature	0.9998	0.0100
Rate to check error status	0.90	1.00
Fast Base Preemp. Z	0.09998	0.00100
Gain [%]	0.00	1.00

Gradient Calibration



Note:

The gradient calibration has to be done BEFORE recording diffusion experiments, as the file *difflist* stored with the data set is calculated from the gradient calibration constant and the entries of *Difframp*.

Further information:

TopSpin -> Help menu -> Software and Application Manuals

1D and 2D Step-by-Step Advanced

Dosy

DOSY CONTIN/DECRA

Diffusion