

Pulse Programming

Elements of Pulse Programs

Trim Pulses

Flip back pulses

Shaped pulses in triple resonance experiments

Trim Pulses

Flip back pulses

Shaped pulses in triple resonance experiments

Purge pulses (spin-lock or trim pulses)

selectively dephase coherences orthogonal to the RF-field (I_x, I_z) while preserving coherences locked along the RF-field (I_y)

experiments: ^1H TOCSY, HSQC-type inverse, triple resonance experiments,...

power levels: **p110** for ^1H TOCSY, others **p11** (same as hard proton pulses!!, but not more than power level corresponding to 10 μs 90 degree pulse!!)


purge pulses: **p17** ($= < 2.5\text{ms}$); **p28** ($= < 1.0\text{ms}$, if p11 set as described above)

Note: make sure that the probe can stand the rather high power level used for the trim pulse. Otherwise the probe is damaged!!

Selective dephasing of solvent signal by RF pulses

Purge pulses (spin-lock- or trim pulses) $p28 = < 1.0\text{ms}$ @ $pl1$ (max. level of $pl1$ corresponding to $10\mu\text{s}$ 90 degree pulse)! use $p28$ only if you have a solvent signal to suppress, try without ($p28 = 0.1\text{ms}$) for inverse probeheads with latest design.

Example: hsqcetgpsi




TOCSY experiment: effect of trim pulses




1 trim pulse

(ppm)




2 trim pulses

(ppm)



Effect of trim pulses when presaturation is applied



Problem:

The two trim pulses are acting as B_1 -gradients.

First TP is defocussing residual water in inhomogeneous regions.

This water has not been presaturated.

Second TP is refocussing that water.

Solution:


remove first trim pulse



TOCSY experiment: effect of trim pulses



first serial file of TOCSY experiment with presaturation



Keeping the water magnetization along the +z axis




Common element in heteronuclear correlation experiments to improve WATERGATE

- **flip-back pulses** throughout the pulse sequence (prior to WATERGATE) to re-align transverse water magnetization along +z axis
- Pulse sequence names: **...fp...** with the following parameters:
 - selective on-resonant flip-back pulse:** **p11** 1 - 2ms, the shorter the broader
 - shape of the flip-back pulse:** **spnam1** sinc or rectangular pulse
 - power of the flip-back pulse:** **sp1** calculate selective 90° pulse in stdisp (might be optimized in gs-mode)

HSQC with water-flip-back pulses & WATERGATE:

hsqcfpf3gp phwg



Trim Pulses


Flip back pulses

Shaped pulses in triple resonance experiments

Shaped pulses in Triple Resonance Experiments




HNCOCA



Time reversed 90° shaped pulses



HNCOCA



hncocagp3d

(p21 ph1):f3

(p13:sp2 ph3):f2

4u

(p14:sp5 ph1):f2

DELTA3

(p14:sp3 ph1):f2

4u

(p14:sp5 ph1):f2

DELTA3

(p13:sp8 ph2):f2

time reversed

4u

30u fq=cnst23(bf ppm):f2

(p13:sp2 ph4):f2

d28

(p14:sp3 ph1):f2

d28

(p13:sp8 ph2):f2

time reversed

d0

(center (p14:sp7 ph1):f2 (p22 ph8):f3)

d0

(p14:sp3 ph1):f2

DELTA4

(p14:sp7 ph1):f2

4u

(p13:sp2 ph9):f2

d28

(p14:sp3 ph1):f2

d28

(p13:sp8 ph10):f2

time reversed

4u

30u fq=cnst21(bf ppm):f2

(p13:sp2 ph2):f2

DELTA3

(p14:sp5 ph1):f2

4u

(p14:sp3 ph1):f2

DELTA3

(p14:sp5 ph1):f2

4u

(p13:sp8 ph1):f2


time reversed

(p21 ph1):f3

Time reversed 90° shaped pulses




Time reversed pulses



Simultaneous 180° shaped pulses



HNCOCA




- (p21 ph1):f3
- (p13:sp2 ph3):f2
- 4u
- (p14:sp5 ph1):f2
- DELTA3
- (p14:sp3 ph1):f2**
- 4u**
- (p14:sp5 ph1):f2**
- DELTA3
- (p13:sp8 ph2):f2
- 4u
- 30u fq=cnst23(bf ppm):f2
- (p13:sp2 ph4):f2
- d28
- (p14:sp3 ph1):f2
- d28
- (p13:sp8 ph2):f2
- d0
- (center (p14:sp7 ph1):f2 (p22 ph8):f3)
- d0

Simultaneous 180° shaped pulses



HNCOCA




- (p21 ph1):f3
- (p13:sp2 ph3):f2
- 4u
- (p14:sp5 ph1):f2
- DELTA3
- (p14:sp3 ph1):f2
- 4u
- (p14:sp5 ph1):f2
- DELTA3
- (p13:sp8 ph2):f2
- 4u
- 30u fq=cnst23(bf ppm):f2
- (p13:sp2 ph4):f2
- d28
- (p14:sp3 ph1):f2
- d28
- (p13:sp8 ph2):f2
- d0
- (center (p14:sp7 ph1):f2 (p22 ph8):f3)
- d0

Single 180° shaped pulses



HNCOCA




(p21 ph1):f3
 (p13:sp2 ph3):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p14:sp3 ph1):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p13:sp8 ph2):f2
 4u
 30u fq=cnst23(bf ppm):f2
 (p13:sp2 ph4):f2
 d28
 (p14:sp3 ph1):f2
 d28
 (p13:sp8 ph2):f2
 d0
 (center (p14:sp7 ph1):f2 (p22 ph8):f3)
 d0

Refocused phase evolution during t_1 : 180° shaped pulses



HNCOCA



(p21 ph1):f3
 (p13:sp2 ph3):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p14:sp3 ph1):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p13:sp8 ph2):f2
 4u
 30u fq=cnst23(bf ppm):f2
 (p13:sp2 ph4):f2
 d28
 (p14:sp3 ph1):f2
 d28
 (p13:sp8 ph2):f2
 d0
 (center (p14:sp7 ph1):f2 (p22 ph8):f3)
 d0
 (p14:sp3 ph1):f2
 DELTA4
 (p14:sp7 ph1):f2

$d0*2 + \text{larger}(p14, p22) - p14$

Shaped pulses in Triple Resonance Experiments



Pairs of 90° pulses: final pulse is time-reversed



Compensation for phase evolution during t_1 : spin-echo



BS compensation for 180° pulses: pair of 180° pulses