



New Features of TopSpin and New Hardware

Presented By Mike Brown
04222012

TopSpin Versions



TopSpin 2.1.pl6

TopSpin 3.0


- Released Spring 2010
- Status: **TS3.0.pl4**

TopSpin 3.1



- Released April 2011 (ENC Conference)
- Status: TS3.1.pl2
- Minor 3.1 pl2 Upgrade

TopSpin 3.1.5

Released April 2012 (ENC Conference)



TopSpin Versions



Student licenses are now available



- For students only (student certificate may be checked)
- **inexpensive**
- limited for 3 years
- **full access to processing, simulation and plotting**
- node locked license type
- **no license transfer possible**
- cannot be upgraded to next major TopSpin version
- **can be ordered via Bruker web shop**
- has to be bought with a credit card or invoice
- **software has to be downloaded**

Student Version



Available for

- Macintosh
- Linux RedHat or CentOS
- Windows 7 or XP

contain:

- TopSpin (processing only)
- FlexLM
- NMR-Sim

You should have a fast connection to download these, The Flexlm license will be e-mailed to you




TopSpin 3.1.5 Developments



Preparing data sets:

edc / rpar

- Solvent selection, access to `getprosol`
- "Use current parameters" to transfer all parameters to new data set
- "Keep parameters" to transfer a personal selection only (configurable)
- Choose experiment via rpar dialog, link to NMRGuide available
- "Show recommended" restricts selection to standards




Selective Experiments are easy to set up and run !



Selective Acquisition Experiments:



- GUI now provides better help when setting up Selective experiments
- Supports 1D NOESY, COSY, TOCSY, Region selective 2D HMBC




AVANCE III HD

New features

- Digital RF synthesizer **SGU3**
- Enhanced Receiver **RXAD/2**
- Enhanced **HPLNA** ^1H with ^{19}F
- **HPLNA XBB** solids preamplifier
- **Enhanced Tune & Match**
- **Preamplifier user interface / Touchscreen**
- **Nanobay 'CryoProbe ready'**
- Support of **NMR Thermometer™**



AVANCE III HD - Enhanced receiver and new synthesizer electronics



HPLNA ^1H Preamplifier

- Enhanced dynamic range
- New ^{19}F capability

Receiver / Synthesizer



- RXAD2: enhanced electronics
- SGU3: 960 MHz digital up-converter

Benefits:


- Higher dynamic range
- Cleaner RF, less phase noise
- Enhanced digital RF and power setting
- Better stability & reproducibility

Sample with high concentration / protonated buffers

AVANCE III HD – Receiver electronics Enhanced Dynamic Range




AVANCE III HD




~30% less noise

x 8k




600 MHz ^1H of 100mg/ml sertraline (Zoloft) in DMSO-d6

AVANCE III HD – Enhanced ‘wobb’ Accurate Tune & Match



31P @ 500 MHz with 200 MHz span



$$r = \frac{U_r}{U_h} = \frac{Z_a - Z_l}{Z_a + Z_l} \quad |r| \leq 1$$

Benefits:

- **Precise tuning and matching**
- What you see is what you get
- Easy T&M with probes
- Factory calibrated preamplifiers

Good to know:

- May be used to check instrument

AVANCE III HD – HPPR

Advanced Preamp User Interface

- Touch Screen display
- Ready to run
- Hot Plug operation with PICS



Now available

'wobble curve' displayed at the probe

- Easy Tune & Match for solids probes
- Optimal impedance matching for probes without ATM™

AVANCE III HD – NMR Thermometer

Where to measure sample temperature?

High Resolution




H₂O: ~ 0.01ppm / K
@600MHz = 6Hz / K

RF Loading 0 to ~2K

0 – 12 Hz shifts

HR-MAS




Other heating effects

MAS spinning speed

up to 5 - 20K

AVANCE III HD – NMR Thermometer

Measure the temperature with the spins



293.0 K
Sensor: 290.7 K



Edit lock parameters for solvent "Naac50".

Lock parameters

Probe name: 5 mm PABBO BB-1H/D Z-GRD Z114607/0007

Probe description:

Lock power:

Loop gain:

Loop time:

Loop filter:

Lock Phase:

Lock power instep:

Temperature lock power:

Signals



Signal	Shift [ppm]	Relative intensity	Type	Description	Delete
1	4.7	1	Lock	<input type="text"/>	<input type="checkbox"/>
2	1.8	1	Temperature	<input type="text"/>	<input type="checkbox"/>

Temperature shift values

Value	Shift [ppm]	Temperature [K]	Delete
1	2.4811	346.16	<input type="checkbox"/>
2	2.5203	341.37	<input type="checkbox"/>
3	2.5588	336.62	<input type="checkbox"/>
4	2.6122	331.8	<input type="checkbox"/>

Sample with thermo-sensitive ^2H compound

AVANCE III HD – NMR Thermometer Topspin User Interface



Benefits NMR Thermometer:

- Measure
 - Maintain
 - Monitor
- } **temperature within the sample**

AVANCE III HD - Ultimate NMR Platform for Life Science and Materials Research



Further improved and new electronics:

- ^1H Preamplifier, Receiver & Digital RF Synthesizer
- High power Broadband LNA Preamplifier with 50 Ohm
- Touchscreen Tune & Match at the probe

Benefits:

- Higher dynamic range and 30% less noise
- Up to 35% SNR gain with the dynamic sample
- Up to 20% increase in sensitivity for X-nuclei Solid State experiments
- State of the art solids spectra without manual adjustments
- Easy and accurate experiment setup

New NMR Thermometer™:

- Accurate sample temperature
- Perfect chemical shifts match
- Measure the temperature with your spins!


CryoProbe Revolution Continues

→ over 1,100 CryoProbes installed worldwide

CryoProbe™ Prodigy TCI 600



- 2.5 times proton sensitivity boost at affordable prices (compared to RT probes)
- 2 times carbon sensitivity
- More than 6-fold throughput increase
- Minimum operating and maintenance costs, long service intervals of 2 years
- Requires no additional infrastructure: easy to site within a small footprint
- Avance III HD 600 spectrometer with Prodigy TCI surpasses conventional 950 MHz RT-probe system in proton sensitivity!
- ~35 ultra-high field 900, 950 and 1 GHz systems now installed by Bruker in customer labs worldwide.











^{15}N -HSQC spectrum of a 2mM lysozyme sample with ^{15}N at natural abundance in 50 min

WineScreener™


Cost-Effective, Automated and High Throughput Solution for Wine Quality Control



- Rapid, cost-effective quantitative targeted and non-targeted statistical analyses of wine
- Based on proven NMR *JuiceScreener*™
- Frequent wine scandals in recent years demonstrated need for new technology that assesses quality and safety within a single experiment
- First product release will offer a targeted analysis applicable to wines produced worldwide
- Non-targeted analysis developed on and applicable to German wines (for now)

Compound	Value	Unit	Official Reference		Wine-Profiling™
			Flag	max. Value	Authentic NMR Database
acetic acid	525	mg/L	○	-	165  1007
benzoic acid	<5	mg/L	●	0 mg/L	<5 mg/L in reference set
citric acid	299	mg/L	●	1000 mg/L	<200  471
ethanol	91.7	g/L	○	-	56.2  108.0
fumaric acid	<5	mg/L	○	-	<5  12
glycerol	7.1	g/L	○	-	4.5  11.1
malic acid	5.5	g/L	○	-	<0.2  8.2
methanol	53	mg/L	●	250 mg/L	6  146
shikimic acid	88	mg/L	○	-	<20  100
sorbic acid	<5	mg/L	●	200 mg/L	<5 mg/L in reference set

WineScreener sample quantification report.



White wine model: 1,3-Propanediol indicates uncontrolled organic acid decomposition

Nitrogen Liquefier BNL

More Efficient, Sustainable and Economic Magnet Operation



- Essentially zero nitrogen boil-off for Ascend™ standard-bore magnets up to 700/54
- **BNL** eliminates nitrogen refills for maximum user convenience
- No interruptions of long-term experiments
- Lower cost of ownership
- Solution to concerns about limited availability or logistics of cryogenics in some emerging markets
- Monitoring Unit ensures easy visual monitoring
- Everyone can benefit from nitrogen refrigeration (CryoProbe™ customers could already benefit from BSNL, using extra cooling capacity of latest generation CryoPlatform™)



Ascend 600, equipped with a **BNL** nitrogen liquefier

CryoProbe™ BBFO 500 MHz



- The new BBFO CryoProbe offers the highest commercially available sensitivity for the broadest range of nuclei.
- The probe is designed for both observe and inverse detection. This first fully broad banded CryoProbe can be automatically tuned over its entire range from ^{15}N to ^{19}F , as well as of course on the proton channel.
- Allows not only all X-observe ^1H decoupled experiments and their inverse detected counterparts, but also $^{19}\text{F}\{^1\text{H}\}/^1\text{H}\{^{19}\text{F}\}$ in full automation.
- In addition it provides a great solution for research applications on less commonly measured nuclei as well as low natural abundance / low sensitivity nuclei.

CryoProbe™ BBFO 500 MHz




Features

- BBFO including $^{19}\text{F}\{^1\text{H}\}$ and $^1\text{H}\{^{19}\text{F}\}$
- ^{15}N - ^{19}F S/N gain: factor $\sim 4^*$
- ^1H S/N gain: factor $\sim 3^*$
- Cold preamplifiers for all nuclei (BB/ ^{19}F / ^1H / ^2H)
- Z-gradient and ATM™ compatible
- $0^\circ\text{C} \rightarrow 80^\circ\text{C}$ sample temperature range, other ranges on request
- Available at 500 MHz, other frequencies on request

* Compared to The Smart Probe

CryoProbe™ BBFO 500 MHz




Left: ¹H-³¹P HMQC of MePt(PR₃)₂(C₄F₆H).

Right: ¹H-¹⁹⁵Pt-HMQC of MePt(PR₃)₂(C₄F₆H).

A new command has been added: tfp



A new command called tfp has been added which does a traf, ft and a pk on a spectra




New Manuals for TopSpin 3.1



- There is a new manual that explains the edlock panel. It gives a method for updating the Field value, using any solvent. It is called **Edlock Guide**
- There is a new manual for edprosol, which explains the inner workings and how to use it for different solvent, and probes. It is called the **Edprosol Manual**
- There is a new manual for diffusion, called **Diffusion**. It gives new information about how to measure diffusion and the new interface for this.
- The **Shapetool** manual has been updated.



Are there any Questions ?



www.bruker-biospin.com