


Pulse programming

Some notes and news



Content




1. Timing for change of rf parameters
2. New mc-statement
3. Shapes and power level settings

Timing for changing rf parameters phase, frequency, amplitude (shapes only)



HNCOCA

hncocagp3d



(p21 ph1):f3
 (p13:sp2 ph3):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p14:sp3 ph1):f2
 4u
 (p14:sp5 ph1):f2
 DELTA3
 (p13:sp8 ph2):f2 time reversed
 4u
 30u fq=cnst23(bf ppm):f2
 (p13:sp2 ph4):f2
 d28
 (p14:sp3 ph1):f2
 d28
 (p13:sp8 ph2):f2 time reversed
 d0
 (center (p14:sp7 ph1):f2 (p22 ph8):f3)
 d0

(p14:sp3 ph1):f2
 DELTA4
 (p14:sp7 ph1):f2
 4u
 (p13:sp2 ph9):f2
 d28
 (p14:sp3 ph1):f2
 d28
 (p13:sp8 ph10):f2 time reversed
 4u
 30u fq=cnst21(bf ppm):f2
 (p13:sp2 ph2):f2
 DELTA3
 (p14:sp5 ph1):f2
 4u
 (p14:sp3 ph1):f2
 DELTA3
 (p14:sp5 ph1):f2
 4u
 (p13:sp8 ph1):f2 time reversed
 (p21 ph1):f3

mc-statement

Introduction



FnMODE parameter =

FQ	phase sensitive 1D
F1QF	magnitude mode QF
F1PH, F2PH	QSEQ, TPPI, States or States-TPPI
F1EA, F2EA	phase sensitive echo-antiecho

- Simplifies pulse programming
 - includes
 - disk write (wr)
 - file pointer incrementation (if)
 - memory initialization (zd)
 - expanded loop structure possible
 - no need to control actions required for phase sensitive 2D experiments (phase or delay incrementation) with delays
- One pulse program can be used for different 2D phase modes

mc-statement

Original implementation, 1D experiments



without mc:

```
1 ze
2 d1
  p1 ph1
  go=2 ph31
  wr #0
exit
```

with mc:

```
1 ze
2 30m
  d1
  p1 ph1
  go=2 ph31
  30m mc #0 to 2 F0(zd)
exit
```

Expanded with mc:

```
define delay MCWRK
define delay MCREST
"MCWRK=0.333333*30"
"MCREST=30m-30m"
```

```
1 ze
2 MCWRK*3
  LBLF0, MCREST
  d1
  p1 ph1
  go=2 ph31
  MCWRK wr #0
  MCWRK ze
  MCWRK zd
  lo to LBLF= times td0
exit
```

mc-statement

Original implementation, 2D experiments



up to XWINNMR 2.6

```
;noesytp
```

```
#include <Avance.incl>
```

```
"d0=3u"
```

```
1 ze
```

```
2 d1
```

```
3 p1 ph1
```

```
  d0
```

```
  p1 ph2
```

```
  d8
```

```
  p1 ph3
```

```
  go=2 ph31
```

```
  d1 wr #0 if #0 ip1 id0 zd
```

```
  lo to 3 times td1
```

```
exit
```

original mc:

```
;noesyph
```

```
#include <Avance.incl>
```

```
"d0=3u"
```

```
1 ze
```

```
2 d1
```

```
3 p1 ph1
```

```
  d0
```

```
  p1 ph2
```

```
  d8
```

```
  p1 ph3
```

```
  go=2 ph31
```

```
  d1 mc #0 to 2 F1PH(ip1, id0)
```

```
exit
```

mc-statement

Original implementation, 2D experiments



F1QF
F1PH
F1EA

phase insensitive
phase sensitive
Echo-Antiecho

QF
QSEQ, States, TPPI, States-TPPI
Echo-Antiecho

The acquisition and processing parameter FnMODE define the mode for F1 (and F2, for 3D-experiments) incrementation

For details, see XWINNMR help -> Other topics -> Writing pulse programs

MC clause	t1 quadrature mode	action after	
		odd increment	even increment
F1PH(ip1, id0) F1PH(ip1, id0)	TPPI States-TPPI	ip1+ id0 ip1	ip1 + id0 again id0
F1PH(rd10 & rd30 & ip4, id0) F1PH(rd10 & rd30 & ip4, id0)	TPPI States-TPPI	rd10+rd30+ip4+ rd10+rd30+ip4	id0 id0

mc-statement

Challenge for non-uniform sampling



Explicit programming of a 2D experiment with non-uniform sampling

define list<loopcounter> t1list=<\$VCLIST> variable counter list for sampling scheme

1 ze

.....

4 d1

"d0=3u+in0*t1list"

calculate value of t_1 delay d0

3m t1list.inc

.....

go=2 ph31 cpd2:f2

d11 do:f2 wr #0 if #0 zd

"cnst30=(t1list%2)*180"

calculate phase increments

"cnst31=cnst30"

3m ip6+cnst30

3m ip9+cnst30

.....

mc-statement

New implementation for non-uniform sampling



	<i>Old</i>	<i>New</i>
Phase	ip1	calph(ph1, +90)
Delay	id10	caldel(d10, +id10)
Loop counter	iu0	calclc(l0, 1)
Phase reset	n.a.	exec(rppall)

Execution of conventional commands:

Incrementation of a frequency list defined by: `define list<frequency stdlist=<$FQ2LIST>`

exec(stdlist.inc)

Lists

calclist(name, incr)

mc-statement

changes for non-uniform sampling

new mc (TopSpin 3.0 +)

increment

decrement

d11 do:f3 mc #0 to 2

F1PH(calph(ph4, +90), caldel(d0, +in0))

F2PH(calph(ph5, +90), caldel(d10, +in10) & caldel(d29, +in29) & caldel(d30, -in30))

mc-statement

Important note



Do not mix old and new mc-statements

Shaped pulses and power level settings



Power level switching before and after shaped pulses

d16 pl0:f1

(p12:sp1 ph2:r):f1

shaped pulse

4u

d12 pl1:f1

p2 ph3

hard pulse