

# RF power handling and more



Typical Pulses for Triple Resonance Inverse 5 mm H-C/N-D CryoProbes™ at 850 MHz

1. PURPOSE ..... 3

2. SCOPE ..... 3

3. REFERENCE TO DOCUMENTS ..... 3

4. IMPLEMENTATION ..... 4

4.1 INTRODUCTION ..... 3

4.2 PULSE AND POWER RECOMMENDATIONS ..... 4

4.3 NOTES ..... 5

Technical Customer Document  
ZTSC0021, Index 00

Bruker BioSpin AG  
BBI0

Typical Pulses for Triple Resonance Inverse 5 mm H-C/N-D CryoProbes™ at 850 MHz

4.2 Pulse and power recommendations

Note: all values require a minimum relaxation delay  $d1$  of 1 s unless differently specified in footnotes <sup>2)</sup>. Important: observe the separate notes on the following pages!

5 mm 850 MHz TCI H-C/N-D CryoProbe™	
<sup>1</sup> H	
hard pulse <sup>3)</sup> (max. length 1ms)	8.0 μs
hard pulse for lossy samples	Power level corresponding to 8.0 μs pulse for non-lossy sample
trim pulse p28 <sup>4)</sup>	1 ms @ 8 μs 2 ms @ 11 μs
TOCZY spin lock <sup>5)</sup>	120 ms @ 20 μs 400 ms @ 35 μs
ROESY spin lock	Up to CW for power level corresponding to a 80 μs pulse
WALTZ65 decoupling during <sup>13</sup> C- detection	Up to CW for power level corresponding to a 80 μs pulse
DIPS12-decoupling in triple resonance	400 ms @ 35 μs
<sup>13</sup> C	
hard pulse <sup>3)</sup> (max. length 1000°)	12.0 μs
trim pulse <sup>4)</sup>	2 ms @ 22 μs 40 ms @ 22 μs
CC spin lock <sup>5)</sup>	150 ms @ 12 μs up to 2 kHz
CPMG REX	250 ms @ 45 μs (> 142 ppm bandwidth)
GARP4 decoupling <sup>7)</sup>	500 ms @ 55 μs
Adiabatic decoupling	Crp48.1.5.20.2 (> 170 ppm bandwidth)
selective pulses <sup>8)</sup>	Q2: 290 μs / Q5: 226 μs Q3: 180 μs CHIRP: 2 ms @ 22 μs
<sup>15</sup> N	
hard pulse <sup>3)</sup> (max. length 1000°)	32.0 μs
GARP4 decoupling <sup>7)</sup>	125 ms @ 120 μs (> 134 ppm bandwidth) 250 ms @ 170 μs (> 95 ppm bandwidth) 1000 ms @ 350 μs (> 45 ppm bandwidth; sofast-HMQC)
CPMG T2 <sup>9)</sup>	250 ms @ 40 μs (see warning <sup>9)</sup> ) 400 ms @ 50 μs
CPMG REX <sup>10)</sup>	50 ms @ 40 μs up to 2 kHz
T1 rho <sup>11)</sup>	100 ms @ 40 μs up to 1 kHz 100 ms @ 3 kHz 200 ms @ 1.5 kHz
<sup>1</sup> H	
hard pulse (max. length 1000°)	100 μs
WALTZ64 decoupling	100 ms @ 250 μs
Z-Gradient	
Absolute max. current $I_{max}$	10 A
Max. overall length $T_{max}$ <sup>12)</sup>	10 ms @ 10 A

Table 1

Technical Customer Document  
ZTSC0021, Index 00

Bruker BioSpin AG  
BBI0

4 / 6

What does “power handling” mean?

1. rf pulse deposit energy (heat) in the rf circuit and rf coils
2. Long pulses trains deposit more heat
3. Examples for long pulse trains:
  - CPD decoupling
  - Long CW irradiation:  $T_{1\rho}$  experiments
  - Long CPMG sequence:  $T_2$ , REX

What might happen?

1. If deposited energy exceed specifications:
  - Excessive sample heating
  - Detuning of the rf circuit
  - Pulse width change: longer pulses which are difficult to control
  - Poor quality spectra

Will the CryoProbe be destroyed?

1. Power check of Topspin: ensures that the pulse voltage is not exceeded
  - Too high voltage **WILL** destroy the probe
  - Applying long pulses just close below the voltage limit are critical


Are there also be specifications for the allowed gradient power?

1. Yes
2. Why?
  - DC current is applied to a coil
  - The gradient coil behaves like a heating coil
  - DC current applied to long = long intense gradient pulses will overheat the wired
  - Wires will burn

## Sofast-HMQC

PC9

R-Snob




# RF power handling



## Sofast-HMQC

1mM Ubiquitin ( $^{13}\text{C}/^{15}\text{N}$ ) in 90% H<sub>2</sub>O / 10% D<sub>2</sub>O

**50 seconds acquisition time !**



Acquisition Parameter:

AQ: 42 ms

D1: 8 ms

$^{15}\text{N}$  decoupling:  $\gamma B_1$  of 892 Hz



# RF power handling



	5 mm 850 MHz TCI H-C/N-D CryoProbe™
<sup>1</sup> H	
hard pulse <sup>3)</sup> (max. length 1ms)	8.0 μs
hard pulse for lossy samples	Power level corresponding to 8.0 μs pulse for non-lossy sample
trim pulse p28 <sup>4)</sup>	1 ms @ 8 μs 2 ms @ 11 μs
TOCSY spin lock <sup>5)</sup>	120 ms @ 20 μs 400 ms @ 35 μs
ROESY spin lock	Up to CW for power level corresponding to a 80 μs pulse
WALTZ65 decoupling during <sup>13</sup> C- detection	Up to CW for power level corresponding to a 80 μs pulse
DIPSII2-decoupling in triple resonance	400 ms @ 35 μs



# RF power handling



5 mm 850 MHz TCI H-C/N-D CryoProbe™	
<sup>1</sup> H	
<sup>13</sup> C	
hard pulse <sup>6)</sup> (max. length 1000°)	12.0 μs
trim pulse <sup>5)</sup>	2 ms @ 22 μs
CC spin lock <sup>5)</sup>	40 ms @ 22 μs
CPMG REX	150 ms @ 12 μs up to 2 kHz
GARP4 decoupling <sup>7)</sup>	250 ms @ 45 μs (> 142 ppm bandwidth)
Adiabatic decoupling	500 ms @ 55 μs Crp48,1.5,20.2 (> 170 ppm bandwidth)
selective pulses <sup>8)</sup>	G4: 290 μs / Q5: 226 μs Q3: 180 μs CHIRP: 2 ms @ 22 μs

# RF power handling



5 mm 850 MHz TCI H-C/N-D CryoProbe™	
<sup>1</sup> H	
<sup>13</sup> C	
<sup>15</sup> N	
hard pulse <sup>6)</sup> (max. length 1000°)	32.0 μs
GARP4 decoupling <sup>7)</sup>	125 ms @ 120 μs (> 134 ppm bandwidth) 250 ms @ 170 μs (> 95 ppm bandwidth) 1000 ms @ 350 μs (> 45 ppm bandwidth; sofast-HMQC)
CPMG T2 <sup>9)</sup>	250 ms @ 40 μs (see warning <sup>9)</sup> ) 400 ms @ 50 μs
CPMG REX <sup>10)</sup>	50 ms @ 40 μs up to 2 kHz 100 ms @ 40 μs up to 1 kHz
T1 rho <sup>11)</sup>	100 ms @ 3 kHz 200 ms @ 1.5 kHz

# Gradient pulses




Z-Gradient		
Absolute max. current <sup>10)</sup>	10 A	10 A
Max. overall length <sup>11)</sup>	10 ms @ 10 A	10 ms @ 10 A

Note: some pulse sequences might have more gradients than you expect!

- Experiment with bipolar gradients used for diffusion/DOSY experiments

# Gradient pulses DOSY

Double STimulated Echo (DSTE) for convection compensation



- LED also
- 3 spoil gradients
- 4 diffusion gradients
- same, but with bipolar gradients:
- LED and 3 spoil gradients (G7, G8, G9)
- bipolar gradients for diffusion (G6), makes 8 diffusion gradients!!!

**NOTE: gradient limitation for CryoProbes: maximum 10ms@10A**

example 1. SMSQ10.100 gradients with 2ms length, 100%:  $5 \times 2\text{ms} = 10\text{ms}$

example 2. SINE.100 gradients with 2ms length, 100%:  $8 \times 2\text{ms} \times 0.67 = 10.72\text{ms}$

because the area of a sine curve is 67% of the area of a square