

Avance 3D / Triple Resonance

Introduction to

3-dimensional and triple resonance

NMR-spectroscopy

with focus on biomolecules

on Avance spectrometers

The information in this course may be altered without notice. BRUKER accepts no responsibility for actions taken as a result of use of this course. BRUKER accepts no liability for any mistakes contained in the course, leading to coincidental damage, whether during installation or operation of the instrument.

Unauthorised reproduction of course contents, without written permission from the publishers, or translation into another language, either in full or in part, is forbidden.

This course manual was written
And desktop published by

Helena Kovacs (Application)

© September 1, 1999, Bruker AG
Fällanden, Switzerland
P/N Z31525
DWG-No. 1232001

INDEX

INTRODUCTION	1	7
Goals of the course 1.1		7
From 2D to 3D 1.2		8
The basic principle		8
Acquisition parameter display		9
Processing parameter display		12
Standard parameter sets for 3D 1.3		13
Manual set-up		13
Automated set-up		13
Choice of the pulse sequence 1.4		13
Strategies for protein NMR studies 1.5		15
Biomolecular sample preparation 1.6		18
PULSE CALIBRATION	2	19
RF routing 2.1		19
¹H observe pulse 2.2		20
¹³C inverse mode 2.3		20
¹⁵N inverse mode 2.4		20
¹³C selective excitation/inversion, inverse mode 2.5		20
²H decoupling pulse 2.6		23
Delays 2.7		24
Summary of the calibrated pulse lengths 2.8		25
Preparatory 1D and 2D experiments 2.9		26
¹H-¹³C TOCSY-HSQC	3	28
Introduction 3.1		28
Pulse Sequence Diagram 3.2		29
Data Acquisition 3.3		30
Setting up the experiment		30
Recording the two 2D planes 3.4		32
Inspection of 2D planes		32
Phase corrections		32
Window functions		33
Spectral Processing 3.5		33
Processing the 3D spectrum		33
AVANCE 3D	Bruker	3

Additional processing parameters	3.6	35
Linear prediction	3.7	36
¹H-¹⁵N NOESY-HSQC	4	37
Introduction	4.1	37
Pulse Sequence Diagram	4.2	38
Data Acquisition	4.3	39
Setting up the experiment		39
The „GS“ Interactive mode to optimize acquisition parameters		40
Spectral Processing	4.4	41
HNHA	5	43
Introduction	5.1	43
Pulse Sequence Diagram	5.2	44
Data Acquisition	5.3	45
„GS“ Interactive mode to optimize acquisition parameters		46
Broadband, adiabatic and selective ¹³C decoupling	5.4	46
Broadband decoupling		46
Adiabatic decoupling		47
Selective decoupling		47
Spectral Processing	5.5	47
HNCA	6	49
Introduction	6.1	49
Pulse Sequence Diagram	6.2	50
Data Acquisition	6.3	51
Adding ² H decoupling		53
Spectral Processing	6.4	53
HN(CO)CA	7	55
Introduction	7.1	55
Pulse Sequence Diagram	7.2	56
Data Acquisition	7.3	57
Creating the ¹³ C frequency list		59
Spectral Processing	7.4	59
HCCH-TOCSY	8	61
Introduction	8.1	61

Pulse Sequence Diagram	8.2	62
Data Acquisition	8.3	63
Spectral Processing	8.4	65
CP-DRIVEN HCCH-TOCSY	9	67
Introduction	9.1	67
Pulse Sequence Diagram	9.2	68
Data Acquisition	9.3	69
Spectral Processing	9.4	70
H(CC)(CO)NH	10	72
Introduction	10.1	72
Pulse Sequence Diagram	10.2	73
Data Acquisition	10.3	74
Creating the ¹ H and ¹³ C frequency lists		76
Spectral Processing	10.4	76
CBCA(CO)NH	11	78
Introduction	11.1	78
Pulse Sequence Diagram	11.2	79
Data Acquisition	11.3	80
Creating the ¹³ C frequency list		82
Spectral Processing	11.4	82
HBHA(CO)NH	12	84
Introduction	12.1	84
Pulse Sequence Diagram	12.2	85
Data Acquisition	12.3	86
Creating the ¹ H and ¹³ C frequency lists		88
Setting up the selective decoupling on the alpha-carbons		88
Spectral Processing	12.4	89
APPENDICES	13	91
Coupling constants in polypeptides	13.1	91
¹³C and ¹H chemical shifts in different residues	13.2	92
Sensitivity of triple resonance experiments	13.3	94
AVANCE 3D	Bruker	5

FIGURES

Fig. 1.	The 3D pulse sequence.	8
Fig. 2.	The 3D spectrum.	9
Fig. 3.	The acquisition parameter display.	10
Fig. 4.	The processing parameter display.	11
Fig. 5.	Choice of the pulse sequence.	13
Fig. 6.	NMR-structure of a 120-residue protein.	18
Fig. 7.	RF-routing.	19
Fig. 8.	Shaped 90° - ^{13}C -pulse calibration (anti-phase signal).	21
Fig. 9.	Shaped 180° - ^{13}C -pulse calibration (in-phase signal).	22
Fig. 10.	Deuterium pulse calibration.	24
Fig. 11.	Uncoupled ^1H 1D-spectrum of a ^{15}N , ^{13}C -labeled protein.	26
Fig. 12.	Folding and aliasing of the ^{13}C -dimension.	26
Fig. 13.	Example of aliasing of the ^{13}C -dimension in the ^{13}C -HSQC spectrum.	26
Fig. 14.	^1H - ^{13}C TOCSY-HSQC.	28
Fig. 15.	The ^1H - ^1H and ^1H - ^{13}C planes of the ^1H - ^{13}C TOCSY-HSQC spectrum.	34
Fig. 16.	^1H - ^{15}N NOESY-HSQC.	37
Fig. 17.	The ^1H - ^1H and ^1H - ^{15}N planes of the ^1H - ^{15}N NOESY-HSQC spectrum.	41
Fig. 18.	HNHA.	43
Fig. 19.	The ^1H - ^1H plane of HNHA.	47
Fig. 20.	Magnetization transfer in the HNCA experiment.	49
Fig. 21.	HNCA.	49
Fig. 22.	The ^1H - ^{13}C and ^1H - ^{15}N planes of the HNCA spectrum.	54
Fig. 23.	Magnetization transfer in HN(CO)CA.	55
Fig. 24.	HN(CO)CA.	55
Fig. 25.	The ^1H - ^{13}C and ^1H - ^{15}N planes of the HN(CO)CA spectrum.	60
Fig. 26.	HCCH-TOCSY.	62
Fig. 27.	The ^1H - ^1H ($l1=3$) and ^1H - ^{13}C planes of the HCCH-TOCSY spectrum.	66
Fig. 28.	CP-driven HCCH-TOCSY.	68
Fig. 29.	The CP-driven HCCH-TOCSY spectrum.	71
Fig. 30.	Magnetization transfer in H(CC)(CO)NH.	72
Fig. 31.	H(CC)(CO)NH.	73
Fig. 32.	The ^1H - ^1H and ^1H - ^{15}N planes of the H(CC)(CO)NH spectrum.	77
Fig. 33.	The magnetization transfer in CBCA(CO)NH.	78
Fig. 34.	CBCA(CO)NH.	79
Fig. 35.	The ^1H - ^{13}C and ^1H - ^{15}N planes of the CBCA(CO)NH spectrum.	83
Fig. 36.	Magnetization transfer in HBHA(CO)NH.	84
Fig. 37.	HBHA(CO)NH.	85
Fig. 38.	The ^1H - ^1H and ^1H - ^{15}N planes of the HBHA(CO)NH spectrum.	89
Fig. 39.	One-bond and two-bond coupling constants in the polypeptide backbone.	91
Fig. 40.	The random coil ^1H and ^{13}C chemical shifts in the natural amino acids.	92
Fig. 41.	The random coil ring ^1H and ^{13}C chemical shifts in the natural amino acids.	93
Fig. 42.	Sensitivity of triple resonance experiments.	94

Introduction

1

Goals of the course

1.1

The present course will provide you with basic knowledge and tools in 3-dimensional NMR spectroscopy, and make you acquainted with a number of 3D-experiments. The course is designed to serve as a general introduction with focus on the implementation on your Bruker instrument. We have achieved our goal if you become confident with setting-up and optimizing the numerous 3D-experiments available in the Bruker standard pulse sequence library and also feel encouraged to embark upon modifying existing ones and writing additional pulse programs for your own specific purposes.

Although the scientific literature abounds in various 3D-experiments, we have restricted the course to a representative selection of them. The experiments are recorded on unlabeled, single-, double- and triple-labeled samples in order to highlight the particularities of the experimental set-up in different cases.

Usually when working with biomolecular samples the experimental sensitivity is of major concern. Although a multitude of experiments are available, the concentration, molecular size and special characteristics of the sample should be taken into consideration when choosing an NMR investigation protocol. It is advisable to proceed step by step, starting with the most sensitive experiments.

The topics of the present course cover the setting-up and processing of 3D experiments, triple resonance spectroscopy, calibration of the necessary pulses including pulses for selective excitation, applying broadband or adiabatic decoupling schemes and deuterium decoupling. Further, heteronuclear cross-polarization and optimization of water suppression are included, as well as the basics of Bruker pulse programming.

We wish you great pleasure and success in using 3D NMR as a tool in your research.

The application team of Bruker Switzerland.


From 2D to 3D

1.2

The basic principle


A 3-dimensional experiment can be regarded as a combination of two 2-dimensional experiments.

Fig. 1. The 3D pulse sequence.



In the 3D experiments there are two evolution times that are incremented, t_1 and t_2 (the acquisition time is called t_3). During the t_1 -evolution time (following the initial 90° - ^1H pulse which creates a $-H_y$ -magnetization) the magnetization is on the proton nuclei and becomes labeled by the ^1H -chemical shift. Thereafter the TOCSY mixing of the ^1H -coherences takes place, followed by an INEPT-transfer of magnetization to the ^{13}C -nuclei. During the t_2 -evolution time, the ^{13}C -coherence, which is anti-phase with respect to the protons (in terms of the product operator H_2C_x , the ^1H -magnetization is aligned along z-axis at this point) becomes labeled by the ^{13}C -frequency. Finally the magnetization of interest is transferred back to the protons for observation during t_3 by a reverse-INEPT scheme.

Fig. 2. The 3D spectrum.



After a Fourier transform of each of the three time domains the spectrum has three frequency domains, that is, F1(^1H), F2(^{13}C) and F3(^1H).

In the F1-F3 dimension the individual planes correspond to TOCSY-spectra at a certain ^{13}C -frequency (note that only ^{13}C -bound protons give rise to signals in this spectrum). In the F2-F3 dimension the individual planes correspond to the $^1\text{H}, ^{13}\text{C}$ -HSQC-planes. You can view the individual 2D-planes in XWINNMR by using the command **display** or scan through the planes by using the **scan** option.

Acquisition parameter display

In order to set-up a 3D-experiment you need to first change **parmode** to 3D. Now the acquisition parameter display, evoked by the command **eda**, consist of three columns which correspond to the F3, F2 and F1 dimensions in this very order.

Fig. 3. The acquisition parameter display.

Parameter	F3	F2	F1	Unit
PULPROG	mleviief2gp3			
AQ_mod	DQD			
TD	2048	64	128	
PARMODE	3D			
NS	4			
DS	128			
D	** Array **			sec
P	** Array **			usec
NDO		4		
INO		0.00003471		sec
ND10	2			
IN10	0.00004418			sec
SW	11.9705	75.0000	12.0000	ppm
SWH	7183.908	11318.125	7201.594	Hz
FIDRES	3.507768	176.845703	56.262451	Hz
FW	125000.00			Hz

Buttons: SAVE, Parameter, Next, CANCEL

NUCLEI: 1H 13C 1H

You need to enter the nuclei, their sweep widths, offsets and the number of points in each dimension. The increments, d0 in F1 and d10 in F2, respectively, are deduced automatically.

There is an additional parameter for each indirect dimension, called nd0 in F1 and nd10 in F2, that you need to set. This parameter indicates the number of times each increment (d0 and d10, respectively) actually occurs in the pulse program. In the above example, the value of nd10 is 2, because each evolution time is split into two parts by a 180°-refocusing pulse on the other nucleus. There is one exception, however; the TPPI-mode requires doubling of the nd-parameter value. Thus in the above example nd0 is set to 4. For the Bruker standard pulse sequences the appropriate nd-values can be found at the end of the pulse programs.

Note #1

As a matter of fact, the choice of denoting the two indirect evolution times by t_1 and t_2 (and the corresponding frequency domains by F1 and F2) is arbitrary and interchangeable. The choice has actually been made by the pulse programmer. In case of uncertainty a quick look at the pulse program (it can be evoked by the **edcpul** command as a text file or by the **ppg** command in visual display) reveals which nucleus evolves during which evolution time.

The following example shows the beginning of the pulse sequence mleviief2gp3d. The delay d0 follows (note that it is split into two by a 180°-carbon pulse) right after the initial 90°-proton pulse. The first pulse brings the proton equilibrium magnetization the transverse plane, where it becomes labeled by the ^1H chemical shift during t_1 -evolution time. Thus F1 is the indirect proton dimension. Consequently F2 must be the carbon-13 frequency domain.

Ex.	d1 pl1:f1 pl2:f2	recycle delay, setpower levels
	(p1ph1)	initial 90°- ^1H pulse
	d0	first half of in0 (initially set equal to d0)
	(p4 ph0):f2	180°- ^{13}C refocusing(decoupling) pulse
	d0	second half of in0 (initially set equal to d0)
	(p1 ph10)	90°- ^1H pulse

Note #2

Although not needed for the correct set-up of the experiment, you may note that the order in which the two evolution times are incremented is also arbitrary, and it can be found out by inspecting the pulse program.

The following example shows the end of the mleviief3gp3d sequence containing the loops for incrementation. Here, d10 (belonging to F2, the carbon-13 frequency domain) occurs first, the meaning of the command id10 is „increment delay d10“. After a complex F2-F3 plane (that is, an HSQC plane) has been recorded, the delay d10 is reset by the command rd10. Then d0 (belonging to F1, the indirect proton domain) is incremented by the command id0. Thus, in this sequence a HSQC-plane is recorded for each proton increment.

Ex.	go=2 ph31 cpd2:f2	acquisition, decoupling, go to 2
	d1 do:f2 wr #0 if #0 zd	stop decoupling, write on disc
	3m igrad EA	echo-antiecho gradient selection
	lo to 3 times 2	go to 3, repeat the loop twice
	d11 id10	increase increment d10
	3m 2*ip3	180° change of phase 3
	3m 2*ip6	180° change of phase 6
	3m 2*ip31	180° change of receiver phase
	lo to 4 times l13	go to 4, repeat l13=td2/2 times => E/A
	10m rd10 ip7	reset in10, increase phase 7
	d11 id0	increase increment in0
	lo to 5 times td1	go to 5, repeat td1 times => TPPI

Processing parameter display

Fig. 4. The processing parameter display.

Parameter	F3	F2	F1	Unit
SI	2048	128	256	
PPARMOD	3D			
SF	600.1300000	150.9027490	600.1300000	MHz
OFFSET	10.685	74.503	10.700	ppm
SR	0.00	0.00	0.00	Hz
HZpPT	3.507768	88.426712	28.132877	Hz
MC2		echo-antiecho	TPPI	
AQORDER	3-2-1			
WDW	QSINE	QSINE	QSINE	
SSB	3.3	2	2	
LB	1.00	0.30	0.30	Hz
GB	0	0.1	0.1	
PH_mod	pk	no	pk	
PKNL	TRUE			
PHC0	-64.000	0.000	55.000	degrees
PHC1	33.000	0.000	-113.000	degrees

SAVE Parameter Next CANCEL

NUCLEI: 1H 13C 1H

The command **edp** evokes a display processing parameters in XWIN-NMR. These are discussed in detail in Chapter 3 (processing of the ^1H , ^{13}C -TOCSY-HSQC experiment). The order of the dimensions is the same as for the acquisition parameters.

Standard parameter sets for 3D

1.3

Manual set-up

XWIN-NMR 2.5 and later versions provide an ever increasing number of ready parameter sets for various 3D experiments. You can simply load a set by typing the command **rpar**. Check the RF-routing through the **edsp** command. Thereafter examine the number of scans, increments and sweep widths in the **eda** display. In particular, you need to set the appropriate pulse lengths and power levels manually. This can conveniently be done through the **ased** listing, which displays the relevant parameters only. The standard Bruker pulse programs contain information about the meaning and appropriate setting of the parameters. Use the **edcpul** command to view the pulse program in a texteditor mode or type **ppg** for a visual display.

Automated set-up

An alternative way of setting up 3D triple resonance experiments is to employ two au-programs available in XWIN-NMR Version 2.5 and later. The command **xau setpulses** evokes a program that will interactively request you to enter all the pulse lengths and power levels necessary for any triple resonance experiment in the Bruker standard pulse sequence library. This au-program will also request you the names of the selective pulses you wish to use, and the proton offset. Note that from the proton offset all the other offsets are derived automatically. The information is stored in a new parameter set called **SETUP3D.1**, where the extension indicates the current probe head number.

Now you can create a data set for your chosen experiment by typing, for instance, **xau setpar3d pulprog <hncagp3d> td1 <64> td2 <32> ns <16>** where the brackets denote your own entries. You will be requested to enter the experiment number where you want to store the parameter set. The rest of the set-up is completed automatically.

Choice of the pulse sequence


1.4

The Bruker standard pulse program and parameter set libraries in XWIN-NMR Version 2.5 contain a number of commonly used 3D-experiments. The diagram below illustrates one possible classification of these. The „double resonance“ experiments require a ^{15}N - or ^{13}C -labeled sample or an unlabeled sample at high concentration. This is essential particularly if nuclear Overhauser (NOE) type of interactions are to be monitored.

The „triple resonance“ experiments require ^{15}N - and ^{13}C -labeling. Often several improved versions of the experiments are found in literature. When different implementations are included in the Bruker standard pulse program library, they are denoted with an increasing extension, for instance, **hncagp3d.1**, **hncagp3d.2** etc. The triple resonance experiments are also

available with ^2H -decoupling. Furthermore, XWIN-NMR Version 2.6 provides the TROSY-implementations of the common triple resonance experiments.

Fig. 5. Choice of the pulse sequence.



In the rightmost column of the diagram above, a selection of pulse sequences is given in order to illustrate some of the different features that have been implemented. The asterisk denotes the sequences that are used in the current 3D manual. A summary of the Bruker pulse sequence nomenclature is found on your spectrometer in the file /u/exp/stan/nmr/lists/pp/Pulprog.info. Below is a key to decipher the names in the diagram.

mlev = MLEV TOCSY mixing
pr = presaturation
ii = inverse / HSQC
et = echo-antiecho-TPPI
si = sensitivity enhanced
gp = gradient pulse syntax
f2 = heteronucleus on channel f2
2h = ^2H decoupled
tr=TROSY

dips = DIPSI-2 TOCSY mixing
19=3-9-19 WATERGATE
i4 = inverse with 4 pulses / HMQC
ea = echo-antiecho
sm = simultaneous
gs = gradient program syntax
f3 = heteronucleus on channel f3
3d = 3-dimensional
noes = NOESY

Strategies for protein NMR studies**1.5**

Table 1 proposes another classification of the common NMR experiments, now according to their applicability on different kinds of protein samples. Depending on the molecular size and the degree of labeling, a strategy is chosen in order to obtain the sequential and sidechain assignment and conformational constraints (NOEs, coupling constants and chemical shift constraints) needed for a structure elucidation. The experimental sensitivity should be considered.

For example, compare the two complementary pairs of experiments, HNCO and HN(CA)CO on one hand and HNCA and HN(CO)CA on the other. The former pair consists of a sensitive and an insensitive experiment, whereas both experiments in the latter pair are reasonably sensitive. Besides, the latter pair yields the alpha-carbon assignments which are likely to be more useful than the carbonyl-carbon assignments from the former pair. This is also why we have chosen to include the HNCA and HN(CO)CA experiments in the current manual. The pulse sequences for HNCO and HNCA are, by the way, identical except for the interchange of alpha-carbon and carbonyl frequencies. The differences in sensitivity are due to differences in the scalar coupling constants (and the corresponding delays) responsible for the magnetization transfer.

Table 1: Strategies for protein studies

Protein/Size	Experiment	Information obtained	Sensitivity
<i>Unlabeled/ less than 50.a.a.</i>	<i>2D Homonuclear</i>		
	COSY, TOCSY	intra-residue assignments	
	NOESY	sequential connectivities NOE distance constraints $^3J_{\text{HN}\alpha}$ coupling constants	
	E.COSY	$^3J_{\text{H}\alpha\beta}$ coupling constants	
<i>^{15}N-labeled/ ~ 50-80 .a.a.</i>	<i>3D Double resonance</i>		
	^{15}N -TOCSY	intra-residue assignments	
*	^{15}N -NOESY	sequential connectivities NOE constraints	
*	^{15}N -HNHA	$^3J_{\text{HN}\alpha}$ coupling constants	
	<i>or 2D HMQC-J</i>	$^3J_{\text{HN}\alpha}$ coupling constants	
	^{15}N -HNHB	$^3J_{\text{H}\alpha\beta}$ coupling constants	

INTRODUCTION

Protein/Size	Experiment	Information obtained	Sensitivity
¹³C, ¹⁵N-labeled/ ~ 80-150 .a.a.	3D Double resonance	NB. Possibly fractionally ²H-labeled	
*	¹⁵ N-NOESY	NOE constraints	
*	¹⁵ N-HNHA	³ J _{HNα} coupling constants	
	¹⁵ N-HNHB	³ J _{Hαβ} coupling constants	
	¹³ C HCCH-COSY	intra-residue assignments	
*	¹³ C HCCH-TOCSY	intra-residue assignments	
	¹³ C NOESY	sidechain NOE constraints	
	3D Triple resonance		
	HNCO	sequential connectivity	100 inter
	HN(CA)CO	sequential connectivity (combine with HNCO)	13/4 inter/intra
*	HNCA	sequential connectivity ¹³ C ^α chemical shift constraints	50/15 intra/inter
*	HN(CO)CA	(combine with HNCA)	71 inter
	HNCAH	sequential connectivity	
	HN(CO)CAH	(combine with HNCAH)	
*	CBCA(CO)NH	sequential connectivity ¹³ C ^α and ¹³ C ^β chemical shifts	13/9 ¹³ C ^α / ¹³ C ^β intra
	CBCANH	for smaller proteins (combine with CBCA(CO)NH)	4/1.7 ¹³ C ^α / ¹³ C ^β intra
	HNCACB	for bigger proteins (combine with CBCA(CO)NH)	1.3/0.5 ¹³ C ^α / ¹³ C ^β intra
*	HBHA(CO)NH	¹ H ^α and ¹ H ^β assignments	13/9 ¹ H ^α / ¹ H ^β intra
*	H(CCCO)NH	sidechain ¹ H assignments	
	(H)CC(CO)NH	sidechain ¹³ C assignments	
¹³C, ¹⁵N, ²H-label. >160 .a.a.	3D Triple resonance with ²H-decoupling		
	CT-HNCA	sequential connectivity	
	HN(CO)CA	(combine with HNCA)	
	CBCA(CO)NH	sequential connectivity ¹³ C ^α / ¹³ C ^β chemical shifts	
	CT-HNCACB	(combine with CBCA(CO)NH)	
	C(CO)NH	sidechain ¹³ C assignments	
	¹⁵ N-HSQC-NOESY- HSQC	sequential and long-range NH-NH NOE constraints	

The experiments presented in this manual are denoted by an asterisk.

Notes:

The sensitivity is given for both the inter- and intra-residual signals where appropriate (the experiments HNCA, HN(CA)CO) and for both the α - and β -position signals where appropriate (the experiments CBCA(CO)NH, CBCANH, HNCACB, HBHA(CO)NH).

Note that in the triple resonance experiments the signals are separated according to their backbone amide frequencies. This enables a „ ^{15}N -strip“ analysis of the spectra.

For ^{15}N , ^{13}C - and fractionally ^2H -labeled proteins 4-dimensional heteronuclear correlated NOESY techniques ($^{15}\text{N}/^{13}\text{C}$, $^{13}\text{C}/^{13}\text{C}$ or $^{15}\text{N}/^{15}\text{N}$) might become useful.

The TROSY-modifications of triple resonance experiments have been developed for large biomolecules with broad lines monitored at high magnetic fields, see Salzmänn *et al.* (1999).

Special experiments have been designed for NMR-investigations of labeled nucleic acids. More information about these is given in the references by Dieckmann & Feigon (1994) and Wijmenga & van Buuren (1998).

References:

Heteronuclear multidimensional NMR experiments for the structure determination of proteins in solution employing pulsed field gradients. M. Sattler, J. Schleucher & C. Griesinger. Prog. in NMR Spectroscopy 34 (1999) 93-158.

Determining the structures of large proteins and protein complexes by NMR. G. M. Clore & A. M. Gronenborn. TIBTECH 16 (1998) 23-24.

TROSY-type triple resonance experiments for sequential NMR assignments of large proteins. M. Salzmänn, G. Wider, K. Pervushin, H. Senn & K. Wüthrich. J. Am. Chem. Soc. 121 (1999) 844-848.

Heteronuclear techniques in NMR studies of RNA and DNA. T. Dieckmann & J. Feigon. Curr. Op. Struct. Biol. 4 (1994) 745-749.

The use of NMR methods for conformational studies of nucleic acids. S. Wijmenga & B. N. M. van Buuren. Progr. NMR Spectr. 32 (1998) 287-387.

Biomolecular sample preparation

1.6

Following points are worthy of consideration when preparing an NMR sample of a protein or a nucleic acid fragment.

- ◆ **Molecular weight** is decisive for the degree of isotope labeling and choice of NMR experiments as indicated in Table 1.
- ◆ **Volume** should be 0.5 ml with 5% D₂O. For the special Shigemi tubes 0.2 ml is sufficient.
- ◆ **Sample concentration** should be around 1 mmol for a complete structure determination. Certainly the higher the better.
- ◆ **Solvent** is preferably water. If the sample is not stable in pure water try to dissolve it in a salt buffer, for instance, 25 mmol PO₄²⁻ and 25 mmol NaCl. The buffer should not contain hydrogen atoms. The total ionic strength of the buffer should not be too high, as this leads to longer pulses.
- ◆ **pH** should be the lowest possible that still does not destabilize the native fold. Observation of exchanging amide protons is possible only up to pH ~7.5 when using optimal techniques for water suppression.
- ◆ **Temperature** should be the highest possible which still does not destabilize the native fold. Due to increased mobility at higher temperatures the NMR-signals are narrowed, which is favourable.
- ◆ **Impurities** that contain hydrogen atoms, even at concentrations of 5-10% of the macromolecular sample, will complicate the spectral analysis, particularly if they are peptides or of low molecular weight.
- ◆ **Growth of fungi** is prevented by adding a minor amount of NaN₃.
- ◆ **If reduced cysteins** are present in the polypeptide sequence, add dithiotreitol (DTT) in order to prevent aggregation through disulfide bridging. Since DTT contains protons, its concentration should maximally be similar to the protein concentration.

Fig. 6. NMR-structure of a 120-residue protein, determined by using the experimental strategy for ¹⁵N, ¹³C-labeled proteins outlined in the current manual.



Pulse Calibration


2

Proper optimization of the pulse lengths is essential in biomolecular NMR where the sample concentrations often are minimal. Besides, in the pulse sequences with many RF pulses the sensitivity losses due to pulse imperfections accumulate.

RF routing

2.1

Fig. 7. RF-routing.



Enter the **edsp** display and select the RF-routing for triple resonance in the following way (note that the ^2H through the 4th channel is optional). Tune and match each channel, moving on from the lowest to the highest frequency.

If the console is equipped with the 2H-TX board: for the above 4th channel routing the connection to the 2H-TX board should be recabled from output-3 on the *first* router board to output-3 on the *second* router board. Change also 2H-TX address in the bsmstool from 3 to 7. This is done by typing in a UNIX window „bsms“, then selecting options 6, 8 and finally changing the address to 7. Select also „save“ and do „cf“ in XWINNMR. With such cabling deuterium observe is possible only with special pulse programs (which use o1=o4).

¹H observe pulse

2.2

The ¹H pulse calibrations should be done on the actual biomolecular sample because they strongly depend on the ionic strength of the solution. If the concentration is sufficient, use the sequence **zgpr** and look at the protein signals. If the protein concentration is too low for accurate pulse calibrations, use pulse sequence **zg** monitoring the water signal. For detailed instructions see the manual for Bruker Avance 1D/2D course.

¹³C inverse mode

2.3

The ¹³C pulse calibrations should be done in the inverse mode (¹³C through the 2nd channel) using the pulse sequence **dec90**. For instance, a sample of 0.1 M CH₃OH-¹³C in DMSO-d₆ can be used. The ¹³C-offset for the methyl group in methanol is 49 ppm, the methyl ¹H-offset is 3.28 ppm and the 1-bond coupling constant is 142 Hz. For detailed instructions see the manual for Bruker Avance 1D/2D course.

¹⁵N inverse mode

2.4

The ¹⁵N pulse calibrations should be done in the inverse mode (¹⁵N through the 3rd channel) using the pulse sequence **dec90f3**. For instance, a sample of 0.1 M ¹⁵N-urea in DMSO-d₆ can be used. The ¹⁵N-offset for the amide groups in urea is 76 ppm, the amide ¹H-offset is 5.4 ppm and the 1-bond coupling constant is 89 Hz. The procedure is analogous to that for ¹³C in 2.3.

¹³C selective excitation/inversion, inverse mode

2.5

The ¹³C pulse calibrations should be done in the inverse mode (¹³C through the 2nd channel) using the pulse sequence **dec90sp** and **dec180sp**. For instance, a sample of 0.1 M ¹³C-CH₃OH in DMSO-d₆ can be used. You also calculate the values by using the „shape tool“ in XWINNMR, **stdisp**.


Selective excitation

The procedure is analogous to the pulse calibration in Chapter 2.3, except that for ^{13}C -excitation a shaped pulse instead of a hard pulse is applied. The pulse sequence starts with an initial proton magnetization that is allowed to dephase under the 1-bond scalar coupling to ^{13}C during $(2 \cdot J_{\text{CH}})^{-1}$. An antiphase signal is obtained, the signal intensity of which depends on the flip angle of the selective ^{13}C -pulse. When the ^{13}C -pulse is a perfect 90° -pulse, the magnetization is completely converted to a multiple quantum coherence and the observed signal reaches its *minimum*.

Table 2: Acquisition parameters for selective ^{13}C 90° pulse calibration.

Parameter	Value	Comments
PULPROG	dec90sp	pulse program
NUC1	1H	nucleus on f1 channel
O1P	3.5 ppm	^1H offset
NUC2	13C	nucleus on f2 channel
O2P	49 ppm	^{13}C offset
NS	1	
DS	0	
CNST2	139	$^1J_{\text{CH}}$
PL2	120 dB	no ^{13}C hard pulses
P13		selective pulse length (409 μs , e.g.)
SP2		power level for selective pulse
SPNAM2	G4.256	name of selective pulse
SPOAL2	0.5	phase alignment of selective pulse
SPOFF2	0.0	offset of selective pulse

Fig 8. Shaped 90° - ^{13}C -pulse calibration (anti-phase signal).



Selective inversion


The pulse sequence starts with an initial proton magnetization that is allowed to dephase under the 1-bond scalar coupling to ^{13}C during $(2 \cdot J_{\text{CH}})^{-1}$. Subsequently a hard 90° -pulse on the ^{13}C -nuclei creates a multiple quantum coherence. This coherence is refocused by the subsequent 180° - ^1H -pulse and a refocussing delay of $(2 \cdot J_{\text{CH}})^{-1}$, yielding an in-phase signal.

If the selective ^{13}C -pulse, applied after the 180° - ^1H -pulse, is a perfect 180° -pulse, the magnetization stays as multiple quantum coherence and the observable signal is *minimal*.

Table 3: Acquisition parameters for selective ^{13}C 180° pulse calibration.

Parameter	Value	Comments
PULPROG	dec180sp	pulse program
NUC1	1H	nucleus on f1 channel
O1P	3.5 ppm	^1H offset
NUC2	13C	nucleus on f2 channel
O2P	49 ppm	^{13}C offset
NS	2	
DS	2	
CNST2	139	$^1J_{\text{CH}}$
PL2	120 dB	no ^{13}C hard pulses
P13		selective pulse length (256 μs , e.g.)
SP3		dB power level for selective pulse
SPNAM3	Q3.256	name of selective pulse
SPOAL3	0.5	phase alignment of selective pulse
SPOFF3	0.0	offset of selective pulse

Fig. 9. Shaped 180° - ^{13}C -pulse calibration (in-phase signal).



²H decoupling pulse

2.6

The use of partially and/or fully ²H-labeled samples enables NMR investigations of larger biomolecules because the major cause of ¹³C line broadening and signal losses is eliminated, that is, the major source of the carbon relaxation (the dipolar coupling to the attached protons) is removed.

The calibration of the deuterium decoupling pulse is done on the ASTM sample consisting of 60% benzene-d₆ in p-dioxane. The following instructions refer to spectrometers equipped with a so called 2H-TX board in the console.


Since almost every probe is equipped with a deuterium lock channel, the only additional hardware required for deuterium decoupling is a free RF-channel and a switch, which allows alternating the lock channel between locking and decoupling during the experiment. This switch is installed in the 2H-TX board. The board also contains a 20.0 W deuterium amplifier.

First do the ²H-recabling as described in the section 2.1 „RF-routing“. It is necessary to change the RF-routing to enable ¹³C observation. In the **edsp** display select 13C through F1 channel. Deuterium should be routed from the F4 channel through the 300.0 W amplifier, which is now replaced by the 20.0 W 2H-TX-amplifier because of the recabling. Save the RF-routing and exit the display. Now read in a pulse sequence called **decp902hf4**. This is the same experiment as used in Section 2.3 for inverse pulse calibration, but modified to activate the 2H switch.

Table 4: Acquisition parameters for deuterium pulse calibration.

Parameter	Value	Comments
PULPROG	decp902hf4	pulse program
NUC1	13C	nucleus on f1 channel
O1P	128 ppm	¹³ C offset
NUC4	2H	nucleus on f4 channel
O4P	7.28 ppm	² H offset
NS	1	
DS	0	
CNST5	25	¹ J _{CD}

Fig. 10. Deuterium pulse calibration.



Since deuterium is a spin 1 nucleus the resulting spectrum contains a triplet instead of the doublet of a spin $\frac{1}{2}$ nucleus. Fig. 10 shows the three cases from the calibration. At the 90° -pulse the outer lines are at minimum and the middle line reaches its maximum.

Delays

2.7

In most 3D experiments, J-couplings (see App. 13.1) are used to transfer magnetization from one nucleus to another one. The corresponding delays have been preset in the standard pulse programs to the optimized values reported in the original literature. You might want to change a value, for example, due to fast T_2 relaxation. This can only be done by modifying the value in the pulse program. We recommend you to start with the set values.

Table 5: Some predefined delays

Coupling const.	delay	length /ms	experiment	remarks
$(3 \cdot J_{CH})^{-1}$	d3	2.2		
$(4 \cdot J_{CH})^{-1}$	d4	1.6-1.8		
$(6 \cdot J_{CH})^{-1}$	d21	1.1		
$(2 \cdot J_{NH})^{-1}$	d21	5.5		
$(4 \cdot J_{C^{\alpha}C^{\beta}})^{-1}$	d22	4.0		
$(2 \cdot J_{C^{\alpha}C^{\beta}})^{-1}$	d22	3.6		2-bond
$(2 \cdot J_{C^{\alpha}C^{\beta}})^{-1}$	d22	4.4		1-bond
$(4 \cdot J_{NC})^{-1}$	d23	12.0		
$(4 \cdot J_{CC})^{-1}$	d24	3.6		
$(4 \cdot J_{NH})^{-1}$	d26	2.3		
T(N)	d21	12.4		
T(C)	d23	12.0		$2T=24ms \leq 1/4J_{CH}$

Summary of the calibrated pulse lengths**2.8**Table 6: ^1H -field: MHz Date:**F1-channel / ^1H**

nucl./shape	length/ms	power/dB	comment
^1H	p1	pl1	90° hard pulse
	p6	pl10	90° mixing pulse (~25 μs)
			90° DIPSI decoupling (~40 μs)
	pcpd1	pl19	90° decoupling & ROESY (~100 μs)
Sinc1.1000	p11	sp1	90° water flip-back (~2ms)

F2-channel / ^{13}C

nucl./shape	length/ms	power/dB	comment
^{13}C	p3	pl2	90° hard pulse
	p9	pl15	90° mixing pulse (~25 μs)
	pcpd2	pl12	90° decoupling (~50 μs)
G4.256	p13	sp2	90° selective on-res.(409 μs) (60ppm)
G4.256	p13	sp4/sp6	off-resonance
G4tr.256	p13	sp8	on-res., time reversed
Q3.256	p14	sp3	180° selective on-res.(256 μs)(60 ppm)
Q3.256	p14	sp5/sp7	off-resonance
Q3.256			180° selective (373 μs) (ca. 40 ppm)
Q3.256	p24	sp9	180° selective (768 μs) (ca. 24 ppm)

F3-channel / ^{15}N

nucl./shape	length/ms	power/dB	comment
^{15}N	p21	pl3	90° hard pulse (~35 μs)
	pcpd3	pl16	90° decoupling (~200 μs)

F4-channel / ^2H

nucl./shape	length/ms	power/dB	comment
^2H	pcpd4	pl17	90° decoupling (200-300 μs)


In the above Table the pulse lengths and power levels are denoted according to the Bruker nomenclature. The standard pulse programs are written in accordance with these conventions. A complete overview of the nomenclature is found on your spectrometer in the file /u/exp/stan/nmr/lists/pp/Param.info.

Preparatory 1D and 2D experiments

2.9

Having inserted your biomolecular sample, set the appropriate temperature, lock on H₂O, tune and match the probe head and perform gradient shimming. Thereafter, set-up a 1D-experiment using the **zgpr** pulse sequence. Optimize the ¹H-offset in the „gs“ mode for the maximal water suppression and determine the ¹H pulse lengths.

Fig.11. Coupled ¹H 1D-spectrum of a ¹⁵N, ¹³C-labeled protein.



REGIONS: *amide* *aromatic* *alpha* *beta* *other-aliphatic* *methyl*

Now record the 2-dimensional ¹³C-HSQC and ¹⁵N-HSQC experiments in order to check your sample and the spectral quality. For the ¹³C-HSQC use the pulse sequence **invietgpsi** and set the coupling constant to 135 Hz. For the ¹⁵N-HSQC use the pulse sequence **invietgpsif3** and set the coupling constant to 92 Hz. Detailed instructions to these two experiments are given in the manuals for the Bruker Avance 1D/2D and GRASP training courses.

The aliphatic- ^{13}C spectral width extends from 0 to 75 ppm and the aromatic- ^{13}C frequencies are located between 105 and 135 ppm. In the order to increase the resolution in the carbon-13 dimension without increasing the number of points - and thus the experiment time - the spectra can be folded (if MC2=TPPI) or aliased (if MC2=States or States-TPPI), see Fig. 12.

Fig. 12 . Folding and aliasing of the ^{13}C -dimension.





Fig. 13. Example of aliasing of the ^{13}C -dimension in the ^{13}C -HSQC spectrum.



In Fig. 13 the ^{13}C -spectral window has been reduced to one fourth and the position of the ^{13}C -offset has been chosen such as to minimize signal overlap in the resulting spectrum. The alpha and methyl carbon frequencies outside the spectral width become aliased. For additional information about ^{13}C -frequencies in polypeptides see Appendix 13.2.

^1H - ^{13}C TOCSY-HSQC


3

Introduction

3.1

Overlap in two-dimensional experiments such as NOESY and TOCSY is readily resolved by spreading the signals to a third dimension according to a heteronuclear frequency. The ^{13}C -edited TOCSY experiment consists of a ^1H - ^1H TOCSY and gradient selected ^1H - ^{13}C HSQC part. The information content of the experiment corresponds to the 2D TOCSY experiment, that is, magnetisation transfer takes place between scalar coupled protons within a so called spin system. The ^1H -signals are spread out to the additional dimension according to the ^{13}C chemical shift in order to alleviate overlap.

The flow of the magnetization is as follows



In a protein each residue forms a separate spin system because magnetization is not transferred over the backbone carbonyl group in this experiment. Hence the ^1H frequencies and the ^{13}C frequencies of their attached carbons can be identified in the ^{13}C -correlated TOCSY experiment.


References:

Experiments for recording pure-absorption heteronuclear correlation spectra using pulsed field gradients. A. L. Davis, J. Keeler, E. D. Laue, D. Moskau. J. Magn. Res. 98 (1992) 207-216.

Protein NMR Spectroscopy. J. Cavanagh, W. J. Fairbrother, A. G. Palmer III, N. J. Skelton. Academic Press Inc. (1996).

Pulse Sequence Diagram

3.2

Fig. 14. ^1H - ^{13}C TOCSY-HSQC.

In the above diagram 90° -pulses are denoted by thin bars and 180° -pulses are denoted by thick bars. The pulse phases are x if not specified. On the bottom of the Figure the TOCSY and HSQC parts of the pulse sequence are indicated, as well as the pulses and delays that achieve the heteronuclear gradient echo (echo-antiecho processing, E/A).

The trim pulses before and after the TOCSY mixing suppress solvent and artifactual coherences perpendicular to the coherence of interest. The amplitude of the first gradient should add up to 3.976 times the amplitude of the last gradient in order to select the ^{13}C -coherences (as the ratio of the magnetogyric ratios $\gamma(^1\text{H})/\gamma(^{13}\text{C})$ is 3.976). For the echo-antiecho (E/A) coherence selection the amplitude of the first gradient is inverted together with the phases of receiver and the two ^{13}C -pulses preceding t_2 . For each t_1 -increment the phase ϕ_1 is incremented.

There is a trade off between the number of points in the indirect dimensions and the number of transients that can be recorded within a reasonable time. Particularly in the indirect ^1H -dimension the digital resolution needs to be sufficient to resolve cross signals. Consequently the phase cycle has to be kept short. Here a minimal phase cycle is employed, two steps for the isotope filter and two steps for the axial peak suppression in t_1 and t_2 .

Setting up the experiment

Sample: 50 mm unlabeled cyclosporin (12 a.a.) in benzene- d_6 .

Experiment time: 16.5 hours.

First record a 1D-spectrum to determine the required ^1H sweep width and to optimize the ^1H set-up. Type **iepxno** to create and enter a new experiment. Invoke the **edsp** display in order to set up the RF-channel routing. The current pulse sequence uses the second channel for ^{13}C -excitation. Depress the off-button underneath the F2-button and select ^{13}C and simply click on the **default setting** button. Click on **save** upon leaving the display. Now change to 3D parameter mode by typing **parmode** and selecting **3D**. The program will ask you whether you want to delete the existing meta etc. files, answer yes. Enter the following values by using **eda** display for defining the three dimensions and the **ased** command which invokes only the parameters that are active in the particular experiment. You need to edit the pulse program **mleviief3gp3d** to use the second channel instead of the third (change all *f3* to *f2* and *cpd3* to *cpd2*).

Table 7: Acquisition Parameters

Parameter	Value	Comments
PULPROG	mleviief2gp3d	pulse program
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	^1H offset
NUC2	13C	nucleus on f2 channel
O2P	37 ppm	^{13}C offset
PL1		high power level f1 channel (^1H)
PL10		power level for TOCSY mixing (^1H)
PL2		high power level f2 channel (^{13}C)
PL12		decoupling power level f2 channel (^{13}C)
PI16	120 dB	for safety only
P1		90° ^1H pulse (f1 channel)
P2	preset to 2*p1	180° ^1H pulse (f1 channel)
P5	preset 0.667*p1	60° pulse for TOCSY mixing (f1)
P6		90° pulse for TOCSY mixing (f1)
P7	preset to 2*p6	180° pulse for TOCSY mixing (f1)
P17	2.5 ms	low power trim pulse (f1)
P28	1u	high power trim pulse (not applied!)
P3		90° ^{13}C pulse (f2 channel)

^1H - ^{13}C TOCSY-HSQC

Parameter	Value	Comments
P4	preset to 2*p3	180° ^{13}C pulse (f2 channel)
D1	1.5 s	recycle delay
D9	70 ms	TOCSY mixing time
SL_time	calculated	actual mixing time, for information only
CNST4	135	J_{CH} coupling constant, $D26=1\text{s}/(\text{CNST4}^4)$
P16	1000 us	gradient length
D16	50 us	gradient recovery
GPNAM1	SINE.100	#1 gradient name
GPNAM2	SINE.100	#2 gradient name
GPZ1	80	#1 gradient amplitude
GPZ2	20.1	#2 gradient amplitude
NS	4	number of scans
DS	128	number of dummy scans (multiple of ns)
CPDPRG2	garp	decoupling scheme f2 channel (^{13}C)
PCPD2		decoupler pulse length f2 channel (^{13}C)
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	12 ppm	sweep width
F2 indirect ^{13}C	*****	(middle column)
TD	64	number of real points
SW	75 ppm	sweep width indirect ^{13}C
ND10	2	no of in10 in pulse program
F1 indirect ^1H	*****	(rightmost column)
TD	128	number of real points
SW	12 ppm	sweep width ^1H
ND0	4	2x no of in0 in pulse program for TPPI

Note that the delays d0 and d10 are set automatically and that the increments are calculated through $\text{in0}=[\text{nd0}*\text{swh}(\text{F1})]^{-1}$ and $\text{in10}=[\text{nd10}*\text{swh}(\text{F2})]^{-1}$. The acquisition time AQ is set through $\text{AQ}=\text{TD}/\text{SW}^2$. Therefore, these parameters are not reported above.

To inspect and edit the pulse program type **edcpul** on the command line. To view the pulse program as the pulse sequence diagram, type **ppg** on the command line. Before starting an experiment, you may check the receiver gain automatically with the **rga** command and the total measurement time with the **expt** command. High receiver gain is necessary in this case.

The individual parameters for each dimension can be called for on the command line by typing, for instance,

td (reports the number of points set in F3)

2 td (reports the number of points set in F2)

1 td (reports the number of points set in F1).

While the experiment is running, the actual status of it can be monitored by the commands

2s td (reports the number of points recorded in F2)

1s td (reports the number of points recorded in F1).

Recording the two 2D planes

3.4

In order to test the experimental set-up and the quality of the spectrum, it is highly recommended to initially record the corresponding two 2D experiments, that is, the ^1H , ^1H -NOESY with only ^{13}C -bound protons selected in the acquisition dimension („13“) and the ^1H , ^{13}C -HSQC plane („23“). In XWIN-NMR create and enter a new experiment with the command **ixpno**. The individual planes can be recorded by using the 3D parameters and setting the number of points collected, **td**, in the inactive dimension to 2.

Inspection of 2D planes

First enter the **edp** display and insert the appropriate processing parameters from Table 6 below. To inspect the 2D planes in a 3D experiment, you can either click on the button **CREATE 2D** in the 3D data set, or just simply type **xfb** on the command line. The program will ask you the following questions.

Select direction? 13 or 23. In the present experiment **13** corresponds to the ^1H - ^1H NOESY and **23** corresponds to the ^1H - ^{13}C HSQC.

Enter slice number? You can inspect the different planes one by one, start with the first one.

Enter 2D procno? - You can choose any number except, of course, processing number 1 where the 3D is stored.

The program will now change to a 2D display and perform the processing using appropriate parameters from the 3D. Phase correct the 2D spectrum if necessary.

Phase corrections

Record the two 2D experiments and process them with the appropriate processing parameters given in Table 6. Phase the spectra and note the correction **phc0** and **phc1** values needed for the indirect dimension in each spectrum. Save the phase corrections for the 3D spectrum by returning, **rep 1**, and entering them in **edp**.

Observe that in the echo-anti-echo type of experiments the required zeroth order phase correction (phc0) in the acquisition dimension of the HSQC-plane differs -90° from the value required for the correct phasing of the ^1H - ^1H -plane. The appropriate value for the F3 dimension in the 3D processing is the value determined for the ^1H - ^{13}C HSQC plane.

In the indirect dimensions phase any deviations are due to phase evolution during the pulses within and flanking the particular evolution time. These effects can be compensated for by manipulating the first increment in the pulse sequence. For instance, the postacquisition phase corrections are phco=90 and phc1=-180 if the pulse sequence contains the statement: „d0=(in0-p22)/2-p1*2/3.14“. Depending on the length of in0 this statement may, however, accrue a negative value. In that case, it should be modified to „d0=in0-p22/2-p1*2/3.14“ giving phco=180 and phc1=-360.

Window functions

Optimize the window functions interactively on 1D slices from the 2D planes by using the „**manual window adjust**“ option available in the **process** menu. which shows the window function and its impact on the fid and the spectrum simultaneously. First type **rser 1** to display and manipulate the first fid. Process it with **ef** and phase the spectrum. Now activate the „**manual window adjust**“ tool and optimize the window function of your choice. Upon returning save it with the option „**save as 2D & return**“. Now perform the fourier transform along the rows only with the command **xf2**. The interferograms are displayed. In the **utilities** menu click on the button „**column**“ and select a column which clearly contains some signal. You will move into 1D processing mode when you upon returning „**store the column as the ~TEMP experiment**“. Now repeat the previous procedure of manual window adjustment. As the number of points in the indirect dimensions of a 3D is limited, the best results are achieved by applying a window function which removes truncation artifacts, such as shifted sine bell, SINE, or squared sine bell, QSINE. The shifting is given as fractions of π , for instance, ssb=2 gives 90° shifted (= cosine window) and ssb=3 gives a 60° shifted sine bell. After having saved the optimized window, start the second fourier transform with the command **xf1** and inspect the result.

Spectral Processing

3.5

Processing the 3D spectrum

You can monitor the measurement while it is running or after its completion by simply typing the command **xfb** to monitor the 2D planes. The program will ask you, which plane you wish to process and under which processing to store it, as described above in Chapter 3.3.

Invoke the **edp** display and set the following parameters.

Table 8: Processing parameters


Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	2k	zero fill
WDW	QSINE	squared sine bell
SSB	3.3	70° shifted
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	no automatic baseline correction
F2 indirect ^{13}C	*****	(middle column)
SI	128 or more	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell
SSB	2	90° shifted
PH_mod	no	phase correction is not needed
ABSG	1	first order base line correction
BC_mod	no	no automatic baseline correction
F1 indirect ^1H	*****	(rightmost column)
SI	512	zero fill
MC2	TPPI	
WDW	QSINE	squared sine bell
SSB	2	90° shifted
PH_mod	pk	phase correction
PHC0		zero order phase correction
PH1		first order phase correction
ABSG	1	first order base line correction
BC_mod	no	no automatic baseline correction

The fourier transforms of the three dimensions are performed by typing the commands **tf3 no**, **tf2 no** and **tf1 no**, respectively (the „no“ here means that the imaginary part of the processed data will not be stored on the disc). Instead, you can write a macro by typing **edmac** and a new name, for instance, **fast3dproc**, which contains the three fourier transform commands, each on a separate line.

Upon completion of the processing the message „3D PROCESSED DATA AVAILABLE“ appears on the screen. Now you can inspect the spectrum by clicking on the **display** button. The program will ask you which countour levels you wish to display, 23, 13 or 12. Answer „yes“ „no“ „no“, and do not store the compressed data, in order to save disc space. It may be necessary

to reset the contour levels by **edlev**, an appropriate starting value is 10 000 000. You can inspect the projections of the spectrum by depressing the **project** and the **x90**, **y90** or **z90** buttons, or you can run the **movie** to monitor the complete cube. In order to inspect the individual planes, click on the button **scan**. Now you can either manually or automatically scan through the planes in the three orthogonal directions. To activate the display depress the appropriate button, **23**, **13** or **12**.

Fig. 15. The ^1H - ^1H and ^1H - ^{13}C planes of the ^1H - ^{13}C TOCSY-HSQC spectrum.



Additional processing parameters

3.6

A baseline correction with a polynomial can be applied to any of the three dimensions with the command **tabs1**, **tabs2** or **tabs3**. The order of the polynomial is given by the ABSG parameter for respective dimension. First order is recommended for simple removal of t_1 -noise. The range where the correction is applied is indicated by the ABSF1 and ABSF2 parameters. For instance, by setting these to 100 and 6 ppm the subsequent **tabs3** command will avoid the baseline correction in the vicinity and to the right of the water signal. In 2D spectra a baseline correction in the acquisition dimension which excludes the water signal within a user defined region, is invoked by the command **abs2.water**.

BC_mod when using digital acquisition should be set to „no“ „no“ „no“ when using digital quadrature detection (DQD). If post-acquisitional water deconvolution is needed, the baseline correction mode in F3 can be set to „**qpol**“ (a polynomial of 5th order is subtracted from the fid) or „**qfil**“ (filtering of the fid according to Marion, Ikura & Bax, JMR 84 (1989) 424-430). For the latter method the effective range of the filter should be set with the parameter BCFW. The parameter COROFFS applies to both methods and gives the correction offset from the center of the spectrum if the observation frequency o1 was not set on the water frequency.

PKNL should be set to **TRUE** for data collected in digital acquisition mode.

FT_mod is only activated for special transforms (performed using the „trf“ command).

If during the processing the command **tdeff** in any dimension is set to a value different from zero, only the specified number of points are used. Thus with the **tdeff** command points can be discarded from the end of the fid. To discard points from the beginning of the fid instead, the parameter **tdoff** defines the number of points with which the fid is left-shifted.

The spectrum is referenced by the command **sref**. The fine adjustment is done automatically with respect to the TMS signal as internal reference. If TMS is not used, type **sr 0**, **1 sr 0** and **2 sr 0** to make sure that the referencing is correct.

AQORDER defines in which order the fids were recorded in the serial data file. A quick look at the end of the pulse program will tell this. If the F1-dimension is incremented first, the inner loop in the pulse program contains the „id0“ command and the outer loop resets this parameter with the command „rd0“, then **AQORDER** is 3-1-2, otherwise it is 3-2-1. Note, however, that the **AQORDER** parameter is no longer essential for the correct processing of the spectrum, but it sets the referencing of the dimensions right.

Linear prediction

3.7

Since the number of points in the indirect dimensions of a 3D spectrum is limited, it is useful to try linear prediction of additional points. Select **ME_mod LPfc** (forward prediction complex). The parameter **NCOEF** represents the number of coefficients used in the calculation. Ideally this parameter should be set to 2-3 times the number of expected peaks. The number of points used is **td**, except if **TDEFF** > 0, in which case **tdeff** points are used for the prediction. The number of predicted points is $2 \cdot \text{SI} - \text{TD}$, thus it replaces the zero-filling. Linear prediction and zero filling can also be combined by setting the parameter **LPBIN** to a value between **TD** and $2 \cdot \text{SI}$. In that case the points from **TD** up to **LPBIN** will be predicted and the remaining points up to $2 \cdot \text{SI}$ are set to zero. In the ^1H - ^{15}N NOESY-HSQC experiment below, the number of real points in the ^{15}N -dimension is 40. This number is extended to 60 by setting **LPBIN** to 60. The size of the final complex data matrix **SI** is 64, thus the points 61-128 are set to zero.

^1H - ^{15}N NOESY-HSQC

4


Introduction

4.1

The ^{15}N -correlated NOESY experiment consists of a ^1H - ^1H NOESY and ^1H - ^{15}N HSQC part. The sensitivity enhanced and gradient selected version of HSQC is implemented. Although the pulse sequence is longer, the relaxation of ^{15}N -nuclei is not usually so fast as to lead to relaxation losses due to the additional delays. Gradient selection also effectively suppresses the solvent signal. The pulse sequence also contains two defocusing gradients while the coherences of interest are stored along the z-axis. Furthermore, the phases of the proton pulses are chosen such as to realign the water magnetization to the +z-axis. These measures are taken in order to avoid saturation transfer from water to the NH-groups which undergo fast exchange. Thus it becomes possible to observe exchanging NH-signals up to neutral pH.

For protein studies the heteronucleus of primary choice is ^{15}N while ^{15}N -labeling is usually feasible and also less costly than ^{13}C -labeling. In the ^{15}N -correlated NOESY spectrum of a protein the signals are dispersed according to the backbone amide frequency of each residue. Thus the NOE connectivities between a backbone amide group and the spacially neighbouring protons, either within the same or in other residues, are observed. In the ^{15}N -HSQC experiment maximal sensitivity enhancement can be expected since all backbone amides are of the same multiplicity, that is, singly protonated.


The flow of the magnetization is as follows



The ^{15}N -correlated NOESY experiment is usually combined with a ^{15}N -correlated TOCSY experiment. The latter yields the proton connectivities within each residue in a protein sequence. The sequential order of the residues can be traced by means of the NOE-connectivities, or through the triple resonance experiments if a ^{15}N , ^{13}C -labeled sample is available. It should be pointed out, however, that for larger molecules the ^1H - ^1H TOCSY transfer becomes less effective. The triple resonance experiment H(CCCO)NH actually confers identical information as ^{15}N -edited TOCSY, although the former relies on the transfer through ^{13}C - ^{13}C scalar couplings.

Pulse Sequence Diagram

4.2

Fig. 16. ^1H - ^{15}N NOESY-HSQC.

In the above diagram 90° -pulses are denoted by thin bars and 180° -pulses are denoted by thick bars. The pulse phases are x if not specified. The phases of the ^1H -pulses are chosen such that the H_2O -magnetization is realigned along the $+z$ -axis. At the bottom the NOESY and HSQC parts of the pulse sequence are indicated, as well as the pulses and delays that achieve the sensitivity enhancement using a heteronuclear gradient echo selection (echo-antiecho). The first and second gradients defocus solvent and artifact coherences in the x - y plane while the magnetization of interest is aligned along the z -axis. Thus the amplitudes of these two gradients are arbitrary and can be optimized. The amplitude of the third gradient should be 9.862 times the amplitude of the fourth gradient in order to select the desired coherences (the ratio of the magnetogyric ratios $\gamma(^1\text{H})/\gamma(^{15}\text{N})=9.862$).

For the sensitivity enhancement in every second HSQC-plane the amplitude of the third gradient is inverted together with the phase ϕ_6 . For each t_2 -increment the phases ϕ_3 , ϕ_4 and the receiver are inverted. For each t_1 -increment the phases ϕ_1 and ϕ_2 are incremented.

References:

Improved Accuracy of NMR Structures by a Modified NOESY-HSQC Experiment. W. Jahnke, M. Baur, G. Gemmecker, H. Kessler. JMR 106 B (1995) 86-88.

Data Acquisition

4.3

Setting up the experiment

Sample: 1mM ^{15}N - or ^{15}N , ^{13}C -labeled (requires ^{13}C -decoupling) polypeptide.

Experiment time: 2 days 20h

First record a 1D-spectrum to determine the required ^1H sweep width and to optimize the ^1H -offset (using the **gs**-tool, see below). Type **ixpno** to create and enter a new experiment. Invoke the **edsp** display in order to set up the RF-channel routing. The current pulse sequence uses the third channel for ^{15}N -excitation. Then change to 3D parameter mode by typing **parmode** and selecting 3D. The program will ask you whether you want to delete the existing meta etc. files, answer „yes“. Enter the parameters for the three dimensions in the **eda** and the rest in the **ased** display.

Table 9: Acquisition Parameters

Parameter	Value	Comments
pulse program	noesiif3gpsi3d	- adding ^{13}C -decoupling, see Section 5.4
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimized with „gs“)
NUC3	15N	nucleus on f3 channel
O3P	117 ppm	offset ^{15}N
PL1		high power level (f1 channel, ^1H)
PL3		high power level (f3 channel, ^{13}C)
PL16		decoupler power level (f3 channel, ^{13}C)
P1		90° ^1H pulse (f1 channel)
P2	automatically set	180° ^1H pulse (f1 channel)
P21		90° ^{15}N pulse (f3 channel)
P22	automatically set	180° ^{15}N pulse (f3 channel)
P28	1u	trim pulse, not used
D1	1.5 s	recycle delay
D8	100 ms	mixing time
CNST4	90	$d24=1\text{s}/(\text{CNST}4^4)$
P16	1000 us	gradient length
D16	50 us	gradient recovery
GPNAM1...4	SINE.100	gradient shape
GPZ1	30	#1 gradient amplitude

Parameter	Value	Comments
GPZ2	50	#2 gradient amplitude
GPZ3	80	#3 gradient amplitude
GPZ4	8.1	#4 gradient amplitude
NS	16	number of scans
DS	32	2*ns*x (because of echo-antiecho)
CPDPRG3	garp	decoupling scheme (^{15}N)
PCPD3		decoupler pulse length
F3 acquisition ^1H	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
F2 indirect ^{15}N	*****	(middle column)
TD	40	number of real points
SW	35 ppm	sweep width indirect ^{15}N
ND10	2	no of in10 in pulse program
F1 indirect ^1H	*****	(rightmost column)
TD	220	number of real points
SW	14 ppm	sweep width indirect ^1H
ND0	2	no of in0 in pulse program

The „GS“ Interactive mode to optimize acquisition parameters

The gs command sets the spectrometer in a mode where acquisition parameters can be adjusted interactively while the first scan of the experiment is repeated. The impact of the adjustment can be observed either on the fid or on the spectrum. Start with the command **gs** and enter the acquisition window by typing **acqu**. Click on the iconified gs-window to open it and modify the parameters. In the DISPLAY menu, invoked through a button on the top bar, you can alternate between „**display time domain**“ or „**display frequency domain**“. The additional option „**phasing**“ allows you to display the spectrum with applying known phase corrections, or as a magnitude spectrum. A condition for the frequency domain display to work properly is a previously processed spectrum. It is worth noting that you can also perform shimming in the gs-mode.

For the current experiment you should optimize the ^1H -offset, **o1**. You can also test changing the amplitude of the first two gradients, gpz1 and gpz2. The desired effect upon changing any of these parameters is the minimization of the solvent signal.


Table 10: Processing parameters

Parameter	Value	Comment
F3 acquisition ^1H	*****	(leftmost column)
SI	2k	zero fill
WDW	QSINE	squared sine bell
SSB	3	70° shifted squared sine bell
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	qpol	water deconvolution
BCFW	1	range for water deconvolution (in ppm)
ABSG	1	first order baseline correction
ABSF1	100 ppm	range for abs3 command
ABSF2	6 ppm	
STSR	0	- display only the left half of
STSI	1k	the spectrum (points 0-1024)
F2 indirect ^{15}N	*****	(middle column)
SI	128	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell
SSB	2	90° shifted squared sine bell
PH_mod	pk	phase correction (only if needed)
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	automatic baseline correction
ABSG	1	First order baseline correction
ABSF1	100 ppm	range for abs2 command
ABSF2	-100 ppm	
STSR	0	range for display
STSI	0	
ME_mod	LPfc	forward complex linear prediction
NCOEF	32	number of coefficients
LPBIN	64	extent of linear prediction
F1 indirect ^1H	*****	(rightmost column)
SI	512	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell

^1H - ^{15}N NOESY-HSQC

Parameter	Value	Comment
SSB	2	90° shifted squared sine bell
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	automatic baseline correction
ABSG	1	first order baseline correction
ABSF1	100 ppm	range for abs1 command
ABSF2	-100 ppm	
STSR	0	range for display
STSI	0	

Fig. 17. The ^1H - ^1H and ^1H - ^{15}N planes of the ^1H - ^{15}N NOESY-HSQC spectrum.




Introduction

5.1

The three-dimensional HNHA experiment is designed to accurately determine three-bond H^N-H^α J-coupling constants. First, problems with overlap are alleviated by spreading the signals to an additional dimension according to the ^{15}N -frequency. Second, the coupling constant is determined from the ratio between the intensities of the diagonal and cross-peak which are readily quantized. The magnetization transfer is of HMQC-type involving the zero- and double quantum coherences. Therefore the signals are not modulated by the one-bond HN-coupling. The ^{15}N -evolution occurs in a constant time fashion and it is incorporated into the periods during which also the magnetization transfer from the H^N to the H^α coherence takes place. This magnetization transfer is proportional to the length of the transfer period (denoted below by 2ϵ) and the size of the H^N-H^α scalar coupling. Finally the coupling constant is deduced from the relation

$$\text{Intensity}(\text{cross peak}) / \text{Intensity}(\text{diagonal peak}) = -\tan^2(2\pi\epsilon J_{HN\alpha})$$

Note that the diagonal and cross peaks in the resulting spectrum have opposite phase. The flow of the magnetization in the HNHA experiment is as follows



In proteins the three-bond H^N-H^α J coupling constants are an important source of information on the secondary structure and improve convergence and accuracy of the structure calculation particularly for α -helical fragments. Accurate determination of the H^N-H^α couplings is, however, complicated by their small size relative to the natural proton line width. A direct measurement is possible only for very small peptides, and also 2D methods combined with curve fitting to the fid become insufficient for proteins over 10 kDa.

References:

Quantitative J Correlation: A new approach for measuring homonuclear three-bond $J(H^N-H^\alpha)$ coupling constants in ^{15}N -enriched proteins. G. W. Vuister & A. Bax. JACS 115 (1993) 7772-7777.

Measurement of H^N-H^α J couplings in calcium-free calmodulin using new 2D and 3D water-flip-back methods. H. Kuboniwa, S. Grzesiek, F. Delaglio & A. Bax. J. Biomol. NMR 4 (1994) 871-878.

Protein NMR Spectroscopy. J. Cavanagh, W. J. Fairbrother, A. G. Palmer III, N. J. Skelton. Academic Press Inc. (1996).

Pulse Sequence Diagram

5.2




Fig. 18. HNHA.

$$\varepsilon = 13.5 \text{ ms}$$

$$\delta = 1/(J_{\text{HN}\alpha} * 2)$$

$$\phi_1 = 2x, 2(-x) + \text{TPPI}(t_1)$$

$$\phi_2 = x + \text{TPPI}(t_2) \quad \phi_3 = 4x, 4y$$

$$\phi_4 = x, y$$

$$\phi_{\text{rec}} = x, 2(-x), x, y, 2(-y), y$$

$$\gamma = \text{gradient} + \text{recovery}$$

$$\tau = 90 \text{ } \mu\text{s}$$

In the above diagram narrow and wide pulses denote 90° - and 180° -flip angles, respectively. The pulse phases are x if not specified. The first heteronuclear pulse creates zero- and double-quantum coherences in a HMQC type fashion. Simultaneous displacement of the first and the third heteronuclear 180° pulse during the subsequent delays, of total duration $4T+t_2$, causes chemical shift labeling of the ^{15}N -nucleus in a constant time manner. The gradient amplitudes should be $G_1+G_2=G_3$, for instance, $G_1=30$, $G_2=45$, $G_3=75$. The last gradient dephases residual magnetization in the x,y -plane while the coherence of interest is align along the z -axis. Its amplitude is arbitrary and should be optimized. The decrement $d30$ is automatically set equal to the increment $d10$ in the pulse program. To inspect the crossing of the last pair of 180° -pulses, the pulse program simulation tool **pulsdisp** found in the WINDOWS menu, can be used. The phases of the ^1H -pulses are chosen such that the H_2O -magnetization is realigned along the $+z$ -axis.

Solvent suppression is achieved by a flip-back method also known as the jump-and-return or 1-1 technique. For this purpose two additional 90° -pulses on-resonance for water, flanking a delay τ , are applied prior to the acquisition. During the delay τ the off-resonance frequencies precess in the x,y -plane. The last 90° -pulse brings water back to the z -axis, thus avoiding solvent saturation, whereas for the off-resonance frequencies the excitation profile becomes frequency dependent with a maximum at offset $\pm 1/4 * \tau$.

Sample: 1mM ^{15}N - or $^{15}\text{N},^{13}\text{C}$ -labeled (requires ^{13}C -decoupling) polypeptide.

Experiment time: 16.5h

Table 11: Acquisition Parameters

Parameter	Value	Comments
F3 acquisition	*****	(leftmost column)
PULPROG	hnha	or hnhagp.x with ^{13}C -decoupling (edit!)
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimized with „gs“)
NUC3	15N	nucleus on f3 channel
O3P	117 ppm	offset ^{15}N
PL1		high power level f1 channel
PL21	(=p1)	1-1 power level f1 channel
PL9	120	if presaturation is used, 70-90 dB
PL3		high power level f3 channel
P1		90° ^1H pulse (f1 channel)
P2	preset to 2*p1	180° ^1H pulse (f1 channel)
P0	(=p1)	first 90° 1-1 pulse (f1 channel)
P28	(=p1)	second 90° 1-1 pulse (f1 channel)
P29	1u	if presaturation is used, presat. pulse
P21		90° ^{15}N pulse (f3 channel)
P22	preset 2*p21	180° ^{15}N pulse (f3 channel)
D1	1.5s	recycle delay minus presaturation pulse
D19	90u	1-1 delay => excitation maximum at O1+/- 4 ppm
D2	preset 5.2 ms	$D2=1s/(CNST2*2) = 1/2 (J_{HN})$
D24	preset 13.5ms	$\leq J \text{ HNH}\alpha$
P16	600u	gradient length
D16	50u	gradient recovery
GPNAME1...4	SINE.50	gradient pulse name
GPZ1	30	#1 gradient amplitude
GPZ2	45	#2 gradient amplitude
GPZ3	75	#3 gradient amplitude
GPZ4	50	#4 gradient amplitude
NS	8	number of scans
DS	16	number of dummy scans
CPDPRG3	garp	decoupling scheme (^{15}N)
PCPD3		^{15}N decoupler pulse length (f3 channel)

Parameter	Value	Comments
PL16		¹⁵ N decoupler power level (f3 channel)
TD	1024	no of points
SW	14 ppm	sweep width
AQ_MOD	DQD	digital quadrature detection
F2 indirect ¹⁵N	*****	(middle column)
TD	40	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	4	no of in10 in pulse program
IN30	=IN10	set CT-decrement in30 equal to in10
F1 indirect ¹H	*****	(rightmost column)
TD	128	number of real points
SW	14 ppm	sweep width ¹ H
ND0	2	no of in0 in pulse program

„GS“ Interactive mode to optimize acquisition parameters

For the current experiment you should optimize the ¹H-offset, **o1**, and **phcor12** which changes the phase of the second jump and return 90°-pulse. This pulse phase should be optimized in very small steps, ~0.01 degrees at a time. The same result can be achieved by optimizing the length of the second ¹H-pulse in steps of 12.5 ns. These parameters are optimized to minimize the solvent signal.

Broadband, adiabatic and selective ¹³C decoupling 5.4

Broadband decoupling

If you are measuring on a ¹³C,¹⁵N-labeled sample you should implement ¹³C-decoupling during the ¹H-evolution and during the acquisition, in order to prevent the splitting of the alpha-proton signal due to coupling to the bound carbon. You can achieve this in two ways.

Insert a 180° ¹³C-pulse in the middle of the ¹H-evolution, centered with respect to the 180° ¹⁵N-pulse. In the beginning of the pulse program type „CEN_HC2=(p22-p4)/2“ and then

```
d0 pl2:f2
(p22 ph0):f3 (CEN_HC2 p4 ph0):f2
d0 pl12:f2
```

Now add composite pulse decoupling during acquisition

```
go=2 ph31 cpd2:f2 cpd3:f3
d1 do:f2 do:f3 wr *0 if *0 zd
```

Pay particular attention to set the power level and the decoupler-off commands (do:f2) correctly.

Thus, for broadband ^{13}C -decoupling the following parameters are needed:

CPDPRG2	garp	decoupling program
PCPD2		decoupling pulse length
PL12		decoupling power level
PL2		hard pulse power level
P4		180° hard pulse

Adiabatic decoupling

You might want to use the adiabatic ^{13}C -decoupling while it enables lower power, and thus, less heating. An appropriate power level is +2dB with respect to the power level determined for ^{13}C -GARP-decoupling. You can create the shaped pulse CHIRP95 optimized for minimal sidebands, with the „shape tool“, **stdisp**. Select the shape „smoothed chirp“, change the total sweep width to 60 000 MHz and the length of pulse to 1500us. Save it under the name CHIRP95. Set the following parameters:

CPDPRG2	p5m4sp180	decoupling program
PCPD2	1.5m	decoupling pulse length
PL12		decoupling power level
SP15	=pl12	note: + 2dB compared to garp
SPNAM15	CHIRP95	adiabatic shaped pulse

Selective decoupling

In some cases it is necessary to use selective ^{13}C -decoupling, for instance on the carbonyls or alpha-carbons (Section 10.3). Set the following parameters:

CPDPRG2	mlev180sp	decoupling program
PCPD2	768u	decoupling pulse length
PL12		decoupling power level
SP15	=pl12	
SPNAM15	Q3.256	selective 180° pulse
SPOFF15		^{13}C -offset for selective decoupling

Spectral Processing

5.5


Table 12: Processing parameters

Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	2k	zero fill
WDW	QSINE	squared sine bell

HNHA

Parameter	Value	Comment
SSB	3	70° shifted
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	qpol	no automatic baseline correction
STSR	0	- display only the left half of
STSI	1k	the spectrum (points 0-1024)
F2 indirect ¹⁵N		<i>(middle column)</i>
SI	128	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	2	90° shifted
PH_mod	no	phase correction
BC_mod	no	no automatic baseline correction
F1 indirect ¹H		<i>(rightmost column)</i>
SI	512	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	2	90° shifted
PH_mod	pk	phase correction
PHC0	180	zero order phase correction
PHC1	-360	first order phase correction
BC_mod	no	no automatic baseline correction

Fig. 19. The ¹H-¹H plane of HNHA. Nb. Crosspeaks are negative.







Fig. 20. Magnetization transfer in the HNCA experiment.

The HNCA experiment provides two types of correlations, namely the correlations from the backbone amide resonances H^N_i and N_i to the alpha-carbon C^α_i within the same residue, and to the alpha-carbon C^α_{i-1} of the preceding residue. This is because the two coupling constants, $^1J_{NiC\alpha i}$ (=11 Hz) and $^2J_{NiC\alpha i-1}$ (=7 Hz) are of similar size. Since only the latter, the sequential connectivity, is observed in the HN(CO)CA experiment, combining information from HNCA and HN(CO)CA yields unambiguous sequential correlations. The HNCA experiment consists of four coherence transfer steps. The flow of the magnetization is the following:



The magnetization transfer in the second and third steps, between ${}^{15}N$ and ${}^{13}C^\alpha$ can take place via a single quantum (INEPT-type) or multiple quantum (HMQC-type) coherence. The former, the single quantum experiment, is presented here. Its advantage is that the ${}^{15}N$ spin part does not undergo transverse relaxation since it stays longitudinal during the ${}^{13}C$ -evolution time.

References:

Improved 3D triple resonance NMR techniques applied to a 31 kDa protein. S. Grzesiek & A. Bax. J. Magn. Res. 96 (1992) 432-440.

Minimization of sensitivity losses due to the use of gradient pulses in triple-resonance NMR of proteins. J. Stonehouse, R. T. Clowes, G. L. Shaw, J. Keeler & E. D. Laue. J. Biomol. NMR 5 (1995) 226-232.

Pulse Sequence Diagram

6.2




Fig. 21. HNCA.

$$\phi_2 = X, X, -X, -X$$

$$\phi_1 = X, -X + \text{States-TPPI } (t_1)$$

$$\phi_3 = Y, Y, -Y, -Y + \text{echo-antiecho } (t_2)$$

$$\phi_{\text{rec}} = X, -X, -X, X$$

$$\delta = 2.3 \text{ ms}$$

$$\epsilon = \text{gradient} + \text{recovery}$$

$$\gamma = 5.5 \text{ ms}$$

$$T = 12.0 \text{ ms}$$

The sequence starts with an INEPT transfer of the magnetization from the backbone ^1H to ^{15}N . A dephasing delay T follows, during which the 1-bond coupling between ^{15}N and $^{13}\text{C}^\alpha$ leads to an antiphase three-spin coherence $\text{H}_z\text{N}_x\text{C}_z^\alpha$. This is converted to $\text{H}_z\text{N}_z\text{C}_y^\alpha$ by the first selective ^{13}C 90° -pulse and labeled by the $^{13}\text{C}^\alpha$ chemical shift during the t_1 -evolution time. It is converted back to the $\text{H}_z\text{N}_y\text{C}_z^\alpha$ antiphase coherence by the time reversed ^{13}C 90° -pulse. This pulse is time reversed in order to accomplish a pure 90° -rotation in the reverse direction, that is, from the transverse plane back to the z -axis.

The three-spin coherence rephases to an ^{15}N -coherence during the delay T and becomes simultaneously labeled by the ^{15}N -chemical shift during the constant time-type of t_2 -evolution. Overlaying the rephasing and evolution delays minimizes transverse relaxation losses. Use of constant time evolution removes the relaxation decay of the magnetization of interest. So called „mirror image“ linear prediction can be used as the signal in the t_2 -dimension becomes a sum of undamped cosinusoidal oscillations of known phase ($=0^\circ$). A reverse INEPT transfer with sensitivity enhancement converts the coherences of interest back to the amide protons for observation.

Protons are decoupled by the DIPSII-2 scheme. Use of broadband decoupling prevents transverse ^{15}N -relaxation, giving further gain in sensitivity compared to the use of 180° -pulses. For reproducible water suppression it is essential to perform the decoupling in a synchronous mode. The ^{13}C -labeled carbonyl carbons are decoupled by a selective 180° -pulse during the $^{13}\text{C}^\alpha$ -evolution.

Sample: 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).

Experiment time: 4h

Table 13: Acquisition Parameters

Parameter	Value	Comments
PULPROG	hncagp3d.2	pulse program
NUC1	^1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimize with „gs“)
NUC2	^{13}C	nucleus on f2 channel
O2P	52 ppm	offset ^{13}C (centered on C^α region)
NUC3	^{15}N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse f2
PL1		high power level f1 channel
PL2	120 dB	high power level f2 channel (not used)
PL3		high power level f3 channel
PL16		power level for GARP decoupling (f3)
PL19		power level for DIPSI-2 decoupling (f1)
SP1		power level, selective H_2O flip-back, optimize in „gs“-mode!
SP2		power level, 90° shaped pulse (C^α on res.) f2
SP3		power level, 180° shaped pulse (C^α on res.) f2
SP5	=sp3	power level, 180° shaped pulse (C' off res.) f2
SP8	=sp2	power level, time reversed 90° (C^α on res.) f2
SPNAM1	Sinc1.1000	selective H_2O flip-back
SPNAM2	G4.256	90° shape: 4 Gaussian cascade
SPNAM3	Q3.256	180° shape: 3 Gaussian cascade
SPNAM5	Q3.256	180° shape: 3 Gaussian cascade
SPNAM8	G4tr.256	90° , time reversed 4 Gaussian cascade
SPOFF1	0	selective H_2O flip-back
SPOFF2	0	
SPOFF3	0	
SPOFF5	18270 (600Mz)	≤ 121 ppm (^{13}C) in Hz (=173-52 ppm)
SPOFF8	0	
P1		90° ^1H pulse f1 channel
P2	preset to 2*p1	180° ^1H pulse f1 channel
P3		90° ^{13}C pulse f2 channel
P4	preset to 2*p3	180° ^{13}C pulse f2 channel

HNCA

Parameter	Value	Comments
P11	2000u	selective H ₂ O flip-back
P13	409u	90° selective pulse f2 (333u at 700MHz)
P14	256u	180° selective pulse f2 channel
P16	1000u	gradient pulse length
P26	=pcpd1	90° DIPSI-2 pulse at pl19, proton
p21		90° ¹⁵ N pulse f3 channel
p22		180° ¹⁵ N pulse f3 channel
PCDP1		decoupling pulse DIPSI-2 channel f1
CPDPRG1	dipsi2	decoupling program channel f1
PCDP3		decoupling pulse GARP channel f3
CPDPRG3	garp	decoupling program channel f3
D1	1.0s	recycle delay
D16	50u	gradient recovery
GPZ1	30	#1 gradient amplitude
GPZ2	80	#2 gradient amplitude
GPZ3	8.1	#3 gradient amplitude
GNAM1	SINE.100	#1 gradient shape
GNAM2	SINE.100	#2 gradient shape
GNAM3	SINE.100	#3 gradient shape
NS	4	number of scans
DS	16	dummy scans (x*2*ns because of E/A)
RG	2k or more	receiver gain
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
NUCLEI	1H	
F2 indirect ¹⁵N	*****	(middle column)
TD	40	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	2	no of in10 in pulse program
NUCLEI	15N	
IN30	= IN10	Constant time: decrement equal to increment
F1 indirect ¹³C	*****	(rightmost column)
TD	64	number of real points
SW	40 ppm	sweep width indirect ¹³ C
ND0	2	no of in0 in pulse program
NUCLEI	13C	

Adding ^2H decoupling

Deuteron decoupling in the aliphatic region is applied during the ^{13}C -evolution time t_1 when the magnetization is on the alpha carbons. The broadband ^1H -decoupling is interrupted during this time period. The proton-carbon couplings are refocused by a 180° ^1H -pulse, instead.

Table 14: ^2H decoupling parameters

Parameter	Value	Comments
PULPROG	hncagp2h3d.2	pulse program
NUC4	^2H	nucleus on f4 channel
O4P	3.0 ppm	offset ^2H (centered on aliphatic protons)
pl4	120 dB	power level for hard pulses f2 channel
pl17		power level for deuterium decoupling
cpdprg4	waltz16	decoupling program
pcpd4		decoupling pulse length ($\sim 300\mu\text{s}$)

Spectral Processing

6.4


Table 15: Processing parameters

Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	2k	zero fill to 2048 complex points
WDW	QSINE	squared sine bell window function, e.g.
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction applied
BC_mod	qpol	water deconvolution
BCFW	0.5-1.0	range for water deconvolution
F2 indirect ^{15}N	*****	(middle column)
SI	128	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell window function
SSB	2	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction
F1 indirect ^{13}C	*****	(rightmost column)

HNCA

Parameter	Value	Comment
SI	256	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell window function
SSB	2	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction

Fig. 22. The ^1H - ^{13}C and ^1H - ^{15}N planes of the HNCA spectrum.



HN(CO)CA

7

Introduction

7.1





Fig. 23. Magnetization transfer in HN(CO)CA.

The HN(CO)CA experiment correlates the backbone amide proton and nitrogen (H^N_i and N_i) to the alpha-carbon of the preceding residue C^{α}_{i-1} . The magnetization transfer between the amide nitrogen and the alpha carbon takes place in two steps, via the intervening carbonyl carbon. The coupling constants are $J_{NC'} = 15$ Hz and $J_{C'C\alpha} = 55$ Hz. Particularly the latter is large, entailing shorter delays for the transfer at the time when the magnetization is on the fast relaxing heteronuclei. Thus the experiment can be used even at large line widths. Note that the second peak in HNCA originates from the direct transfer between the same nuclei due to the weaker two-bond coupling. HN(CO)CA consists of in total six coherence transfer steps:



Note that in both the HNCA and HN(CO)CA experiments the ${}^{13}\text{C}$ -evolution time $t_{1,\text{max}}$ should be kept shorter than $1/(2J_{\text{C}\alpha\text{C}\beta})$, in practice at around 10 ms, in order not to resolve this J-splitting of the signals in the carbon-13 dimension. Similarly, in order to avoid resolving the one-bond ${}^1J_{\text{NC}'}$ coupling in the ${}^{15}\text{N}$ dimension, the total time of transverse ${}^{15}\text{N}$ magnetization ($t_{2,\text{max}} + 2\delta + 2T$) should be kept smaller than $1/J_{\text{NC}'}$. To facilitate the analysis of the HNCA and HN(CO)CA spectra the values for the sweep widths and offset in these two experiments should be identical.

References:

An efficient 3D NMR technique for correlating the proton and ${}^{15}\text{N}$ backbone amide resonances with the α -carbon of the preceding residue in uniformly ${}^{15}\text{N}/{}^{13}\text{C}$ enriched proteins. A. Bax, M. Ikura. J.Biomol. NMR. 1 (1991) 99-104.

Pulse Sequence Diagram

7.2




Fig. 24. HN(CO)CA.

 $\phi_2 = x, -x$ +States-TPPI (t_1) $\phi_3 = x, x, -x, -x$ $\phi_4 = -y, -y, y, y$ +echo-antiecho(t_2) $\phi_{rec} = x, -x, -x, x$ $\phi_1 = x, -x$ $\delta = 2.3$ ms $\beta = 4.0$ ms $\gamma = 5.5$ ms $T = 12.0$ ms $\epsilon =$ gradient+recovery

Above the arrows indicate the points of shifting the ^{13}C offset and the striped pulses are off-resonance. All carbon pulses are selective.

The sequence starts with an INEPT transfer of the magnetization from the backbone ^1H to the ^{15}N . This is followed by a delay T , during which the antiphase term is refocused and simultaneously the 1-bond coupling between ^{15}N and carbonyl- ^{13}C creates an antiphase coherence $N_y C'_z$. It is converted to $N_z C'_y$ by the pair of ^{13}C and ^{15}N 90° -pulses, and subsequently correlated with the $^{13}\text{C}^\alpha$ spins in an INEPT-transfer during the delays β , relying on the 1-bond coupling between the carbonyl and alpha carbons. The coherence is labeled by the $^{13}\text{C}^\alpha$ chemical shift during the t_1 -evolution time, and subsequently converted back to $N_y C'_z$. This rephases to an ^{15}N -coherence during the delay T and becomes simultaneously labeled by the ^{15}N -chemical shift during the constant time-type of t_2 -evolution. Overlaying of the rephasing and evolution delays minimizes transverse relaxation losses. A reverse INEPT transfer with sensitivity enhancement converts the coherence of interest back to the amide protons for observation.

Protons are decoupled by broadband DIPSII-2 scheme as in the HNCA experiment. The ^{13}C -labeled carbonyl carbons are decoupled by a selective 180° -pulse during the $^{13}\text{C}^\alpha$ -evolution. The ^{13}C -labeled carbonyl and alpha carbons are decoupled by selective 180° -pulses during the ^{15}N -CT-evolution.

Sample: 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).

Experiment time: 4h

Table 16: Acquisition Parameters

Parameter	Value	Comments
PULPROG	hncocagp3d.2	pulse program
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimize with „gs“)
NUC2	^{13}C	nucleus on f2 channel
O2P	52 ppm	offset ^{13}C (centered on C^α region)
NUC3	^{15}N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse f2
PL1		high power level f1 channel
PL2	120 dB	high power level f2 channel (not used)
PL3		high power level f3 channel
PL16		power level for GARP decoupling, f3 channel
PL19		power level for DIPSI-2 decoupling f1 channel
SP1		power level, selective H_2O flip-back optimize in „gs“-mode!
SP2		power level, 90° shaped pulse (on res.) f2
SP3		power level, 180° shaped pulse (on res.) f2
SP5	=sp3	power level, 180° shaped pulse (C^α off res.) f2
SP7	=sp3	power level, 180° shaped pulse (C' off res.) f2
SP8	=sp2	power level, time reversed 90° (on res.) f2
SPNAM1	Sinc1.1000	selective H_2O flip-back
SPNAM2	G4.256	90° shape: 4 Gaussian cascade
SPNAM3	Q3.256	180° shape: 3 Gaussian cascade
SPNAM5	Q3.256	180° shape: 3 Gaussian cascade
SPNAM7	Q3.256	180° shape: 3 Gaussian cascade
SPNAM8	G4tr.256	90° shape: time reversed 4 Gaussian cascade
SPOFF1	0	selective H_2O flip-back offset
SPOFF2	0	
SPOFF3	0	
SPOFF5	-18270 Hz	<= minus 121 ppm (^{13}C) in Hz (=173-52 ppm)

HN(CO)CA

Parameter	Value	Comments
SPOFF7	+18270 Hz	<= plus 121 ppm (¹³ C) in Hz (=173-52 ppm)
SPOFF8	0	
P1		90° ¹ H pulse f1 channel
P2	preset to 2*p1	180° ¹ H pulse f1 channel
P3		90° ¹³ C pulse f2 channel
P4	preset to 2*p3	180° ¹³ C pulse f2 channel
P11	1000-2000u	selective H ₂ O flip-back if needed!
P13	409u	90° selective pulse f2 (333μs on 700MHz)
P14	256u	180° selective pulse f2
P16	1000u	gradient pulse length
P26	=pcpd1	90° DIPSI-2 pulse at pl19
P21		90° ¹⁵ N pulse f3 channel
P22	preset to 2*p21	180° ¹⁵ N pulse f3 channel
PCDP1		decoupling pulse DIPSI-2 channel f1 (¹ H)
CPDPRG1	dipsi2	decoupling program channel f1 (¹ H)
PCDP3		decoupling pulse GARP channel f3 (¹⁵ N)
CPDPRG3	garp	decoupling program channel f3 (¹⁵ N)
D1	1.0s	recycle delay
D16	50u	gradient recovery
GPZ1	30	#1 gradient amplitude
GPZ2	80	#2 gradient amplitude
GPZ3	8.1	#3 gradient amplitude
GPNAM1	SINE.100	#1 gradient shape
GPNAM2	SINE.100	#2 gradient shape
GPNAM3	SINE.100	#3 gradient shape
NS	4	number of scans
DS	16	dummy scans (x*2*ns because of E/A)
FQ2LIST	hncoca13c	see the instruction below
RG	2048 or more	receiver gain
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
NUCLEI	¹ H	
F2 indirect ¹⁵N	*****	(middle column)
TD	40	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	2	no of in10 in pulse program

Parameter	Value	Comments
NUCLEI	¹⁵ N	
IN30	= IN10	Constant time: size of the increment and decrement should be equal
F1 indirect ¹³ C	*****	(rightmost column)
TD	64	number of real points
SW	40 ppm	sweep width indirect ¹³ C
ND0	2	no of in0 in pulse program
NUCLEI	¹³ C	

Creating the ¹³C frequency list

In this experiment the ¹³C frequency is initially at the carbonyl-carbon frequency of 173 ppm. When the magnetization has been transfer from ¹³C' to ¹³C^α the frequency is changed to the value of 52 ppm in the middle of the alpha-carbon region. For this purpose you need to create the frequency list. Edit a file /u/exp/stan/nmr/lists/f1/hncoca13c with the following contents:

```
O 150.902749
26108.18
8008.18
```

The first line contains the capital letter „O“ and the basic frequency of ¹³C on your instrument. The second line the carbonyl offset in Hz (given by the o2 parameter). The third line corresponds to the alpha-carbon offset used for the spoff5, in the above example 26108.18-18100=8008.18 Hz. Note that in this case the o2p value only sets the frequency in the spectral display.

Spectral Processing

7.4


Table 17: Processing parameters

Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	2k	zero fill to 2048 complex points
WDW	QSINE	squared sine bell window function, e.g.
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	qpol	water deconvolution

HN(CO)CA

Parameter	Value	Comment
BCFW	0.5-1.0	range for water deconvolution, in ppm
F2 indirect ¹⁵ N	*****	(middle column)
SI	128	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell window function
SSB	2	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction
REVERSE	TRUE	
F1 indirect ¹³ C	*****	(rightmost column)
SI	256	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell window function
SSB	2	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction

Fig. 25. The ¹H-¹³C and ¹H-¹⁵N planes of the HN(CO)CA spectrum.




HCCH-TOCSY

8

Introduction

8.1

The HCCH-TOCSY experiment renders the ^1H - ^1H correlations within a spin system (the side chain protons of each residue in a polypeptide, e. g.) while separating the signals in the third dimension according to the ^{13}C frequency. Note however, that here the TOCSY-transfer is accomplished through a ^{13}C - ^{13}C mixing, instead of ^1H - ^1H mixing as in the 2D TOCSY experiment. For slowly tumbling macromolecules the 2D homonuclear experiment has lower sensitivity because, first, the natural ^1H linewidth ($=1/\pi T_2$) becomes larger than the ^1H - ^1H couplings and, second, ^{13}C labeling reduces the proton transverse relaxation time T_2 . The HCCH-TOCSY experiment relies on the strong one-bond ^1H - ^{13}C (125-250Hz) and ^{13}C - ^{13}C (35-55Hz) couplings for the magnetization transfer. The flow of the magnetization is as follows



A complementary experiment to HCCH-TOCSY is the HCCH-COSY experiment, which is based on transfer over one ^{13}C - ^{13}C bond and correlates neighbouring ^1H frequencies only. Both experiments are performed on ^{13}C -enriched samples. Thus, for protein studies, the HCCH-TOCSY and HCCH-COSY experiments are chosen to obtain side chain ^1H and ^{13}C assignments. To monitor through-space connectivities within and between side chains a three dimensional ^{13}C -correlated NOESY experiment can be recorded. Concerning the spectral analysis, the three experiments yield data in a similar format, which however, differs from that of the triple resonance experiments which separate the signals according to the backbone amide frequencies.


References:

A gradient-enhanced HCCH-TOCSY experiment for recording sidechain ^1H and ^{13}C correlations in H_2O samples of proteins. L. E. Kay, G.-Y. Xu, A. U. Singer, D. R. Muhandiram, J. D. Forman-Kay. JMR 101 B (1993) 333-337.

^1H - ^1H correlation via isotropic mixing of ^{13}C magnetization, a new three-dimensional approach for assigning ^1H and ^{13}C spectra of ^{13}C -enriched proteins. A. Bax, M. Clore, A. M. Gronenborn. J. Magn. Res. 88 (1990) 425-431.

Pulse Sequence Diagram

8.2



$$\begin{aligned} \phi_1 &= x, -x + \text{States-TPPI} (t_1) & \phi_2 &= 2(x), 2(-x) + \text{States-TPPI} (t_2) & \phi_3 &= 4(x), 4(-x) \\ \phi_{\text{rec}} &= x, -x, -x, x & \Delta &= 1/(4 * J_{\text{HC}}) & \epsilon &= 1/(6 * J_{\text{HC}}) & 2 * \gamma &= 1/(6 * J_{\text{HC}}) \end{aligned}$$

Fig. 26. HCCH-TOCSY.

In the above diagram 90°-pulses are denoted by thin bars and 180°-pulses are denoted by thick bars. At the carbonyl frequency (¹³C') two 180° shaped pulses are applied. The pulse phases are x if not specified.

The pulse sequence starts with proton magnetization. After the ¹H-evolution time t_1 , the magnetization is transferred to the carbon-nuclei through the INEPT scheme. The antiphase ¹³C-magnetization is refocused during the intervals $2 * \gamma$ (the length of which is a compromise for the different multiplicities), and labeled by the ¹³C-chemical shift during the t_2 -evolution time. The selective 180° pulses at the carbonyl frequency are applied in order to refocus C^α-C' coupling. The trim pulses at the beginning and the end of the mixing defocus all in-phase ¹³C-magnetization that is not parallel to the effective field.

The length of the isotropic mixing is regulated by the loop parameter l1. Assuming a ¹³C-mixing pulse of 25us, the l1-value of 1 gives a mixing time of 5.4ms (resulting in transfer over one bond), 2 gives 10.9 ms (transfer up to two bonds), 3 gives 16.3 ms (transfer up to three or more bonds). Finally, a reverse INEPT sequence brings the carbon magnetization back to the attached protons for detection. Water suppression is achieved by two z-filters and a pair of proton spin lock pulses.

Sample: 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).

Experiment time: 16h

Table 18: Acquisition Parameters

Parameter	Value	Comments
PULPROG	hcchdigp3d	pulse program
NUC1	^1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimized with „gs“)
NUC2	^{13}C	nucleus on f2 channel
O2P	43 ppm	offset ^{13}C
NUC3	^{15}N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse
PL1		^1H high power level f1 channel
PL2		^{13}C high power level f2 channel
PL3		^{15}N high power level f3 channel
PL11		power level for ^1H trim pulses (=pl1+3dB). Edit pulse sequence!
PL12		power level GARP decoupling f2
PL15		power level C-C cross polarization f2
SP5		180° off-resonance shaped pulse for carbony decoupling (C')
SPNAM5	Q3.256	shape: Gaussian cascade
SPOFF5		= 130 ppm (^{13}C) in Hz(=173-43 ppm)
P1		90° ^1H pulse f1 channel
P2		180° ^1H pulse f1 channel
P3		90° ^{13}C pulse f2 channel
P4		180° ^{13}C pulse f2 channel
P9	25u	low power pulse (cross polarization)
P14	256u	180° off-resonance shaped pulse for carbonyl decoupling (C')
P17	1m	^1H trim pulse f1 channel
P27	2m	^{13}C trim pulse f2 channel
P22		180° ^{15}N pulse f3 channel
PCDP2		decoupling pulse on ^{13}C f2
CPDPRG2	garp	decoupling program on ^{13}C f2
D1	1.5 s	recycle delay

HCCH-TOCSY

Parameter	Value	Comments
D4	preset 1.6m	= 1/4 (J_{CH})
D21	preset 1.1m	= 1/6 (J_{CH})
D23	preset 475u	2*d23~ 1/6 (J_{CH})
D16	50u	gradient recovery
P16	500u	length of #1 gradient
P19	2m	length of #2 gradient
P29	300u	length of #3 gradient
P30	5m	length of #4 gradient
P31	4.4m	length of second #4 gradient
GPZ1	16	#1 gradient amplitude
GPZ2	16	#2 gradient amplitude
GPZ3	30	#3 gradient amplitude
GPZ4	60	#4 gradient amplitude
GPNAM1	SINE.50	#1 gradient shape
GPNAM2	SINE.32	#2 gradient shape
GPNAM3	SINE.100	#3 gradient shape
GPNAM4	SINE.100	#3 gradient shape
NS	4	number of scans
DS	128	number of dummy scans
l1	1, 2 or 3	DIPSI mixing time, see text. MAX 3!
F3 1H acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
NUCLEI	1H	
F2 indirect ^{13}C	*****	(middle column)
TD	64	number of real points
SW	75 ppm	sweep width indirect ^{13}C
ND10	2	no of in10 in pulse program
NUCLEI	^{13}C	
F1 indirect 1H	*****	(rightmost column)
TD	128	number of real points
SW	14 ppm	sweep width indirect 1H
ND0	2	no of in0 in pulse program
NUCLEI	1H	

Spectral Processing


8.4

Table 19: Processing parameters


Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	4k	zero fill
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	no automatic baseline correction
F2 indirect ¹³C	*****	(middle column)
SI	256	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	no automatic baseline correction
F1 indirect ¹H	*****	(rightmost column)
SI	1k	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	no	phase correction
BC_mod	no	no automatic baseline correction

HCCH-TOCSY

Fig. 27. The ^1H - ^1H (I1=3) and ^1H - ^{13}C planes of the HCCH-TOCSY spectrum.



The ^1H - ^1H (I1=2) and ^1H - ^1H (I1=1) planes of the HCCH-TOCSY spectrum.




CP-driven HCCH-TOCSY

9

Introduction

9.1

The heteronuclear cross polarization (CP) driven HCCH-TOCSY experiment renders identical spectrum as the INEPT-driven HCCH-TOCSY presented above, albeit with higher sensitivity. Thus the ^1H - ^1H correlations within a spin system are observed while dispersing the signals in the third dimension according to the ^{13}C frequency. The flow of the magnetization is



The initial proton magnetization is labeled by the ^1H -frequency during t_1 -evolution. Heteronuclear cross polarization is employed to transfer the magnetization to the attached carbon-nuclei. To avoid phase distortions the heteronuclear mixing should be performed with a pulse train consisting of 180° phase shifts such as the DIPSI scheme. The resulting ^{13}C -magnetization is frequency labeled during t_2 , subsequently transferred to other carbons within the spin system through homonuclear mixing, and eventually converted back to proton magnetization for detection.

The reason for the sensitivity gains obtained by heteronuclear cross polarization (CP) compared to the INEPT transfer with free precession delays, is the minimization of the residence time of the ^{13}C -magnetization in the transverse plane. This becomes important particularly for larger proteins where the ^{13}C T_2 relaxation is relatively fast.

Water suppression in this sequence is achieved by two strong gradients in the following way. The first gradient dephases the water magnetization in the x,y-plane which has been created by the first ^1H pulse and spinlocked by the CP. The water magnetization that has returned to the z-axis through relaxation, is flipped to the transverse plane by the second 90° ^1H pulse and dephased by the second gradient.


References:

HCCH-TOCSY spectroscopy of ^{13}C -labeled proteins in H_2O using heteronuclear cross-polarization and pulsed-field gradients. H. Wang, E. R. P. Zuiderweg. *J. Biomol. NMR* 5 (1992) 207-211.

Sensitivity improvement in 2D and 3D HCCH-spectroscopy using heteronuclear cross-polarization. A. Majumdar, H. Wang, R. C. Morshauer, E. R. P. Zuiderweg. *J. Biomol. NMR* 3 (1993) 387-397.

Pulse Sequence Diagram

9.2



$$\phi_1 = x, -x + \text{States-TPPI} (t_1)$$

$$\phi_2 = 2(y), 2(-y) + \text{States-TPPI} (t_2)$$

$$\phi_{\text{rec}} = x, -x, -x, x$$

Fig. 28. CP-driven HCCH-TOCSY.

In the above diagram 90° -pulses are denoted by thin bars and 180° -pulses are denoted by thick bars. The pulse phases are x if not specified.

In principle, this sequence accomplishes coherence selection with a single transient.

The optimal transfer time for the heteronuclear cross polarization is the inverse of the coupling constant, that is, $1/J_{\text{CH}}$ between ^1H - ^{13}C . The value of 6 ms that is used here, is an average value for the different carbon multiplicities (the one-bond J_{CH} coupling constant is 125 Hz in a CH-group, 156 Hz in CH_2 and 249 Hz in CH_3), see Appendix 12.1. For effective polarization transfer the simultaneous RF fields applied to the two unlike nuclei need to fulfill the Hartman-Hahn condition, $\gamma_{\text{H}}B_{1\text{H}} = \gamma_{\text{C}}B_{1\text{C}}$, where γ_{H} and γ_{C} are the magnetogyric ratios and $B_{1\text{H}}$ and $B_{1\text{C}}$ are the field amplitudes. For the present purpose, $\gamma B_1/2\pi = 10$ kHz field suffices, requiring a 25 μs mixing pulse on both ^1H and ^{13}C nuclei.

The length of the ^{13}C - ^{13}C mixing is set through the loop parameter l1. Assuming a ^{13}C mixing pulse of 25 μs , the l1-value of 1 gives a mixing time of 5.4 ms (resulting in transfer over one bond), l1=2 corresponds to 10.9 ms (transfer up to two bonds) and l1=3 corresponds to 16.3 ms (transfer up to three or more bonds).

Data Acquisition**9.3****Sample:** 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).**Experiment time:** 16h**Table 20: Acquisition Parameters**

Parameter	Value	Comments
PULPROG	hcchcpdi3d	pulse program
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimized with „gs“)
NUC2	13C	nucleus on f2 channel
O2P	43 ppm	offset ^{13}C
PL0	120 dB	power level preceding shaped pulse
PL1		^1H high power level f1 channel
PL2		^{13}C high power level f2 channel
PL12		^{13}C power level GARP decoupling f2
PL15 (not pl16!)		power level ^{13}C - ^{13}C cross polarization
SP15		180° off-resonance shaped pulse (C') f
SPNAM15	Q3.256	shape: 3 Gaussian cascade
SPOFF15		≤ 130 ppm (^{13}C) in Hz (=173-43 ppm)
P1		90° ^1H pulse f1 channel
P2	2*p1 preset	180° ^1H pulse f1 channel
P3		90° ^{13}C pulse f2 channel
P4	2*p3 preset	180° ^{13}C pulse f2 channel
P9	25u	low power pulse (cross polarization) on both ^1H and ^{13}C
P27	1m	^{13}C trim pulse f2 channel
PCDP1	256u	select. off-resonance decoupling pulse
CPDPRG1	mlevsp180	select. off-resonance decoupling (C') f2
PCPD2		decoupling pulse on ^{13}C , f2 channel
CPDPRG2	garp	broadband decoupling on ^{13}C
D1	1.5 s	recycle delay
D16	50u	gradient recovery
D18	3m, preset	length of #1 gradient
D19	1m, preset	length of #2 gradient
GPZ1	80	#1 gradient amplitude
GPZ2	80	#2 gradient amplitude
GPNAM1	SINE.100	#1 gradient shape

Parameter	Value	Comments
GPnam2	SINE.100	#2 gradient shape
NS	4	number of scans
DS	128	number of dummy scans
I1	1	^1H - ^{13}C DIPSImixing time, see text
I2	1, 2 or 3	^{13}C - ^{13}C DIPSImixing time, see text
F3 ^1H acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	12 ppm	sweep width
NUCLEI	^1H	
F2 indirect ^{13}C	*****	(middle column)
TD	64	number of real points
SW	75 ppm	sweep width indirect ^{13}C
ND10	2	no of in10 in pulse program
NUCLEI	^{13}C	
F1 indirect ^1H	*****	(rightmost column)
TD	600	number of real points
SW	12 ppm	sweep width indirect ^1H
ND0	2	no of in0 in pulse program
NUCLEI	^1H	

Spectral Processing

9.4

Table 21: Processing parameters

Parameter	Value	Comment
F3 acquisition	*****	(leftmost column)
SI	2k	zero fill
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	no automatic baseline correction

CP-driven HCCH-TOCSY

Parameter	Value	Comment
F2 indirect ¹³C	*****	(middle column)
SI	256	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	pk	phase correction
PHC0	0	zero order phase correction
PHC1	0	first order phase correction
F1 indirect ¹H	*****	(rightmost column)
SI	1k	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell
SSB	3	60° shifted
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction

Fig. 29. The CP-driven HCCH-TOCSY spectrum. (See Fig. 27. HCCH-TOCSY)

H(CC)(CO)NH

10

Introduction

10.1





Fig. 30. The magnetization transfer in H(CC)(CO)NH.

The H(CC)(CO)NH experiment correlates the backbone amide proton and nitrogen frequencies of one residue with the side chain proton frequencies of the preceding residue. Thus the spectral layout is similar to the ^{15}N -TOCSY experiment and it can be conveniently analyzed together with the ^{15}N -NOESY HSQC and HAHB(CO)NH experiments. The experimental sensitivity is, however, lower. The flow of the magnetization is the following



References:

Heteronuclear multidimensional NMR experiments for the structure determination of proteins in solution employing pulsed field gradients. M. Sattler, J. Schleucher & C. Griesinger. Prog. NMR Spectr. 34 (1999) 93-158.

Correlation of Backbone amide and aliphatic side-chain resonances in $^{13}\text{C}/^{15}\text{N}$ -enriched proteins by isotropic mixing of ^{13}C magnetization. S. Grzesiek, J. Anglister, A. Bax. J. Magn. Res. B 101 (1993) 114-119.

Pulse Sequence Diagram

10.2




Fig. 31. H(CC)(CO)NH.

The arrow indicates where the ^{13}C -frequency is moved from C^α to C' and the ^1H -offset from the water frequency to the amide region. The carbonyl 180° -pulse denoted by BS compensates for the Bloch-Siegert phase shift caused by the preceding pulse. The striped carbon pulses are off-resonance.

The sequence starts with proton magnetization that is labeled by the ^1H -chemical shift labeled during the t_1 -evolution period. Subsequently the magnetization is transferred to the attached carbons through a heteronuclear cross polarization of 6 ms. Simultaneously the magnetization transfer between the coupled carbon nuclei through homonuclear mixing begins. Its total length is preset in the pulse program to 12 ms. This is followed by the delay 2β in order to relay magnetization to the alpha-carbons. During the long delay 2ζ the C^α -coherence is relayed further to the carbonyl-carbons and at the same time defocusing with respect to the amide- ^{15}N takes place. The magnetization is converted to a ^{15}N -coherence by the first ^{15}N - 90° pulse and labeled by the nitrogen chemical shift in a constant time fashion. Finally the magnetization is refocused to an amide proton coherence for detection.

Sample: 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).


Experiment time: 1 day 20h

Table 22: Acquisition Parameters

Parameter	Value	Comments
PULPROG	hcccconhgp3d1	pulse program
NUC1	^1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (or 3.0 ppm, see note below)
NUC2	^{13}C	nucleus on f2 channel
O2P	40 ppm	offset ^{13}C (centered on aliphatic ^{13}C region)
NUC3	^{15}N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse f2
PL1		high power level f1 channel (^1H)
PL2	120 dB	high power level f2 channel (not used)
PL3		high power level f3 channel (^{15}N)
PL10		low power level for cross-polarization (^1H)
PL15		low power level for cross-polarization (^{13}C)
PL16		power level for GARP dec. f3 channel (^{15}N)
PL19		power level for DIPSI-2 decoupling f1 (^1H)
SP2		power level, 90° shaped pulse (^{13}C on res.)
SP3		power level, 180° shaped pulse (^{13}C on res.)
SP5	=sp3	power level, 180° shaped pulse (C' off res.)
SP7	=sp3	power level, 180° shaped pulse (C^α off res.)
SP8	=sp2	power level, time reversed 90° (^{13}C on res.)
SP9		power level, very selective 180° (C^α on res.)
SPNAM2	G4.256	90° shape: 4 Gaussian cascade
SPNAM3	Q3.256	180° shape: 3 Gaussian cascade
SPNAM5	Q3.256	180° shape: 3 Gaussian cascade
SPNAM7	Q3.256	180° shape: 3 Gaussian cascade
SPNAM8	G4tr.256	90° shape: time reversed 4 Gaussian casc.
SPNAM9	Q3.256	180° shape: 3 Gaussian cascade
SPOFF2	0	
SPOFF3	0	
SPOFF5	20000 Hz	= plus 133 ppm (^{13}C) in Hz (=173-40 ppm)
SPOFF7	-20000 Hz	= minus 133 ppm (^{13}C) in Hz (=173-40 ppm)
SPOFF8	0	

H(CC)(CO)NH

Parameter	Value	Comments
SPOFF9	0	
P1		90° ¹ H pulse f1 channel
P2		180° ¹ H pulse f1 channel
P6		90° ¹ H cross-polarization pulse f1 channel
P9	=p6	90° ¹³ C cross-polarization pulse f2 channel
P13	409u	90° selective pulse f2 (333μs at 700 MHZ)
P14	256u	180° selective pulse f2 channel
P16	500u	gradient pulse length
P17	2000u	trim pulse at low power level
P21		90° ¹⁵ N pulse f3 channel
P22		180° ¹⁵ N pulse f3 channel
P24	768u	180° strongly selective pulse f2 channel
PCDP1		decoupling pulse length DIPSI-2 channel f1
CPDPRG1	dipsi2	decoupling program channel f1 (¹ H)
PCDP3		decoupling pulse length GARP channel f3
CPDPRG3	garp	decoupling program channel f3 (¹³ C)
D1	1.2 s	recycle delay
D16	50u	gradient recovery
GPZ1	80	#1 gradient amplitude
GPZ2	8.1	#2 gradient amplitude
GPNAM1	SINE.100	#1 gradient shape
GPNAM2	SINE.100	#2 gradient shape
NS	32	number of scans
DS	64	dummy scans (x*2*ns because of E/A)
FQ1LIST1	hccconh1h	see the instruction below
FQ2LIST1	hccconh13c	see the instruction below
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
NUCLEI	¹ H	
F2 indirect ¹⁵N	*****	(middle column)
TD	40	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	2	no of in10 in pulse program
NUCLEI	¹⁵ N	
IN30	= IN10	Constant time: size of the increment and decrement should be equal



Parameter	Value	Comments
F1 indirect ^1H	*****	<i>(rightmost column)</i>
TD	160	number of real points
SW	14 ppm	sweep width indirect ^1H
ND0	2	no of in0 in pulse program
NUCLEI	^1H	

Creating the ^1H and ^{13}C frequency lists


It is recommended to set the proton frequency at the water frequency for the first part of the experiment, and shift it in the middle of amide region for the second part. The former action is taken in order to achieve better water suppression and the latter for higher sensitivity. The ^{13}C frequency should be set in the middle of the aliphatic region during the cross polarization. It should be moved to the carbonyl region when the magnetization is transferred to these nuclei. How to create a frequency list, see Section 7.3. The frequency lists are stored in the directory /u/exp/stan/nmr/lists/f1.

Spectral Processing

10.4


Table 23: Processing parameters

Parameter	Value	Comment
F3 acquisition		<i>(leftmost column)</i>
SI	2k	zero fill to 4096 complex points
WDW	QSINE	squared sine bell window function, e.g.
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	qpol	water deconvolution
BCFW	1.0	range for water deconvolution
STSR	0	display only the left half of spectrum
STSI	1k	-, -
F2 indirect ^{15}N	*****	<i>(middle column)</i>
SI	512	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_MOD	no	no phase correction is needed
REVERSE	true	



Parameter	Value	Comment
F1 indirect ^1H	*****	(rightmost column)
SI	1k	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	no phase correction is needed

Fig. 32. The ^1H - ^1H and ^1H - ^{15}N planes of the H(CC)(CO)NH spectrum.



CBCA(CO)NH

11

Introduction

11.1





Fig. 33. The magnetization transfer in CBCA(CO)NH.

The CBCA(CO)NH experiment correlates the backbone amide proton and nitrogen frequencies of one residue with the alpha- and beta-carbon frequencies of the preceding residue. Thus the spectrum contains the same information as HN(CO)CA. The CBCA(CO)NH experiment is less sensitive, but provides in addition the $^{13}\text{C}^{\beta}$ frequencies which are instrumental for identification of residue type and secondary structure, and correlation of side chains to the sequential assignment. The flow of the magnetization is following




References:

Gradient-enhanced triple-resonance three-dimensional NMR experiments with improved sensitivity. D. R. Muhandiram, L. E. Kay. J. Magn. Res. B 103 (1994) 203-216.

Correlating backbone amide and side chain resonances in larger proteins by multiple relayed triple resonance NMR. S. Grzesiek, A. Bax. J. Am. Chem. Soc. 114 (1992) 6291-6293.

Pulse Sequence Diagram

11.2



$$\phi_1 = x, -x + \text{States-TPPI} (t_1)$$

$$\phi_3 = 2x, 2(-x)$$

$$\phi_4 = 2(y), 2(-y)$$

$$\phi_2 = -4x, 4(-x) + \text{echo-antiecho}(t_2)$$

$$\phi_{\text{rec}} = 2(x, -x, -x, x)$$

$$\alpha = 1.7 \text{ ms}$$

$$\beta = 1.1 \text{ ms}$$

$$\zeta = 3.6 \text{ ms}$$

$$\delta = 4.4 \text{ ms}$$

$$\epsilon = 12.4 \text{ ms}$$

$$\nu = 5.5 \text{ ms}$$

$$\gamma = 2.3 \text{ ms}$$

$$\eta = \text{gradient+recovery}$$

$$T_C = 3.6 \text{ ms}$$

$$T = 12.4 \text{ ms}$$

$$\text{arrow} = {}^{13}\text{C-frequence moves from } C^\alpha \text{ to } C'$$

Fig. 34. CBCA(CO)NH.

The sequence starts with an INEPT enhancement of the magnetization of the aliphatic carbons. The ${}^{13}\text{C}$ magnetization undergoes chemical shift t_1 -evolution during the constant time period $2 \cdot T_C$ (in order to avoid line broadening due to ${}^{13}\text{C}$ - ${}^{13}\text{C}$ couplings). The total length of the constant time is $\sim 1/(4J_{CC}) = 7.2 \text{ ms}$. Scalar coupling to the carbonyls is removed by a selective 180° pulse. Note that the carbonyl 180° pulses denoted by BS are applied in order to remove the phase shifts caused by the preceding 180° -pulse. The third 90° carbon pulse converts ${}^{13}\text{C}^\beta$ magnetization to ${}^{13}\text{C}^\alpha$ magnetization. The existing ${}^{13}\text{C}^\alpha$ magnetization is not affected. During the subsequent delay of 2ζ the ${}^{13}\text{C}^\alpha$ -magnetization is conferred to the carbonyls. During the delay δ the carbonyl magnetization refocuses with respect to its attached ${}^{13}\text{C}^\alpha$ spin, and during the delay 2ϵ it is transferred to the adjacent ${}^{15}\text{N}$ spin. During the second constant time period $2 \cdot T_N$ ($=24.8 \text{ ms}$) the antiphase ${}^{15}\text{N}$ magnetization rephases with respect to the carbonyl and becomes labeled by the ${}^{15}\text{N}$ -chemical shift. The ${}^{15}\text{N}$ -magnetization dephases with respect to the attached proton during the final fraction ν , when ${}^1\text{H}$ decoupling is turned off. At the end the magnetization is transferred to the observable proton coherences through a sensitivity enhanced reverse INEPT sequence. Coherent ${}^1\text{H}$ decoupling is interrupted during the gradients.

Sample: 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).

Experiment time: 15h

Table 24: Acquisition Parameters

Parameter	Value	Comments
PULPROG	cbcaconhgp3d.2	pulse program
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (optimize with „gs“)
NUC2	13C	nucleus on f2 channel
O2P	43 ppm	offset ^{13}C (centered on $\text{C}^\alpha/\text{C}^\beta$ region)
NUC3	15N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse f2
PL1		high power level f1 channel
PL2	120 dB	high power level f2 channel (not used)
PL3		high power level f3 channel
PL16		power level for GARP decoupling f3 channel
PL19		power level for DIPSI-2 decoupling ^1H
SP2		power level, 90° shaped pulse (^{13}C on res.)
SP3		power level, 180° shaped pulse (^{13}C on res.)
SP5	=sp3	power level, 180° shaped pulse (C' off res.)
SP7	=sp3	power level, 180° shaped pulse (C^α off res.)
SP8	=sp2	power level, time reversed 90° (^{13}C on res.)
SPNAM2	G4.256	90° shape: 4 Gaussian cascade
SPNAM3	Q3.256	180° shape: 3 Gaussian cascade
SPNAM5	Q3.256	180° shape: 3 Gaussian cascade
SPNAM7	Q3.256	180° shape: 3 Gaussian cascade
SPNAM8	G4tr.256	90° shape: time reversed 4 Gaussian casc.
SPOFF2	0	
SPOFF3	0	
SPOFF5	20000 Hz	\leq plus 130 ppm (^{13}C) in Hz (=173-43 ppm)
SPOFF7	-20000 Hz	= minus 130 ppm (^{13}C) in Hz (=173-43 ppm)
SPOFF8	0	
P1		90° ^1H pulse f1 channel
P2	preset to 2*p1	180° ^1H pulse f1 channel

CBCA(CO)NH

Parameter	Value	Comments
P3		90° ¹³ C pulse f2 channel
P4	preset to 2*p3	180° ¹³ C pulse f2 channel
P13	409u	90° selective pulse f2 (333μs at 700 MHz)
P14	256u	180° selective pulse f2
P16	500u	#1 gradient pulse length
P19	1000u	#2 gradient pulse length
P21		90° ¹⁵ N pulse f3 channel
P22	preset to 2*p21	180° ¹⁵ N pulse f3 channel
P29	250u	#3 gradient pulse length
P30	1500u	#4 gradient pulse length
P31	2500u	#5 gradient pulse length
PCDP1		decoupling pulse length DIPSI-2 channel f1
CPDPRG1	dipsi2	decoupling program channel f1
PCDP3		decoupling pulse length GARP channel f3
CPDPRG3	garp	decoupling program channel f3
D1	1.0 s	recycle delay
D16	50u	gradient recovery
GPZ1	30	#1 gradient amplitude
GPZ2	80	#2 gradient amplitude
GPZ3	8.1	#3 gradient amplitude
GNAM1	SINE.100	#1 gradient shape
GNAM 2	SINE.100	#2 gradient shape
GNAM 3	SINE.100	#3 gradient shape
NS	16	number of scans
DS	32	dummy scans (x*2*ns because of E/A)
FQ2LIST1	cbcaconh13c	see the instruction below
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points
SW	14 ppm	sweep width
NUCLEI	¹ H	
F2 indirect ¹⁵N	*****	(middle column)
TD	36	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	2	no of in10 in pulse program
NUCLEI	¹⁵ N	
IN30	= IN10	constant time: decrement should be equal to increment

Parameter	Value	Comments
F1 indirect ^{13}C	*****	<i>(rightmost column)</i>
TD	128	number of real points
SW	75 ppm	sweep width indirect ^{13}C
ND0	2	no of in0 in pulse program
NUCLEI	^{13}C	
IN20	= IN0	constant time decrement (= increment)

Creating the ^{13}C frequency list

In this experiment the ^{13}C frequency is initially in the middle of the alpha- and beta-carbon regions, at 43 ppm. When the magnetization has been transferred from $^{13}\text{C}^{\alpha}$ to C' the frequency is also changed to the carbonyl-value of 173 ppm. For this purpose you need to create the frequency list „cbcaconh13c“. See Section 7.3. Note that the frequency lists are stored in the directory /u/exp/stan/nmr/lists/f1.

Spectral Processing

11.4


Table 25: Processing parameters

Parameter	Value	Comment
F3 acquisition	*****	<i>(leftmost column)</i>
SI	2k	zero fill to 4096 complex points
WDW	QSINE	squared sine bell window function, e. g.
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	qpol	Removal of water residual water signal
BCFW	1.0	+/- ppm around water is deconvoluted
STSR	0	Display only the left half of spectrum
STSI	1k	- „ -
F2 indirect ^{15}N	*****	<i>(middle column)</i>
SI	256	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$

CBCA(CO)NH

Parameter	Value	Comment
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction
REVERSE	TRUE	
F1 indirect ^{13}C	*****	<i>(rightmost column)</i>
SI	1k	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction not needed

Fig. 35. The ^1H - ^{13}C and ^1H - ^{15}N planes of the CBCA(CO)NH spectrum.



HBHA(CO)NH

12

Introduction

12.1





Fig. 36. The magnetization transfer in HBHA(CO)NH.

The HAHB(CO)NH experiment correlates the backbone amide proton and nitrogen frequencies of one residue with the alpha- and beta-proton frequencies of the preceding residue. Thus it is similar to the CBCA(CO)NH experiment and it is also equally sensitive. The flow of the magnetization is the following



References:

An efficient triple-resonance experiment for proton-directed sequential backbone assignment of medium-sized proteins. A. C. Wang, P. J. Lodi, J. Qin, G. W. Vuister, A. M. Gronenborn, G. M. Clore. *J. Magn. Res. B* 105 (1994) 196-198.

Amino acid type determination in the sequential assignment procedure of uniformly ¹³C/¹⁵N-enriched proteins. S. Grzesiek, A. Bax. *J. Biomol. NMR* 3 (1993) 185-204.

Pulse Sequence Diagram

12.2




Fig. 37. HBHA(CO)NH.

The arrow indicates where the ^{13}C -frequency is moved from C^α to C' and the ^1H -offset from aliphatics to the amide region or water frequency.

The carbonyl 180° -pulse denoted by BS compensates for the phase shift caused by the preceding 180° -pulse. Striped pulses are off-resonance.

The sequence starts with proton magnetization that is labeled by the ^1H -chemical shift labeled during a semi-constant time evolution period. Simultaneously the magnetization becomes antiphase with respect to the attached $^{13}\text{C}^\alpha$ and $^{13}\text{C}^\beta$. It is refocused during delay α . Transfer from $^{13}\text{C}^\beta$ to $^{13}\text{C}^\alpha$ takes place during the delay 2β . From here on the sequence is identical to the CBCA(CO)NH experiment.

Data Acquisition**12.3****Sample:** 2 mM ^{15}N , ^{13}C -labeled Ribonuclease T1, (104 a.a.).**Experiment time:** 19h**Table 26: Acquisition Parameters**

Parameter	Value	Comments
PULPROG	hbhaconhgp3d.2	pulse program
NUC1	1H	nucleus on f1 channel
O1P	4.7 ppm	offset ^1H (or 3.0 ppm, see note below)
NUC2	13C	nucleus on f2 channel
O2P	43 ppm	offset ^{13}C (centered on $\text{C}^\alpha/\text{C}^\beta$ region)
NUC3	15N	nucleus on f3 channel
O3P	116 ppm	offset ^{15}N
PL0	120 dB	power level preceding shaped pulse f2
PL1		high power level f1 channel
PL2	120 dB	high power level f2 channel (not used)
PL3		high power level f3 channel
PL12	=sp15	power level for decoupling (C^α off res.) f2
PL16		power level for GARP decoupling f3 channel
PL19		power level for DIPSI-2 decoupling f1
SP2		power level, 90° shaped pulse (C^α on res.) f2
SP3		power level, 180° shaped pulse (C^α on res.)
SP5	=sp3	power level, 180° shaped pulse (C' off res.)
SP7	=sp3	power level, 180° shaped pulse (C^α off res.)
SP8	=sp2	power level, time reversed 90° (C^α on res.)
SP15		power level, select. decoupling (C^α off res.)
SPNAM2	G4.256	90° shape: 4 Gaussian cascade
SPNAM3	Q3.256	180° shape: 3 Gaussian cascade
SPNAM5	Q3.256	180° shape: 3 Gaussian cascade
SPNAM7	Q3.256	180° shape: 3 Gaussian cascade
SPNAM8	G4tr.256	90° shape: time reversed 4 Gaussian casc.
SPNAM15	Q3.256	180° shape for selective decoupling f2
SPOFF2	0	
SPOFF3	0	
SPOFF5	19500 Hz	= plus 130 ppm (^{13}C) in Hz (=173-43 ppm)
SPOFF7	-19500 Hz	= minus 130 ppm (^{13}C) in Hz (=173-43 ppm)
SPOFF8	0	
SPOFF15	-18000 Hz	= minus 120 ppm (^{13}C) in Hz (=173-53 ppm)

HBHA(CO)NH

Parameter	Value	Comments
P1		90° ¹ H pulse f1 channel
P2	preset to 2*p1	180° ¹ H pulse f1 channel
P3		90° ¹³ C pulse f2 channel
P4	preset to 2*p3	180° ¹³ C pulse f2 channel
P13	409u	90° selective pulse f2 (333μs at 700 MHZ)
P14	256u	180° selective pulse f2
P16	500u	#1 gradient pulse length
P19	1000u	#2 gradient pulse length
P21		90° ¹⁵ N pulse f3 channel
P22		180° ¹⁵ N pulse f3 channel
P29	250u	#5 gradient pulse length
P30	1500u	#3 gradient pulse length
P31	2500u	#4 gradient pulse length
PCDP1		decoupling pulse length DIPSI-2 channel f1
CPDPRG1	dipsi2	decoupling program channel f1
PCDP2	768u	selective off-resonance decoupling pulse f2
CPDPRG2	mlevsp180	selective off-resonance decoupling (C ^α) f2
PCDP3		decoupling pulse length GARP channel f3
CPDPRG3	garp	decoupling program channel f3
D1	1.5 s	recycle delay
D16	50u	gradient recovery
GPZ1	3	#1 gradient amplitude
GPZ2	2	#2 gradient amplitude
GPZ3	30	#3 gradient amplitude
GPZ4	60	#4 gradient amplitude
GPZ5	60.75	#5 gradient amplitude
GNAM1	SINE.50	#1 gradient shape
GNAM2	SINE.100	#2 gradient shape
GNAM3	SINE.100	#3 gradient shape
GNAM4	SINE.100	#4 gradient shape
GNAM5	SINE.32	#5 gradient shape
NS	16	number of scans
DS	32	dummy scans (x*2*ns because of E/A)
FQ1LIST1	hbhaconh1h	see the instruction below
FQ2LIST1	cbcaconh13c	see the instruction below
F3 acquisition	*****	(leftmost column)
AQ_MOD	DQD	digital quadrature detection
TD	2048	no of points

Parameter	Value	Comments
SW	14 ppm	sweep width
NUCLEI	¹ H	
F2 indirect ¹⁵ N	*****	<i>(middle column)</i>
TD	36	number of real points
SW	35 ppm	sweep width indirect ¹⁵ N
ND10	2	no of in10 in pulse program
NUCLEI	¹⁵ N	
IN30	= IN10	Constant time: size of the increment and decrement should be equal
F1 indirect ¹ H	*****	<i>(rightmost column)</i>
TD	256	number of real points
SW	12 ppm	sweep width indirect ¹³ C
ND0	2	no of in0 in pulse program
NUCLEI	¹ H	
IN20		= in0-in29
IN29		=d29 / I3

Creating the ¹H and ¹³C frequency lists

It is recommended that the proton frequency is set in the middle of aliphatic region. It is shifted to the water frequency or the amide region for the second part, the former results in better water suppression and the latter in higher sensitivity. For the ¹³C frequency shift you can use the list „cbcaconh13c“ from the previous experiment. The lists are stored in the directory /u/exp/stan/nmr/lists/f1.

Setting up the selective decoupling on the alpha-carbons


Selective decoupling is applied on the alpha-carbon region (centered at 52 ppm) during the constant time evolution of the ¹⁵N-magnetization, in order to ensure that the scalar coupling between ¹³C^α and amide ¹⁵N does not evolve. The frequency offset is given through the parameter SPOFF15 which should be given a *negative* value, in Hz, corresponding to the difference 173-52 ppm=121 ppm of ¹³C at the field you are working at. Use the **mlevsp180** decoupling scheme (CPDPRG2) with a **Q3.256** pulse (SPNAM15) and length of **768 us** (PCPD2). Two parameters, SP15 and pl12, should be set to the appropriate power level value.

Table 27: Processing parameters

Parameter	Value	Comment
F3 acquisition		<i>(leftmost column)</i>
SI	4k	zero fill to 4096 complex points
WDW	QSINE	squared sine bell window function, e.g.
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction applied
PHC0		zero order phase correction
PHC1		first order phase correction
BC_mod	no	no automatic baseline correction
F2 indirect ¹⁵N	*****	<i>(middle column)</i>
SI	256	zero fill
MC2	echo-antiecho	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	pk	phase correction
PHC0		zero order phase correction
PHC1		first order phase correction
REVERSE	TRUE	
F1 indirect ¹³C	*****	<i>(rightmost column)</i>
SI	1k	zero fill
MC2	States-TPPI	
WDW	QSINE	squared sine bell window function
SSB	3	shifting of the sine bell, $\pi/3=60^\circ$
PH_mod	no	phase correction not needed

HBHA(CO)NH

Fig. 38. The ^1H - ^1H and ^1H - ^{15}N planes of the HBHA(CO)NH spectrum.



Appendices

13

Coupling constants in polypeptides

13.1




Fig. 39. One-bond and two-bond coupling constants in the polypeptide backbone.

The coupling constants used for magnetization transfer in the most common experiments are indicated above. On the left coupling constants involving protons are shown and on the right the ones involving heteronuclei. Note that the one-bond J_{CH} in aromatic rings is around 160 Hz. Usually a compromise value of 135 Hz is used. The dashed lines and corresponding framed values indicate some relevant coupling constants over two bonds.

Reference:

Heteronuclear multidimensional NMR experiments for the structure determination of proteins in solution employing pulsed field gradients. M. Sattler, J. Schleucher & C. Griesinger. Prog. in NMR Spectroscopy 34 (1999) 93-158.


^{13}C and ^1H chemical shifts in different residues**13.2**

Fig. 40. The random coil ^1H and ^{13}C chemical shifts in the natural amino acids.

Compared to the above indicated aliphatic random coil values, the alpha-carbons shift approximately +3 ppm in α -helices and -2 ppm in β -sheets. The corresponding values for the beta-carbons are: -1 ppm in α -helices, +3 ppm in β -sheets. The alpha-carbon shifts about -2 ppm in any residue preceding a proline residue.

Above the values for a reduced cysteine are indicated, whereas in a disulfide bridged cysteine the alpha-carbon is found at 52.5 ppm and the beta-carbon 39.4 ppm. For the aspartic acid and glutamic acid values of the non-protonated sidechain carboxyl acid are reported. The values for the protonated form are 36 ppm for beta-carbon in aspartic acid and 32 ppm for gamma-carbon in glutamic acid. Proline in the trans-conformation has alpha-carbon at 63.2 ppm (in the cis-form 62.4 ppm) and beta-carbon at 31.1 ppm (the cis-form 33.1 ppm).




Fig. 41. The random coil ring ^1H and ^{13}C chemical shifts of histidine, tyrosine, phenylalanine and tryptophan.

Note:

Tryptophan: W2= δ 1, W4= ϵ 3, W5= ζ 3, W6= η 2, W7= ζ 2

References:

Carbon-13-NMR of Peptides and Proteins. O. W. Howarth & D. M. J. Lilley (1978) Prog. NMR Spectr. 12, 1-40.

NMR of proteins and nucleic acids. Wüthrich, K. (1986) Wiley, New York.


Sensitivity of triple resonance experiments**13.3**

Fig. 42. The first increment of some triple resonance experiments. The sample is 2mM double-labeled ribonuclease T1, number of scans 8. The spectra were recorded on a CryoProbe system at 600MHz. The numbers indicate the signal to noise ration for the largest peak. The HNCA is not comparable due to a very long first ^{13}C -increment.